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en Instructios for use / SELF-ETCH AND SELF-ADHESIVE DUAL RESIN CEMENT

PATIENT POPULATION

Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and there are no exclusions for specific patient groups..

INDICATIONS

- Permanent cementation of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges.
- Metal alloys adhesive opaquer.
- Porcelain repair.
- Proven Color Stability over time.

PROPERTIES

- Proclinic Expert Dual Resin Cement is a self-etch, self-adhesive dual-cure, radiopaque permanent resin cement and does not require etching, priming or bonding.
- Proclinic Expert Dual Resin Cement classification according to ISO 4049 is Type 2, Class 3.
- Proclinic Expert Dual Resin Cement has a radiopacity value above 2.5 mm Aluminium.)Radiopacity is defined as the ability of a material to be visible in x-ray photographs.)
- Proclinic Expert Dual Resin Cement forms ideal bond with metal alloys, porcelain, ceramic, dentin, enamel, amalgam.
- Proclinic Expert Dual Resin Cement has a chemical affinity for metal.
- Proclinic Expert Dual Resin Cement is a dual-cure cement – it has a self curing mechanism which can also be light cured for an immediate set and extra strength at the margins.
- Light curing Proclinic Expert Dual Resin Cement gives restorations immediate stability and easy clean up of excess cement.

- Proclinic Expert Dual Resin Cement features high retention and total margin integrity.
- Proclinic Expert Dual Resin Cement automix syringe saves application time and guarantees consistent mix.
- Proclinic Expert Dual Resin Cement is available in 3 shades: A2, White and Translucent. To see shade guide visit www.bjmlabs.com.
- Proclinic Expert Dual Resin Cement has been on the market since 2008 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CONTRAINDICATIONS

- Not to be used on exposed pulp.
- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

CAUTION!

- Proclinic Expert Dual Resin Cement contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact. If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.
- Proclinic Expert Dual Resin Cement contains materials which: sensitive to light / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.
- *Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. You may also contact our EU Authorised Representative: Obelis S.A. (contact details below).*

COMPOSITION

Essential ingredients – BASE	% wt.	Essential ingredients – CATALYST	% wt.
Bisphenol a-glycidyl methacrylate	<40%	Bisphenol a-glycidyl methacrylate	<40%
Triethylenglycoldimethacrylate	<20%	Triethylenglycoldimethacrylate	<20%
Polymerization activators	<1%	2-Hydroxyethylmethacrylate	<10%
Photoinitiators	<1%	4-Methacryloxyethyltrimellitic anhydride	<5%
2-Hydroxyethylmethacrylate	<10%	Benzoyl Peroxide	<1%
Glass filler	<40%	Glass filler	<50%
Fumed silica	<5%	Fumed silica	<5%

FILLER CONTENT

The filler system comprises glass fillers and modified silica:

40-50 weight-% or 10-20 volume-%.

The particle size of the fillers range: <6 um.

INSTRUCTIONS FOR USE

Note: Tips are for single patient use only to prevent cross-contamination between patients.

1. Remove the provisional restoration. Note: Use only temporary cements that are eugenol free. The use of eugenol containing temporary materials can inhibit the polymerization process of the Proclinic Expert Dual Resin Cement.
2. Rinse thoroughly.
3. Air dry the preparation (do not desiccate).
4. Try in the restoration to ensure proper fit.
5. Automix (When using a dual-barrel syringe with automix syringe mix tip): Bleed the syringe before installation of the automix syringe mix tip. Mount a new, automix syringe mix tip onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. Hand Mix (When using a dual-barrel syringe with Detachable Manual Adapter): a) Bleed the syringe before installation of the Manual Adapter. Mount a new manual adaptor onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. b) Using a clean spatula or mixing

stick, mix the two components into a homogenous paste.

6. Dispensing the material can be done as follows:
 - For Inlay / Onlay - Mount the supplied intra oral tip onto the mixing tip and dispense the cement directly, covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all sides.
 - For Crown - Dispense the cement directly into the restoration, covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all margins.
 - For Post - Mount the supplied XXF intra oral tip onto the mixing tip. Syringe the cement into the canal. Seat the post and vibrate slightly avoiding possibility of trapped air. Allow the cement to slowly flow from the canal space.
7. When the restoration is properly seated, remove all excess cement. Excess cement is best removed in its gel state, after reaching the material working time, with a scaler or explorer. After removal of excess cement, light cure all surfaces including margins for 20 seconds. For the material time indications, see table below.
8. Time indications: The working and setting times depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00
Light curing in a light-wave range of 390 – 470 nm. Compatible with all curing units available on the market	00:20

9. Tightly re-cap syringe after use.
10. Use finishing strips for proximal areas. Polish margins with polishing discs, cups, or points.
11. Check occlusion.

MIXING OPTIONS

Note: working and setting times complies to the mentioned above range without dependence on mixing options.

Mixing options	Ratio
Automix - 1 Automix syringe	1:1 - Equal amount of Base and Catalyst
Hand Mix - 1 Automix syringe with Detachable Manual Adapter	1:1 - Equal amount of Base and Catalyst

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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OBELIS SWISS GmbH
B.J.M. Laboratories Ltd.
PROCLINIC S.A.U.