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GB

Definition

Zinc oxide and synthetic materials temporary filling cement, self hardening, ready to use. No expansion when setting, due to its high quality synthetic materials. Eugenol-free to not hinder polymerization of resins and composites. Contains fluoride for a cariostatic action.

Safety notes

Warnings

- Elсотemp contains Zinc Sulfate which may seriously irritate eyes.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek for medical advice.
- Harmful if swallowed. Don't give anything to drink. Do not induce vomiting and seek for medical advice.

Precautions

- Wear suitable protective eyewear, mask, clothing and gloves.

Contraindication

Elсотemp is contraindicated in patients with a history of allergic reaction to Zinc Sulfate or any other of the components.

Instruction for use

Syringe form

- Remove cap protecting syringe end. Because of the thickness of the material, slowly turn the plunger of the syringe clockwise, to express the exact amount of cement required. Re-cap the syringe immediately after use.

Jar form

- Remove the desired amount of cement with a spatula. Tightly close the jar after use.

Application

- Apply the cement in the cavity, with a spatula, hardening start immediately with saliva. Remove it with an ultrasonic device or with a probe.

Storage & shelf-life

Store at a temperature between 12°C and 25°C. When stored under such conditions, Elсотemp has a shelf-life of 48 months. Store in a dry, cool and well-ventilated place, away from humidity.

Presentations & references

Three consistencies :
grey(soft) - white(normal) - pink(hard)
TPS-* – Box of 3 x 28g dispensing syringe
* : B (white) / R (pink)

TPP-* – Box of 6 x 28g jar
* : B (white) / G (grey) / R (pink)
TPP-1 – 1 x 28g jar
* : B (white) / G (grey) / R (pink)

Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions

and legal information are printed.

- Storage at a too high temperature could cause to age prematurely and soften the product.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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