
DETAX

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


REF 02611

INTENDED USE

Bite registration silicone for dental prosthetics

INDICATION:

- ▶ Bite registrations are used in the context of C&B technology as well as in orthodontics for registering and intraoral recording of the occlusion situation.
- ▶ greenbite apple and greenbite colour are suitable for applications that require a high degree of final hardness.
- ▶ Detaseal® bite, milkbite and glassbite are used for registration that requires a slightly lower level of Shore hardness.
- ▶ greenbite colour and milkbite facilitate fast bite registration as less time is required for intraoral setting.
- ▶ Detaseal® bite is used for bite registration that requires slightly longer processing time.
- ▶ milkbite is directly scannable and can be detected and registered visually without any contrast powder.

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist

PRODUCT DESCRIPTION

Detaseal® bite, greenbite apple, greenbite colour, milkbite and glassbite are addition-curing silicones that are categorized as Class 1 medical products. The silicones have a base-to-catalyst mixing ratio of 1:1 for precise bite registration.

MIXING AND DOSING

Extrude material from the cartridge with Automix2-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 3). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed homogeneously (Fig. 4). Only then dosing can be performed individually. Leave mixing cannula on cartridge after use, thus sealing the cartridge.

INSTRUCTION FOR USE

Extrude material from cartridge with mixing cannula or with attached contouring tip onto teeth (Fig. 5). At terminal occlusion application from vestibular is possible. Let registration harden in terminal occlusion and remove from the mouth.

DISINFECTION

Subsequent disinfection can be carried out with suitable commercially available disinfecting solutions according to the manufacturer's instructions.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained personnel.
- ▶ Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- ▶ Where movement is impeded or blocked, discard the cartridge. Do not apply force.
- ▶ Do not intermix with condensation curing silicones.
- ▶ Avoid staining of clothing. Hardened impression materials cannot be removed.
- ▶ Do not leave any residues of the impression material in the patient's mouth, and prevent inadvertent pressing of material into areas outside the oral cavity through mouth-antrum connections or malformations of the palate.
- ▶ Wear suitable personal protective equipment during the finishing of the cured material.
- ▶ Refer to the corresponding safety data sheet for hazard and safety information.

NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Read and understand the safety data sheet!

For users and / or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority of the Member State in which the user and/or patient is established.


CONTRAINDICATION


Do not use this material in case of known allergies to any of the ingredients, or contact allergies. Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

▶ See figures on page 75

Technical Data	Detaseal® bite	greenbite apple	greenbite colour	milkbite	glassbite
DIN 13903	x	x	x	x	
Mixed volume	50 ml (cartridge)				
Color of mixing cannula	 REF 02707	 REF 02705	 REF 02705	 REF 02707	 REF 02705
Colour code Base/Catalyst	yellowish green / white	green / yellow	blue / yellow*	creme / white	transparent
Working time**	approx. 45 sec.	approx. 30 sec.	approx. 30 sec.	approx. 30 sec.	approx. 30 sec.
Setting time**	approx. 2:15 min.	approx. 1:30 min.	approx. 1:20 min.	approx. 1:20 min.	approx. 1:30 min.
Time in mouth	approx. 1:30 min.	approx. 60 sec.	approx. 50 sec.	approx. 50 sec.	approx. 60 sec.
Final hardness	90 Shore A	45 Shore D	45 Shore D	90 Shore A	80 Shore A
Linear dimensional change	approx. 0.2 %				

Storage:

Application: At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity

Mixing time: omitted (Automix2-system)

*Note that the change in color (green to yellow) is an indicator of the setting status.

**from beginning of mixing

Increased temperatures accelerate, decreased temperatures retard a. m. times.

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GmbH & Co. KG