OptraMatrix

Instructions for Use Gebrauchsinformation Mode d'emploi Istruzioni d'uso Instrucciones de uso Instruções de Uso Bruksanvisning Brugsanvisning Καyttööhjeet Bruksanvisning Productinformatie Οδηγίες Χρήσεως

- Partially thinned circular matrix band
- Gezielt ausgedünntes Zirkulärmatrixband
- Matrice circulaire à zones de contact sélectives plus fines, et perforée sur la face opposée
- Matrice circolare a banda sottile
- Banda matriz circular parcialmente adelgazada
- Banda parcialmente afilada para matriz circular
- Rundformade matrixband med uttunnad sektion
- Cirkulært matricebånd med rude
- Osittain ohennettu hampaan ympäri vietävä matriisinauha
- Ekstra tynt sirkulært matrisebånd
- Gedeeltelijk verdunde ronde matrixband
- Κυκλικό τεχνητό τοίχωμα μερικώς μικρού πάχους

For dental use only.
Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Germany
Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

590875/0704/25/WE3/G

ivoclar
vivadent
clinical

English

Instructions for Use

Description

OptraMatrix has been specially designed to facilitate the contouring of direct multi-surface restorations in the posterior region (Class II). The selectively thinned contact areas of the matrix (thinned to 10 micrometer) and the window for the intact contact on the uninvolved side (in two-surface restorations) assist in creating precise, anatomically shaped contact points of consistent tightness. The window and the selectively thinned contact area ensure that only a thin layer of foil separates the adjacent teeth during the restorative procedure. The need for forceful separation of the teeth is thus eliminated.

Indication

Suitable for single use in direct Class II restorations in conjunction with composite filling materials.

Shapes

Two different shapes are available for two-surface restorations. The shapes are selected according to the position of the proximal box and the quadrant in question. The universal MOD shape is used for preparations that involve three surfaces.

On the reverse side of the matrix band, the following schemes can be found:

Legend

Quadrant I	Quadrant II
Quadrant IV	Quadrant III
Or,	
Teeth (1 – 8)	Teeth (9 – 16)
Teeth (25 – 32)	Teeth (17 – 24)
M – mesial	D – distal

1. OptraMatrix Molar MOD / Universal

three surfaces – all quadrants
(The gingival extensions may be trimmed using scissors, if required)

2. OptraMatrix Molar


Or,

I M	II D
IV D	III M
(1 – 8) M	(9 – 16) D
(25 – 32) D	(17 – 24) M

(The gingival extensions may be trimmed using scissors, if required)

3. OptraMatrix Molar


Or,

I D	II M
IV M	III D
(1 – 8) D	(9 – 16) M
(25 – 32) M	(17 – 24) D

(The gingival extensions may be trimmed using scissors, if required)

4. OptraMatrix Premolar MOD/Universal

three surfaces – all quadrants
(May also be suitable for small preparations in molars, or in pedo applications)

5. OptraMatrix Premolar


Or,

I M	II D
IV D	III M
(1 – 8) M	(9 – 16) D
(25 – 32) D	(17 – 24) M

(May also be suitable for small preparations in molars, or in pedo applications)

6. OptraMatrix Premolar

Or,

I D	II M
IV M	III D
(1 – 8) D	(9 – 16) M
(25 – 32) M	(17 – 24) D

(May also be suitable for small preparations in molars, or in pedo applications)

Application:

1. Select the appropriate matrix band according to the position of the preparation, as described above. Insert the matrix in a conventional Tofflemire® matrix retainer. Make sure that the inserted matrix is slightly tapered towards the root surface and the matrix retainer faces the mesial side.
2. Make sure that the thinned (matte finish) section of the matrix faces outwards towards the adjacent tooth. The thinned (matte finish) matrix surface must not face inwards towards the preparation (see picture). This is essential to achieving a smooth interproximal restoration surface.

Example of a matrix band that has been inserted correctly – the matte finished, thinned section of the matrix faces outwards.

3. For a two-surface Class II restoration, place the matrix in such a way that the window is seated around the intact contact on the uninvolved side (intact proximal area). The window enables the matrix band to be placed without separating the teeth in the contact area and thereby facilitates the creation of consistently tight contact points.

The window of the matrix fits around the intact contact point. As a result, the tooth is not forced out of position in the contact area (lower left first molar – restoration of a two-surface mesial cavity; buccal view).

4. Carefully tighten the matrix. Caution: Excessive tightening of the matrix may cause a tear in the thinned section.

5. It is advisable to use interdental wedges for a reliable seal of the gingival margin. The thinned metal band eliminates the need for forceful separation of the teeth.
6. Complete the procedure according to the restorative material manufacturer's instructions.
7. Remove the interdental wedges and undo the matrix retainer before removing the matrix. Be aware of the fact that there is a risk of gingival injury if the matrix is torn while removing it.

Notes

- OptraMatrix is suitable for single use only.
- If a flowable composite is used as a first layer to line the preparation and subsequently a contact point instrument (e.g. OptraContact from Ivoclar Vivadent) is employed, care should be taken to apply the flowable composite only to the preparation and not to the matrix band. If the flowable material comes in contact with the matrix band, it may adhere to the matrix band the band from being pressed against the adjacent tooth with the contact point instrument.

**Store out of the reach of children!
For use in dentistry only!**

**Date information prepared: 05/2004
Revision 0**

Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

Tofflemire is not a registered trademark of
Ivoclar Vivadent, Inc.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Dentigenix Inc.