Multilink®

Multilink®Primer A/B

Metal/Zirconia Primer

Instructions for Use

Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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- Self-curing resin based dental luting material
Self-curing and self-etching primer
- Selbsthärtendes zahnärztliches Befestigungscement
Selbsthärtender und selbstätzender Primer
- Composite autopolymérisant pour le collage des restaurations dentaires
Primer automordançant, autopolymérisable
- Cemento composito autoindurente
Primer autoindurente e automordenzante
- Material de fijación dental en base a resina autopolimerizable
Primer autopolimerizable y autograbante
- Composto de cimentação autopolimerizável
Primer autocondicionante e autopolimerizável

Brugsanvisning

Brugsanvisning

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Brugsanvisning

Productinformatie

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For dental use only.
Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.
Made in Liechtenstein
Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

Complies with / entspricht:
ISO 4049, ANSI/ADA Spec. No. 27

C E 0123

ivoclar
vivadent®
clinical

English

Instructions for Use

Description

Multilink is a self-curing and self-etching luting composite system for the adhesive luting of indirect restorations made of metal, metal-ceramic, all-ceramic and composite. Multilink has been developed for applications where light curing is not possible or difficult and for indications where high retention and bond strength are desired. Multilink is applied together with the Multilink Primer A/B. For this purpose, the two Multilink Primer components A and B are mixed in a 1:1 ratio directly before the application. Subsequently, the mixture is applied to the dentin and the enamel.

Multilink features a reliable and fast self-cure with excellent mechanical properties and a long-term hydrolysis-resistant matrix due to its pure composition of composites. Together with the self-curing primer, very good bond strength values and a good marginal adaptation is achieved. The Metal/Zirconia Primer is recommended as a coupling reagent to achieve a stronger bond to precious and non-precious alloys, as well as all-ceramics made of zirconium and aluminium oxide.

Shades

Multilink is available in three shades with different degrees of translucency:

- Transparent (high translucency)
- Yellow (high translucency)
- Opaque (low translucency)

Working and setting time

The working and setting times depend on the ambient temperature. The following times apply after the two pastes have been mixed (mixing time approx. 20 sec):

	Room temperature 23°C ± 1°C	Intraoral approx. 37°C ± 1°C
Working time (with- (out mixing time)	180 ± 30 sec	120 ± 30 sec
Curing time (without working time)	300 ± 30 sec	180 ± 30 sec

Mixing ratio/Dosage

The single pastes are dispensed from the Multilink double-push syringe and mixed in a 1:1 ratio. Multilink Primer A and Multilink Primer B are mixed in a 1:1 ratio (e.g. 1 drop of Primer A + 1 drop of Primer B) or a multiple of it.

Composition

Multilink

The monomer matrix is composed of ethoxylated Bis-EMA, UDMA, Bis-GMA, HEMA. The inorganic fillers are barium glass, ytterbium trifluoride, spheroid mixed oxide. The particle size is 0.25 – 3.0 µm. The mean particle size is 0.9 µm. The total volume of inorganic fillers is 39.7 %.

Multilink Primer A and B

Multilink Primer A contains an aqueous solution of initiators. Multilink Primer B contains HEMA and phosphonic acid and acrylic acid monomers.

Metal/Zirconia Primer

The Metal/Zirconia Primer contains phosphonic acid acrylate and methacrylate cross linking agents in an organic solution.

Indication

Multilink and Multilink Primer A/B are used for the permanent cementation of indirect restorations where a strong bond is desired:

- Metal and metal-ceramics (inlays, onlays, crowns, short-span bridges)
- Adhesive bridges (Maryland bridges)
- All-ceramics, especially opaque zirconium oxide ceramics, as well as composites (inlays, onlays, crowns, metal-free adhesive bridges)
- Root posts made of metal, ceramic and fibre-reinforced composites

Contraindication

The use of Multilink and Multilink Primer A/B is contraindicated

- If a dry working field cannot be established or the stipulated working procedures cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Multilink and Multilink Primer A/B.
- No direct application on pulpa-near dentin and on the open pulp or irritations of the pulp may occur.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to single components have been reported.


Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization.

Consequently, application of materials containing these substances is to be avoided.

Application

Overview of the individual steps for luting with Multilink and Multilink Primer A/B:

1. Removal of temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or the prepared die with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray. Subsequently, dry with water- and oil-free air.

Note:

Only use eugenol-free temporary cements, as eugenol may inhibit the polymerization of Multilink. If System®.inlay/System®.onlay is used for inlays as temporary filling material, a temporary luting cement is almost superfluous.

2. Try-in of the restoration and drying

Next, the shade, fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic materials before they are permanently cemented as there is a risk of fracture. In individual cases, do not check the occlusion at this stage. If necessary, make corrections with fine diamonds at medium speed and with slight pressure. Polish ground surfaces.

When applying the adhesive luting technique with composites, absolute drying of the operatory area - preferably with a rubber dam, alternatively with cotton rolls and saliva ejector - is indispensable.

3. Pre-treatment of the restoration

According to the instructions of the manufacturer. Otherwise:

3.1. Non-precious/precious alloys, all-ceramics made of zirconium and aluminium oxide

a) Sandblasting (micro-mechanical bond) + Metal/Zirconia Primer

- Sandblast the restoration surface (parameters according to the instructions of the manufacturer)
- Clean the restoration in an ultrasonic unit for about 1 minute.
- Rinse with water spray.
- Dry the restoration
- Apply the Metal/Zirconia Primer: Apply a thin layer of the primer with a micro-brush and leave to react for 180 seconds. Subsequently, dry with water- and oil-free air.

Alternative non-precious/precious alloys:

b) Silicitation (mechanical-chemical bond, especially recommended for high-gold alloys)

The bond strength can be optimized through silicitation and subsequent silanization with a usual bonding silicane (e.g. with Monobond-S). Pre-treatment according to the instructions of the manufacturer.

3.2. All-ceramics based on silicate (e.g. IPS Empress®)

- Etching with IPS® Ceramic etching gel (hydrofluoric acid) according to the instructions of the manufacturer
- Thoroughly rinse and dry the restoration
- Silanize the inner surfaces of the restoration e.g. with Monobond S: apply, let react for 60 sec and dry with water- and oil-free air.

3.3. Lab-fabricated composite restorations (the following notes apply for Tetric® Ceram)

- Sandblast the restoration surface
- Clean the restoration in an ultrasonic unit for about 1 minute.
- Rinse with water spray.
- Dry the restoration
- Silanize the inner surface of the restoration e.g. with Monobond S: apply, let react for 60 sec and dry with water- and oil-free air.

4. Mixing of Multilink Primer A and Multilink Primer B

Mix the two Multilink Primer liquids A and B in a 1:1 mixing ratio (e.g. 1 drop of Primer A and 1 drop of Primer B). The mixed Primer A/B is self-curing and does not need to be protected against light, but it has to be processed within 10 minutes.

5. Application of the mixed Multilink Primer A/B on the enamel and the dentin

Apply the mixed Multilink Primer A/B with a micro-brush on the whole tooth surface (cavity/die), starting from the enamel and brush with slight pressure for 15 sec. A reaction time of 30 sec is recommended on the enamel, and 15 sec on the dentin. The applied primer is subsequently dried with water- and oil-free air. As the primer is solely self-curing, no light curing is necessary!

6. Mixing of the Multilink cement and application on the restoration

Dispense Multilink from the double-push syringe and mix the two pastes in a 1:1 ratio on the mixing pad directly before the restoration is placed.

Note

- Multilink should be applied soon after mixing and the restoration seated into place. As soon as Multilink gets into contact with the mixed Multilink Primer A/B, the self-curing reaction is accelerated so that the curing time is shortened.

Apply the mixed Multilink on the inner surface of the restoration. The application of Multilink directly on the die or in the cavity, which has been pre-treated with Multilink Primer A/B, is not indicated as this would lead to a considerable acceleration of the curing process and thus may cause improper fit.

7. Placement of the restoration and removal of excessive cement

Seat the restoration in place and fix/hold. Remove excesses immediately with a micro-brush/brush/pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). Due to the reaction between Multilink and Multilink Primer A/B a high bond strength and degree of cure is achieved only a few minutes after the placement of the restoration.

Note

As any other composite, Multilink is subject to oxygen inhibition. This means that the surface layer (approx. 100 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. Cover restoration margins with glycerine gel/air block (e.g. Liquid Strip) after removal of excess and rinse off after complete polymerization.

8. Finishing of the restoration

- Apply finishing and polishing strips in the proximal regions.
- Check occlusion and functions and make corrections if necessary.
- Polish restoration margins with silicone polishers (Politip-F®, Astropol®) or discs.

Specific instructions for the cementation of endodontic posts:

1. For the cementation of endodontic posts, the root canal and occlusal surface of the preparation are ideally coated with the mixed Multilink Primer A/B using the thin (violet) micro-brushes, and left to react for about 15 seconds. Remove excess material from the root canal using paper tips.
2. Coat the selected and prepared (according to the instructions of the post manufacturer) endodontic post with the mixed Multilink cement.
Important: Do not lute the Multilink cement into the root canal that has been wetted with the Multilink Primer A/B. Otherwise, premature curing might occur which possibly prevents the post from being placed into its required position.
3. The endodontic post is placed and tightly held, thus producing excess cement.
4. The excess Multilink cement can be distributed across the treated preparation surfaces for total coverage.
5. The core build-up material (eg MultiCore®) can then be directly applied onto the Multilink which thus serves as bonding agent, and cured according to the instructions of the manufacturer (light-curing).

Warning

Unpolymerized Multilink pastes and the Multilink Primer A/B liquids may cause slight irritations. Avoid contact with the skin, mucous membrane or eyes. If Multilink comes into contact with the eyes, immediately rinse with water and, if required, seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water.

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Multilink and Multilink Primer A/B after expiry date.
- Cool storage (2–8 °C/36–46 °F) for Multilink and Multilink Primer A/B.
- Multilink double syringe and primer bottles should be closed immediately after use.
- Shelf life: see expiry date

Keep material out of the reach of children!

For use in dentistry only!

Date information prepared:

07/2005

Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.