parkell®

DIGITEST³

Pulp Vitality Tester

Includes:

- (1) Digitest 3 Pulp Vitality Tester
- (4) Autoclavable Stainless-Steel Tooth Probes
- (1) Lead Wire
- (1) Ground Clip
- (1) High-Output Alkaline 9-Volt Battery

DA..... 9	FR..... 33	NO..... 57
DE..... 17	IT..... 41	SL..... 65
ES..... 25	NL..... 49	SV..... 73

This manual, in whole or in part, should not be considered a substitute for formal training. Appropriate dental professional education is strongly recommended prior to using this device clinically.

Parkell Technical Support Service:

Monday – Friday, 8:30 AM to 5:00 PM EST
+1-800-243-7446 • techsupport@parkell.com

- A. Stimulus Adjustment / Control Button
 B. Digital Display
 C. Ground Clip
 D. Lead Wire

- E. Posterior Autoclavable Probe
 F. Anterior Autoclavable Probe
 G. Precision Labial Autoclavable Probe
 H. Precision Lingual Autoclavable Probe

Device Description

The Digitest® 3 Pulp Vitality Tester is a hand-held, battery-powered dental diagnostic device that identifies a living tooth nerve by stimulating it with a weak electric current. When the operator depresses the button, the strength of the electrical stimulus automatically increases at one of three preset rates. The unique waveform is designed to trigger a patient response in a vital nerve with a minimal amount of discomfort.

Intended Use/Indications

The Digitest 3 Pulp Vitality Tester is intended to be used as a diagnostic instrument to assist in the determination of the vitality of the dental pulp. It is indicated for use on vital and non-vital adult human teeth.

Contraindications

This Digitest 3 Pulp Vitality Tester is contraindicated for use on a patient or by an operator wearing a cardiac pacemaker or any other intracorporeal electronic device (internal defibrillator, insulin pump, etc.).

Warning

- Do not modify this device. Modification may violate safety codes, endanger the patient and the operator, and void the warranty.
- This device should only be used by licensed dental professionals qualified in the use of the unit.
- Read and understand all instruction manuals before using the device.
- Portable RF communications equipment (including peripherals such as antenna cables

and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Digitest 3 unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment may result.

- Use of this equipment adjacent to, or stacked with other equipment, should be avoided because it may result in improper operation.
- This device is to be operated with Parkell accessories only. Use of accessories other than those specified, or provided by Parkell Inc., may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment or improper operation.
- Equipment not suitable for use in the presence of flammable or explosive gases. Use of dental nitrous oxide/oxygen analgesia is acceptable.

Conformance to Standards

- Parkell's quality system is certified to ISO 13485, and this device conforms with IEC 60601-1, IEC 60601-1-2, CAN/CSA-C22.2 No. 60601-1 and IEC 60601-2-40.

Cleaning and Infection Control of the Digitest 3

- **DO NOT AUTOCLAVE THE DIGITEST 3 POWER UNIT, AS THIS WILL CAUSE DAMAGE TO IT.**
- Autoclaving and disinfecting do not remove accumulated debris. Before autoclaving or disinfecting accessories: Rinse the accessories under warm running water for 30 seconds to remove any external or internal soil or debris. Using a soft soapy cleaning brush to assist in the cleaning, if necessary. Use non-ammoniated detergent or dishwashing soap. Do not use ammoniated cleansers or

disinfectants. Rinse the item again under warm running water for 30 seconds to remove any residual soap and blot dry with a dry lint-free towel.

- The Power Unit should be protected using properly-sized disposable plastic barrier sleeves.
- Since the Digitest 3 incorporates sophisticated electronic circuitry, it should not be directly sprayed with or soaked in disinfectant. The device may be disinfected by wiping it with a lint-free towel moistened with an EPA-approved disinfectant containing not more than 30% alcohol, alcohol-free surface disinfectant wipes, or disinfecting wipes containing 3-5% Hydrogen Peroxide. Any disinfectant used on the lens and the housing should eventually be wiped off as per the instructions given on disinfectant labeling to remove chemical residue. If the lens becomes hazy or scratched,

NOVUS No. 2- Fine Scratch Remover can be used to remove the damage as per the manufacturer's instructions.

- Follow the surface disinfecting protocol specified by the disinfectant manufacturer.
- For further information, consult the Centers for Disease Control website at www.CDC.gov, referencing the most recent version of the "Guidelines for Infection Control in Dental Health-Care Settings".
- The Tooth Probes, Ground Clip and Lead Wire may be sterilized in any conventional steam autoclave (130°-134° C for 15 min for Gravity Displacement units, or 4 minutes for Pre-vacuum Cycle units, followed by a 15-minute minimum cool down interval). Comply with the autoclave manufacturer's recommended procedure, if different. These items can be autoclaved for 250 cycles before requiring replacement.

Specifications

Protection Against Electric Shock:	Type BF applied part
Protection Against Ingress of Liquids:	Pulp Tester - IPX0 (Ordinary)
Mode of Operation of Equipment:	Continuous
Operating Conditions:	15-40°C, 10-80% RH (non-condensing)
Transport and Storage Conditions:	10-40°C, 10-80% RH (non-condensing)
Stimulator Output:	The voltage output at the electrode can be characterized as discretely pulsed voltage packets that gradually increase in amplitude as the display count is incremented. The output voltage is purely AC and contains no DC component. See "Electrical Specifications" section

Electrical Specifications

Maximum DC Component Amplitude:	None
Maximum AC Component Amplitude:	500 V peak to peak
Peak Output Current Maximum:	250uA
Pulse Durations (Packet Width):	High (3 bars): 60.5ms Medium (2 bars): 100ms Low (1 bar): 124ms
Pulse Repetition Period:	High (3 bars): 140ms Medium (2 bars): 228ms Low (1 bar): 284ms
Average Voltage Increment:	7 V

The tooth enamel of a patient is simulated using a 2 Megaohm resistive load. Using this load as a reference, the absolute maximum voltage that can appear at the electrode is 500V at a maximum output current of 250uA.

Service and Maintenance

- Do not open the Digitest 3 Power Unit, except to change the battery. There are no user-serviceable parts inside. Internal repairs are to be made only by authorized Parkell personnel, by returning the unit to the service address at the end of these instructions.
- Avoid dropping the Power Unit or subjecting it to physical shock.
- Battery should be removed if unit is to be stored unused for more than 30 days.
- To prolong battery life, the device automatically turns off after 14 seconds of inactivity.
- When the battery is low, the letters “LO” will blink on the screen when the device is initially turned on. If the battery is low, replace it as explained in the section “Installing/Replacing the Battery”.
- If a Tooth Probe becomes loose in the mount of the Power Unit, it may be tightened by carefully inserting a flat screwdriver blade into the split metal base of the probe and gently expanding the sides apart with a twist of the tool (*Figure 1*).

Fig. 1

- Discard and replace any Tooth Probe if the metal or the insulation becomes damaged during use.

Installing / Replacing the Battery

- The Digitest 3 Pulp Vitality Tester comes with a high-output, 9-Volt alkaline battery.
- To install the battery, remove screw (*Figure 2*) and carefully lift off the plastic battery door by gently pulling it away from the Power Unit (*Figure 3*). Withdraw the old battery from the battery enclosure, and unsnap it from the wire battery connector (*Figure 4*). Discard in accordance with applicable laws.
- The fresh battery should be snapped onto the connector and placed in the Power Unit so the top end attached to the connector is placed into the opening first (*Figure 5*). Slight pressure on the bottom of the battery will help it to slide home (*Figure 6*). Replace the plastic battery door by engaging the tabs on the bottom first (*Figure 7*), and tighten the screw gently (*Figure 8*).
- The replacement battery may be a Duracell® MN1604, Panasonic® 6AM6 or equivalent. Although rechargeable 9-Volt batteries may be utilized, they will require more frequent charging because of the considerable voltage required to perform the pulp testing procedure.

Digitest 3 Setup Prior to Performing Pulp Vitality Testing

- To accurately diagnose the condition of the tooth, pulp testing must be performed on a patient who has not received local anaesthesia or gas analgesia. Make sure to use the Digitest 3 **BEFORE** you perform these procedures.
- Fully explain to the patient what you will be doing prior to the pulp test, so they are not surprised during the test. Using the words

Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6

Fig. 7

Fig. 8

“feeling” or “discomfort” are always preferable to the word “pain”.

- Standard infection control protocol must be followed by the operator during pulp testing, by wearing disposable rubber, vinyl or nitrile gloves.
- The tooth to be tested and the adjacent teeth must be clean and dry. Interproximal embrasures should be made free of any impacted food debris, saliva, water or plaque by flossing, scaling with an instrument and air-drying before testing.
- Electrically isolate the tooth being tested from neighboring teeth or metallic restorations by wrapping it with a clear mylar strip, such as those used in fabricating Class 3 composite restorations.
- A dry cotton roll should be placed in the buccal vestibule to isolate the tooth from the lip and cheek.
- A plastic mirror should be used to keep the tongue off the tooth during testing.
- There are two alternatives for completing the circuit using the autoclavable Ground Clip and lead wire.

1. Method 1: Over the Lip

- Plug one end of the wire into the socket on the bottom of the Power Unit (*Figure 9*) and insert the Ground Clip (*Figure 10*) into the socket at the other end of the wire.
- Place a small amount of toothpaste on the patient’s lip (to enhance electrical conductivity), then place the Ground Clip over the patient’s toothpaste-coated lip and confirm there is good contact between the lip and the clip.

2. Method 2: Finger Hold

- Connect the wire and Ground Clip to the Power Unit as per Method 1.
- Place a small amount of toothpaste on the Ground Clip and have the patient securely hold the metal clip between the thumb and forefinger of either hand.

- For both methods, advise the patient what to expect during the test and why the test is being performed. Tell the patient that when the vitality test is underway, they should indicate, either verbally or by carefully raising their hand (for test utilizing the “Finger Hold” of the Ground Clip, they should raise their free hand), at the first sign of an awareness of the stimulus in the tooth. As soon as they make you aware, you will stop the test by releasing pressure on the button and removing the probe from contact with the tooth surface.
- For the Finger Hold Method, the doctor has the option of giving the patient some control over the testing procedure. This is done by telling the patient they can make the sensation in the tooth stop if they let go of the Ground Clip. Letting go of the Ground Clip will not turn off the Digitest 3 and will not stop the numbers on the digital screen from continuing to count higher. It will only stop the patient from feeling the stimulus.
- **IMPORTANT: This option should be offered on a limited basis, and only to patients who may be hesitant to undergo the procedure due to concerns of feeling extreme discomfort. Patients offered the “Drop Clip” alternative may anticipate discomfort and drop the clip prematurely, thereby compromising the accuracy of the test results.**
- Insert Tooth Probe into the open end of the Power Unit (*Figure 11*). To enhance electrical conductivity and contact between tooth and the probe, apply a small amount of toothpaste to the metal tip.
- Metal or ceramic surfaces cannot be used as touch points for the Tooth Probes. Teeth selected for pulp testing must possess enough exposed enamel or dentin to allow the probe to make contact without touching the gingiva or a metal restoration. In certain cases, this may require the use of a Precision Probe (included).

Fig. 12


Fig. 13

Vitality Testing with the Digitest 3 Pulp Tester

1. To activate the unit, press and hold the start button (Figure 12) for a half second, and then release the button. When the button is depressed, the display will show one row of horizontal bars if the stimulus rate is set for "SLOW", two rows for "MEDIUM", or three rows for "FAST" (Figure 13).
2. If the desired stimulus rate mode was displayed, the vitality test may proceed, as detailed in Step 4.
3. If you wish to change the stimulus rate mode, simply press the button and the mode will change after the button is released. **NOTE: The unit will remember the last stimulus rate mode setting, even if the device is powered down.**
4. Place the toothpaste-covered tip of the Tooth Probe on the middle of the labial or lingual surface of the tooth. Avoid soft tissue and restorations such as crowns, amalgams or composites.
5. Depress and hold the button, and the display number will rise, indicating that a gentle stimulus is being applied to the tooth.
6. When the patient indicates that they feel the stimulus, pressure on the button should be released. Releasing pressure on the button will stop the test and freeze the digital screen. The screen will remain illuminated for approximately 14 seconds once pressure is released, which is enough time to remove the probe from contact with the tooth and record the number at which the patient felt sensation.
7. The maximum stimulus reading is 64. Even if there is no response at this level, there is still the possibility that the tooth is vital. No sensation at this number may simply indicate that the tooth is non-responsive at the time of the test, possibly from trauma. However, since teeth have been known to recover from traumatic injury many days after presenting with a "non-vital" reading, follow-up testing is almost always indicated after any initial readings. If this read-

ing persists over several visits, it is reasonable to assume that the tooth is non-vital. However, this conclusion should always be confirmed by another accepted endodontic testing method.

8. To confirm the diagnosis, a corresponding control tooth in the same arch should be tested. Molars should be matched to molars, premolars to premolars, cuspids to cuspids, and incisors to incisors. If this is not possible because teeth are absent, endodontically treated, or have full coverage restorations, a similar tooth in the opposite arch should be used.

- **Testing Tip: Do not combine the "Over the Lip" and "Finger Hold" methods for testing. For each patient, whichever method is chosen should remain consistent for the entire testing process on that patient. If the decision is made to change methods mid-test, the entire test will need to be repeated.**

Clinical Observations

- It is not possible to prepare a "table of normal values" for pulp tester readings, because THERE IS NO "NORMAL" IN PULP TESTING. Rather, the clinician should perform sequential comparisons between the subject tooth and the control tooth at consecutive office visits, observing how the readings are changing as time progresses. By utilizing electric pulp testing, along with all available diagnostic information, it is often possible to predict where the tooth's vitality is heading. This allows the clinician to make informed decisions as to whether endodontic therapy is appropriate, or whether it is prudent to simply watch and wait.
- There are general anatomic trends in pulp vitality readings. Posterior teeth generally require greater stimulus than anteriors, probably because of the greater thickness of enamel and dentin in posterior teeth. Enamel requires a greater stimulus than dentin or cementum, because of the higher percentage of non-conductive mineral, and the lower percentage of water. Similar cross-arch teeth will have similar thresholds to stimulus.

- The stimulus threshold may also be affected by such factors as the age, gender, previous pain history of the patient, pulp chamber size, trauma, pathology and use of prescription and non-prescription medications, or illicit drug use.
- The Parkell Digitest 3 is to be used in conjunction other diagnostic tests such as x-rays, temperature tests, percussion, etc. in order to confirm tooth pulp vitality.


Troubleshooting



Display does not light when button is depressed.
<ul style="list-style-type: none"> • Battery weak or dead – needs replacement. • Unit damaged – needs service.
<i>Tooth Probe loose or rotating in Power Unit.</i>
<ul style="list-style-type: none"> • Tooth Probe base too loose—adjust base with screwdriver (See <i>Figure 1</i>).
<i>Vitality test showing no response at “64” reading, while control tooth is normal.</i>
<ul style="list-style-type: none"> • Tooth is non-vital. • Tooth Probe not adequately contacting tooth – apply toothpaste to probe and reapply to tooth. • Ground Clip Lead Wire not connected to Power Unit. • Patient not holding “Ground Clip” tightly, or without toothpaste on it. • Consider applying small amount of toothpaste on “Ground Clip” to enhance conductivity before placing it on the lip or holding it in hand.
<i>Vitality test showing immediate response at very low reading.</i>
<ul style="list-style-type: none"> • Tooth is hyperemic – recommend endodontic treatment. • If problem persists with all teeth tested, unit needs service.

Warranty and Terms of Use

For full Warranty and Terms of Use information, please see www.parkell.com. Parkell’s Quality System is certified to ISO 13485. If you have any questions, please email our Technical Support Service at techsupport@parkell.com

Explanation of Symbols Used

	Professional use only
	Medical Device
	Catalogue / stock number
	Unique Device Identifier
	Serial number
	Follow instructions for use
	Do not use if package is damaged. This symbol is on packaging.
	Non-sterile
	Temperature limit
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	Keep dry

	Do not dispose this product into the ordinary municipal waste or garbage system
	Package contents
	Pulp Tester
	Autoclavable Probe(s)
	Lead Wire
	Ground Clip
	9V Battery
	Manufacturer
	Importer
	Translation

parkell®

This precision dental device was designed, manufactured and is serviced in the United States of America by:

 **Parkell, Inc.,**
300 Executive Dr., Edgewood, NY 11717 USA
+1-800-243-7446 • www.parkell.com

 **UKRP: Topdental (Products) Ltd,**
12 Ryefield Way, Silsden, West Yorkshire BD20, 0EF England UK

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