FR

COMPOSITION

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PROPRIÉTÉS

Biodentine™ is a bioactive dentine substitute from the «Active Biosilicate Technology™» innovation. Biodentine™ has mechanical properties similar to the active biosilicate technology™ in innovation. Biodentine™ is mainly contains high-purity, monomer-free mineral ingredients and is highly biocompatible. Biodentine™ creates the optimal conditions for maintenance of pulp vitality, by providing a very light seal on the dentine surface. It therefore guarantees the **absence of postoperative sensitivity** and the **longevity of restorations in vital teeth**. Biodentine™ also promotes reactionary dentine genesis and dentine bridges for unique pulp healing properties

INDICATIONS

- Permanent dentine restoration under composites or inlay/onlay
- Temporary dentine-enamel restoration.
- Restoration of deep and/or large coronal carious cavities (sandwich technique).
- The repair of root canal radicular lesions.
- Pulp capping.
- Root-end filling in endodontic surgery (retrograde filling).

IN THE TOOL

- Repair of root perforations.
- Repair of furcation perforations.
- Repair of perforating internal resorptions.
- Repair of external resorption
- Root-end filling in endodontic surgery (retrograde filling).

CONTRAINDICATION

Allergy to one of the ingredients.

LIMITS OF USE

- Restoration of a large loss of tooth substance subjected to high stresses.
- Esthetic restoration of anterior teeth.
- Treatment of teeth with irreversible pulpitis.

SIDE EFFECTS

- No known side effects.

INSTRUCTIONS FOR USE (FOR EACH INDICATION)

- Take a capsule and gently tap it on a hard surface to loosen the powder.
- Open a capsule and place the powder on the white capillary holder.
- Detach a single-dose container of liquid and gently tap on the sealed cap to force all the liquid down the container.

1) Twist cap to open. Be careful that no drop of liquid falls out of the single dose container.

- Pour 5 drops from the single-dose container into the capsule.
- Close the capsule. Place the capsule on a mixing device, such as Technomic, Tac 400 (Lineatoc), Siliamat, CapMix, Rotomix, Ultramat etc., at a speed of 4000 – 4200 rotations/min.
- Mix for 30 seconds.

- Open the capsule and check the material's consistency. If a thicker consistency is preferred, wait for 30 sec to 1 min before checking again. Do not exceed the working time.
- Collect Biodentine™ with the instrument supplied in the box. Depending on the desired application, you may handle Biodentine™ with an amalgam carrier, a spatula or a Root Canal Messing Gun. Rapidly rinse and clean the instruments to remove any residual material.

RESTORATION

Assess pulp vitality by the usual tests: Biodentine™ is not indicated for the treatment of teeth with irreversible pulpitis

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
- Adapt a matrix around the tooth if a wall is missing.
- Insert Biodentine™ in the cavity, avoiding trapped air bubbles. Condense the material and ensure good adaptation to the cavity walls and margins.
- Model the surface of the restoration.
- Wait until the end of the setting time of the material before removing the matrix.
- Optimize the mechanical properties of the material and facilitate removal of the matrix, a varnish can be applied onto the surface of the restoration.
- Check occlusion.
- Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material (composite, inlay/onlay, amalgam). The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restoration techniques, and particularly with all types of bonding systems.

PULP CAPPING

Assess pulp vitality by the usual tests: Biodentine™ is not indicated for the treatment of teeth with irreversible pulpitis

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
- Adapt a matrix around the tooth if a wall is missing.
- If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™.
- Prepare Biodentine™ as indicated above (*Biodentine™ mixing instructions*).
- Place Biodentine™ directly on the exposed pulp without condensation. Ensure good adaptation to the cavity walls and margins.
- Check occlusion.
- Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material (composite, inlay/onlay, amalgam). The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restorations, and particularly with all types of bonding systems.

1) Isolate the tooth with a rubber dam.

2) Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.

3) Adapt a matrix around the tooth if a wall is missing.

4) If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™.

5) Prepare Biodentine™ as indicated above (*Biodentine™ mixing instructions*).

6) Place Biodentine™ directly on the exposed pulp without condensation. Ensure good adaptation to the cavity walls and margins.

7) Check occlusion.

8) Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material (composite, inlay/onlay, amalgam).

9) The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restorations, and particularly with all types of bonding systems.

10) Complete root canal treatment at the next visit according to current recommendations.

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