(DE)	Gebrauchsanweisung	3
(GB)	Instructions for use	11
(FR)	Mode d'emploi	18
(ES)	Instrucciones de uso.....	26
(IT)	Istruzioni per l'uso	34
(PT)	Instruções de uso	42
(NL)	Gebruiksaanwijzing	50
(SE)	Bruksanvisning	58
(DK)	Brugervejledning	65
(NO)	Bruksanvisning	72
(FI)	Käyttöohjeet	79

MD

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Medical device

Dispositif médical

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Medisch hulpmiddel

Medicinteknisk produkt

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i|BOND[®] Etch 35 Gebrauchsanweisung (DE)

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Stand: 2021-12

Medical Device – for use by dental health care professionals only.

Intended purpose

iBond Etch Gel and Fluid are etching agents based on orthophosphoric acid. They are used to etch enamel and condition dentine in the context of adhesive filling therapy.

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Composition

iBond Etch 35 Gel contains: water, phosphoric acid, silicon dioxide.

iBond Etch 35 Fluid contains: water, phosphoric acid, polyvinyl alcohol.

Indications

Enamel etching and enamel/dentine etching in the context of the total etch technique prior to the adhesive bonding of:

- Composite and compomer restorations
- Laboratory-fabricated, indirect restorations (e.g. inlays, crowns, bridges, veneers)
- Sealants

Contraindications

Application to dentine close to the pulp (less than 1 mm). In case of known or suspected allergy to components of the product, its use is contraindicated.

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

Product application

Prior to treatment, carefully clean tooth with fluoride-free polishing paste. The use of a rubber dam is recommended. Visually inspect for damage before use. Damaged products must not be used.

1. After preparation in accordance with adhesive filling therapy, rinse the cavity and dry with oil-free air.
2. Starting from the bevelled enamel areas, apply iBond Etch to the entire cavity surface, including the dentine. Ensure that the enamel is conditioned for at least 20s, a maximum of 30s and the dentine for a maximum of 15s. If desired, only apply iBond Etch to the enamel and let set for a maximum of 30 s. Etch unprepared enamel, e.g. in case of fissure sealing, for at least 30 s and no more than 60 s. In case of fluoridated enamel, etching for no more than 60 s may also be required.

Caution: In case of deep cavities, we recommend covering the dentine near the pulp with a calcium hydroxide preparation and a base (e.g. glass ionomer cement) prior to etching.

3. Use a water sprayer to completely rinse off iBond Etch for 20 s under suction. The cavity must be completely free of any remnants.
4. Dry the cavity with oil-free air (do not dry out dentine!). The etched enamel surface should appear chalky white. In case of contamination of the tooth surface (e.g. with saliva), the etching must be repeated (for a maximum of 10 s dwell time).
5. Perform further treatment with adhesive (e.g. iBond Total Etch) and composite (e.g. VENUS®, CHARISMA®) in accordance with manufacturer directions.

iBond Etch Gel may be applied using a brush or cannula, and iBond Etch Fluid may be applied using a brush. Both brushes and cannulas are for one-time use only.

Release the plunger after applying the gel. When using syringes in mouth, use common hygienic protective sleeves to prevent contamination. Do not re-use after contamination. After use, close iBond Etch Gel and iBond Etch Fluid with the appropriate plug or cap.

Interactions

Liner and base materials, such as calcium hydroxide, may be dissolved by the acid, depending on their type and composition.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Warnings and safety instructions

iBond Etch 35 Gel and Fluid 35 contain 35% phosphoric acid and can cause chemical burns. May cause skin and eye irritation. Avoid contact with skin. Skin contact may cause sensitization. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. Eye protection is recommended

for the patient. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist.

The etching agent must not be swallowed by the patient, otherwise incompatibilities may occur.

Special advice

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com

To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Storage conditions

Store below 25°C (77°F). Protect from sunlight. Store in a dry place. Do not use after expiry date.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number in all correspondence about the product. Article number **REF** and batch number **LOT**: See note on bottle/packaging.

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Kulzer GmbH
Kulzer, LLC