
Document Detail

Type: 1450-DFU
Document No.: 1450-GUTTA CONDENSOR_DFU[8]
Title: GUTTA CONDENSOR_0514_DFU
Comment F19 03 08.X
Status: CURRENT
Effective Date: 30-May-2014

GUTTA - CONDENSOR

Ref A0242 - A 0244

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
Visit our website: www.dentsplymallefer.com

FOR DENTAL USE ONLY

GB DIRECTIONS FOR USE

- INDICATIONS FOR USE**
These instruments are to be used only in a clinical or hospital environment, by qualified users.
Application field: The Gutta-Condenser is used to mechanically condense gutta percha in the root canal, after instrumentation, cleaning and disinfection are finished.
- CONTRAINDICATIONS**
 - The Gutta-Condenser must not be used in curved canals, as in this clinical situation the risk of separation of the instrument is high.
 - The Gutta-Condenser is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause: difficulty breathing, swollen face or eyes, hives or rash. Should any of these symptoms occur, the patient should be advised to contact his/her dental professional immediately.
- WARNINGS**
 - The use of the Gutta-Condenser is not risk-free. The main danger consists in a possible penetration of the gutta percha into the periapical region. The following reasons may be responsible for such a case: the use of a too-small or too-fine master cone, a too-small Gutta-Condenser, excessive speed or axial pressure, or too-deep insertion of the Gutta-Condenser in the root canal. Even a small amount of gutta percha in the periapical region may cause severe clinical reactions. There is no resorption to be expected from the periapical macrophages.
 - For multicuspид mandibular teeth, an overfilling reaching the inferior dental nerve may provoke irreversible damage, such as dyesthesia, paraesthesia, anaesthesia, or paralysis with a possible labiodynism.
- PRECAUTIONS**
 - Never use the Gutta-Condenser in counter-clockwise rotation (danger of getting the stuck in the canal or in the bone).
 - Do not penetrate the tip of the instrument into the canal to a distance of less than 2 mm from the apex.
 - Do not oppose strong resistance when the Gutta-Condenser is backing out of the canal
- ADVERSE REACTIONS**
 - For multicuspид mandibular teeth, an overfilling reaching the inferior dental nerve may provoke irreversible damage, such as dyesthesia, paraesthesia, anaesthesia, or paralysis with a possible labiodynism.
 - The Gutta-Condenser is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause: difficulty breathing, swollen face or eyes, hives or rash. Should any of these symptoms occur, the patient should be advised to contact his/her dental professional immediately.
- STEP-BY-STEP INSTRUCTIONS**
See scheme above
The Gutta-Condenser is a stainless steel instrument. The rotation of this instrument against the Gutta-Percha plastifies this material by heat produced by friction. In this plastic state the Gutta-Percha is advanced apically and condensed.
In order to obtain sufficient heat to plastify the Gutta-Percha, the Gutta-Condenser must be used at a low speed with a high torque instrument capable of rotating at a low speed 8'000 min⁻¹ with a high torque. For a complete filling of the root canal system, bring the instrument to a distance of 2 mm from the apex of the canal which has already been prepared. The penetration depth of the Gutta-Condenser is checked by using the graduations marked on the shaft of the instrument or by using a silicone-stop.
- Establish the working length and apical stop of the root canal.
2. Insert the gutta percha point and the Gutta-Condenser into the canal, and press the Gutta-Condenser against the gutta percha.
3. Descend the Gutta-Condenser into the canal, rotating in a clock-wise direction at 8'000 min⁻¹. Do not penetrate the Gutta-Condenser to a distance of less than 2 mm from the apex.
4. Back the Gutta-Condenser out of the canal, without opposing strong resistance.
- DISINFECTION, CLEANING AND STERILIZATION**
- FOREWORD**
For hygiene and sanitary safety purposes, all instruments must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.
Instruments which are marked as "sterile" do not require any specific treatment before the first use, but have to follow this procedure for all subsequent use if not labelled as "single use".
- AREA OF APPLICATION**
Disinfection and sterilization before each use (except for the first use of sterile instruments) and reprocessing procedures concerning:

DISINFECTION AND STERILIZATION

- A. Device**
A1. Instruments:
Cutting instruments, (hand and engine driven) such as:
• Endodontic instruments (files, broaches, reamers, enlargers, endodontic burs, ultrasonic inserts);
• Rotary cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs).
Root canal filling instruments (Pluggers, spreaders, compactors);
Hand instruments, clamps and Rubber Dam accessories.
- A2. Accessories:**
Supports, kits, instrument organisers and other accessories
- DISINFECTION ONLY**
B. Filling material and calculable plastic posts: Only chemical disinfection (no sterilisation)
Gutta percha, Obturators, Unidip and Mosser Calcinate plastic posts.
- PRESENT PROCEDURE NOT APPLICABLE**
- C. Exclusion**
• Equipment such as Motors, Apex locators and other devices with reprocessing procedures included in the individual Direction for Use.
• MTA, Glyde, TopSeal, Paper Points, Rubber Dam.

- III - GENERAL RECOMMENDATION**
1 - Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFI of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use antirust disinfecting and cleaning agents.
2 - For your own safety, please wear personal protective equipment (gloves, glasses, mask).
3 - The user is responsible for the sterilization or disinfection of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterilization.
4 - It is safest for the practitioner to use our instruments only once. Should our instruments be reused, we recommend to always carefully inspect them before use: the appearance of defects such as cracks, deformations (bent, unwound), corrosion, loss of color coding or marking, are indications that the devices are not able to fulfil the intended use with the required safety level and must therefore be discarded. In any case, we recommend not to exceed the following maximum number of uses for our root canal shaping instruments:

| Type of canal | Stainless Steel instruments with a diameter ≤ ISO 015 | Stainless Steel instruments with a diameter = ISO 015 | NiTi instruments |
|--|---|---|------------------|
| extremely curved (>30°) or S-shaped canals | 1 canal max | 2 canals max | 2 canals max |
| moderately curved canals (10° to 30°) | 1 canal max | 4 canals max | 4 canals max |
| slightly curved (<10°) or straight canals | 1 canal max | 8 canals max | 8 canals max |

- Single use marked devices are not approved for re-use.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Tungsten carbide burs, plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H₂O₂) solution.
- Only the active part of the NiTi instruments should be immersed in a NaOCl solution at NOT more than 5%.
- Do not use acid (pH < 6) or alkaline (pH > 8) solutions with aluminium devices. These types of devices are degraded in presence of caustic soda solutions with mercury salt.
- The washer-disinfector is not recommended for devices made of aluminium, tungsten carbide or carbon steel.

IV - STEP-BY-STEP PROCEDURE

| Operation | Operating mode | Warning | Following uses | |
|--|--|---|----------------|-----|
| | | | First use | |
| 1. Disassembling | - Disassemble the device, if required | - Silicone stops have to be removed. | | X |
| 2. Pre-Disinfection | - Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible. | - Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor. - Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments (See general recommendations). - For visible impurities observed on instruments a pre-cleaning is recommended by brushing them manually with soft material. | | X |
| 3. Rinsing | - Abundant rinsing (at least 1 min) | - Use quality water in accordance with local regulations. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the cleaning. | | X |
| 4a. Automated Cleaning with washer-disinfector | - Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (No value > 3000 or, at least 5 min at 90 °C). | - Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer disinfecter use kits, supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly. | X | X |
| OR | | | | |
| 4b. Manual Cleaning and assisted by an ultrasonic device | - Place the devices in a kit, support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable. | - No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor. | X | X |
| 5. Rinsing | - Abundant rinsing (at least 1 min) | - Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air. | | X X |
| 6. Inspection | - Inspect devices and sort out those with defects. - Assemble the devices (stops) | - Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any defect as described in the General Recommendation above. - Protect carbon steel bur with corrosion inhibitor before packaging. - For Contra Angle : lubricate the device with an adequate spray before packaging | X | X X |
| 7. Packaging | - Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in "Sterilisation pouches". | - Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607. | X | X X |

| | | | | | |
|------------------|---|---|---|---|---|
| 8. Sterilization | - Steam sterilisation at: 134 °C / 273°F during 18 min | - The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use only autoclaves that are matching the requirements of EN 13060, EN 285. - Use a validated sterilisation procedure according ISO 17865 - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure - Control the efficiency (packaging integrity, no humidity, colour change of sterilisation indicators, physico-chemical integrators, digital records of cycles parameters). - Traceability of procedure records | X | X | X |
| 9. Storage | - Keep devices in sterilization packaging in a dry and clean environment. | - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period). | X | X | X |

B. Filling material and calculable plastic posts

| 1. | Operation | Operating mode | Warning |
|----|--------------|--|---|
| | Disinfection | - Immerse the obturation devices in NaOCl (2.5 % at least) during 5 min. at ambient temperature. | - Do not use disinfecting solutions containing Phenol or any products which are not compatible with the treated filling material (See general recommendations). |

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