DE ZINKPHOSPHATZEMENT

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EN ZINC PHOSPHATE CEMENT

INDICATIONS

Cementation in inlays, onlays, crowns, bridges and orthodontic bands.

If the cavity is deep enough to compromise the pulp is recommended to use a layer of calcium hydroxide before the use of this cement.

COMPOSITION

Zinc oxide and Aluminum phosphate hydrated


INSTRUCTIONS FOR USE

- (A) Use a cool mixing slab (72°F, 22°C) and a non-corrosive spatula, for mixing ZINC PHOSPHATE CEMENT.
- (B) Use a ratio of 1.3 to 1.4 gr. of POWDER to 0.85 g. of LIQUID.
- (C) Divide POWDER into four (4) equal portions. Divide one (1) of these quarter portions into two (2) parts.
- (D) ① Spatulate 1/8 of the portions into the entire amount of the liquid for 10 seconds.
② Add the other 1/8 of the portions and mix for 10 seconds.
③ ④ Add consecutively two (2) of the quarter portions, spatulating 10 seconds after each addition.
⑤ Add the final quarter portion and mix the entire amount for 20 seconds. This will give 1 minute total mixing time.

MIXING TECHNIQUE

SETTING - TIME CONTROL

- (A) Powder-Liquid ratio. Thick mixes set faster than thin mixes.
- (B) The faster the powder is added to the liquid, the quicker the set.
- (C) Prolonged mixing tends to lengthen setting time.

IMPORTANT

To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

To eliminate condensation of water vapor on the slab, raise the temperature of the slab slightly, if 72°F (22°C) temperature of slab is below the dewpoint.

Be certain all the powder is thoroughly incorporated into the liquid in each successive step of the mixing.

WARNING

- In case of any change on the original function of the product, stop using it and contact the manufacturer.

For Professional Dental use only.

WARNINGS

Intention of use Cementation in inlays, onlays, crowns, bridges and orthodontic bands.

Composition: Zinc oxide and phosphoric acid

The zinc phosphate cement is not irritable to the skin. In case of contact with the eyes, to wash immediately with abundant water during at least 15 minutes and if the inconveniences persist consult with a doctor. In case of contact with skin wash the area with soap and water. In case of inhalation clean the nose with water and medical assistance is required. Before the use of these products always use rubber dam.

Warnings: Some patients can present allergy to the zinc phosphate or any of its components. In this case suspend its use. Do not storage in temperatures above 24°C. Close the packages immediately after to dispense the product.

Only for use professional dental.

Contraindication: The material must not be used when the techniques application described are not possible.

Warranty proclinic warranties the product against any fault of manufacture or of the material, proclinic does not take responsibility for the use of the product out of the intention of use specified or out of the instructions of use, or out of the precautions, the warnings, and indicated contraindications. Being the only obligation is the replacement of the defective product.

Information Post-production:

If the use of these product on a patient or any one else some anomaly occurs please send detailed report to MEDENTAL INTERNATIONAL. medental_inc@hotmail.es

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