





Soft relining material

HAND MIXING SYSTEM

Instructions for use

- DE** Gebrauchsanweisung.
- EN** Instructions for use.
- ES** Instrucciones de uso.
- FR** Mode d'emploi.
- IT** Istruzioni per l'uso.
- PL** Instrukcja postępowania.
- PT** Instruções de uso.

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
LONG-TERM SOFT RELINING MATERIAL FOR CHAIRSIDE USE, A-SILICONE-BASED, COLD CURING

1. PREPARATION OF THE DENTURE

- Removal of existing relinings (fig. 1): The old, existing relining has to be removed completely. Grind off the denture base approx. 1–2 mm, depending on the situation. Incorporate in the areas of transition a shoulder of approx. 2–3 mm in height with a depth of minimum 1 mm. Round off the borders and clean the denture thoroughly.
- Degreasing (fig. 2): Degrease by means of pure medical alcohol and let dry.
- Applying the adhesive for adhesion (fig. 3): Apply the adhesive uniformly and completely onto the grinded denture areas to which the material has to adhere. Avoid contact with eyes! Wear suitable eye protection! Let the adhesive dry for approx. 1 minute. A continuous film of adhesive is essential for satisfactory reliner adhesion. After use close the bottle immediately. Saliva, humidity or monomer liquids prevent the adhesion of Proclinc Expert Soft Relining Material to the denture base.

2. DIRECT CHAIRSIDE RELINING METHOD

- Coating the denture: Prepare denture as described under point 1. Mix homogeneously equal lengths of Proclinc Expert Soft Relining Material base material and catalyst for 30 seconds. Apply bubble-free material with spatula onto the vestibular, lingual and dorsal borders. Then coat completely the remaining denture with a layer of minimum 2 mm. All denture areas in contact with the gingiva have to be coated.
- Functional movements (fig. 4): Seat the coated denture into the patient's mouth and let him close in normal occlusion. Setting is finished after approx. 5 minutes. Then remove the relined denture and trim and polish as described under point 4.

3. INDIRECT METHOD IN THE LABORATORY

Invest denture to be relined into a flask with

counter. For easier removal of the denture, prior to pouring the counter coat teeth with silicone. Remove denture from the flask and prepare as described under point 1. Coat the model with an usual alginate based insulation. Mix homogeneously equal lengths of Proclinc Expert Soft Relining Material base material and catalyst for 30 seconds. Apply bubble-free material with spatula uniformly onto the denture area to be relined (the layer of Proclinc Expert Soft Relining Material should be approx. 1 mm more than the base material has been grinded off). Close flask with counter. Press with light pressure for approx. 30 minutes or place flask in a pressure pot at 40–45 °C/104–113 °F for approx. 10 minutes in the water bath. Then deflask the relined denture and trim and polish as described under point 4.

4. FINAL TRIMMING AND VARNISHING

- Smoothing of uneven spots: Remove excess material with a sharp scalpel or fine scissors (fig. 5). Then smooth any rough areas or ridges with suitable grinding instruments (fig. 6).
- Polish (fig. 7): Gloss varnish is applied only for smoothing of the trimmed functional borders. Avoid contact with eyes and skin! Wear suitable eye protection and gloves! Open the aluminium cap with a spatula, squeeze both liquids onto a mixing pad in a mixing ratio of 1:1, spatulate to a homogeneous mixture and brush on the mixed liquid. Let dry at room temperature for approx. 5 minutes.

5. RECOMMENDATIONS FOR DAILY CARE OF RELINED DENTURE

Patients should be advised to keep dentures clean by brushing several times daily with a soft-bristled toothbrush and neutral soap and rinsing thoroughly. Dentures should not be soaked in cleaning baths for more than 10 minutes.

IMPORTANT WORKING HINTS

- Do not interchange tube closures.

- Absolutely firm bonding is obtained only by applying adhesive to roughened surfaces.
- Adhesive must cover entire area of contact completely.
- For proper denture adhesion, all surfaces of the denture in contact with gingiva must be coated evenly with Proclinic Expert Soft Relining Material.
- Setting times must be strictly observed.
- Please follow the instructions of the safety data sheets.
- Adhesive: avoid contact with eyes, mucous membranes, and skin. Do not inhale vapors.
- Gloss varnish: avoid contact of catalyst with eyes and skin, this can lead to irritation. In case of accidental contact, rinse eyes immediately with running water for at least 15 minutes, consult an ophthalmologist.

INDICATIONS FOR USE

- Suited for all polymethylmethacrylate based synthetics
- Fast and easy chairside preparation of permanently soft relinings; chairside or indirect use in the dental laboratory
- Correction of pressure areas in the upper or lower jaw
- Post-damming of the A-line in case of adhesion problems or with flabby ridges
- To steady the position of upper and lower dentures by improving adhesion
- Elastic covering of sharp ridged alveolar processes
- For covering implant regions and as overdentures

TECHNICAL DATA

DIN EN ISO 10139 – 2 Type A

- Mixed volume:
2 x 50 ml (tubes)
Mixing ratio: 1:1
- Colour code:
Base: dark red
Catalyst: clear blue
- Mixing time: manually approx. 30 sec.
- Working time: approx. 1 min. 30 sec.*
- Setting time: approx. 7 min.*
- Time in the mouth: approx. 5 min.
- Recovery from deformation: > 99.5%
- Linear dimensional change: < 0.2%
- Application: at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- Storage:

* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Further information

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

For use by trained specialist staff.

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DETAX GmbH & Co. KG