Prior to use, carefully read the instructions for use.

EN

GC Fuji ORTHO™ LC CAPSULE

LIGHT-CURED ORTHODONTIC BONDING ADHESIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- 1. Bonding metal brackets and attachments.
- 2. Bonding ceramic brackets.
- 3. Bonding acrylic appliances.
- 4. Band cementation only when extra band retention is desired.

CONTRAINDICATIONS

- 1. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- 2. Bonding polycarbonate brackets.
- 3. Cementation of stainless steel crown-retained appliances (e.g., Herbst Appliances).

DIRECTIONS FOR USE

- 1. ENAMEL PREPARATION
 - a) Using a prophy cup or a brush, clean the bonding surfaces of the teeth with plain (non-fluoridated) pumice and water.
 - b) Rinse thoroughly with water.
 - c) Using a sponge or a cotton pellet, apply GC ORTHO CONDITIONER to the bonding surfaces of the teeth for 20 seconds. Rinse thoroughly. Although successful results have been obtained with no enamel pretreatment, conditioning of the enamel bonding surfaces will increase the adhesive's bond strength. It is therefore recommended, especially for new users of GC Fuji ORTHO LC.


Etch the enamel bonding surfaces according to the etchant manufacturer's instructions. Rinse thoroughly.

The enamel bonding surfaces must be moist. An overly dry or desiccated enamel surface will adversely affect the bond strength. The optimum level of surface moisture can be obtained by wiping the bonding surfaces of the teeth with a moistened cotton roll immediately after the bracket bonding. If the bonding surfaces become desiccated during the bonding procedure, rehydrate by wiping them with the moistened cotton roll before bonding with GC Fuji ORTHO LC.

- 2. CAPSULE ACTIVATION AND MIXING
 - a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1). Hold the capsule and push the plunger until the end of the plunger is flush with the main body of the capsule (Fig. 2).

- b) Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 3). The capsule is now activated.
- c) Promptly place the capsule into a capsule mixer (or amalgamator). Mix for 10 seconds at high speed (Fig. 3). The working time is approximately 3 minutes from start of mix at 23°C (73°F). Higher temperatures will shorten working time and lower temperatures will extend it. Refrigeration of capsules is therefore recommended until immediately prior to activation and use.

- 3. EXTRUDING THE ADHESIVE
 - a) Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
 - b) Extrude the material by pulling the lever 4 or 5 times (Fig. 5). Coat the bonding surface of the bracket or appliance completely with adhesive.

- 4. BONDING PROCEDURE
 - A. Brackets
 - a) Coat the bonding surface of the bracket completely with the mixed adhesive.
 - b) Position the coated bracket on the tooth.
 - c) Press the bracket firmly against the enamel surface. Using an explorer or a scaler, remove the excess adhesive at this time. Press the bracket again. Continue placing additional brackets. Place all brackets in a quadrant or in the full arch. More than one mix of adhesive may be necessary.

- 5. DEBONDING PROCEDURE
 - A. Brackets
 - a) Cut the bracket with a ligature cutter and, while supporting the tooth with fingers from the lingual side, give a sharp clockwise rotation to the bracket. If resistance is met, try again twisting in the opposite direction after desiccating (using an air syringe) the tooth around the bracket. Remove any residual adhesive with a scaler or rotary instrument.



- B. Appliances
 - a) To remove bonded appliances, flex the appliance in multiple sites with forceps. If this does not remove the appliance, section the appliance, desiccate the adhesive, and remove the appliance.

- STORAGE
 - Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F).
 - (Shelf life : 2 years from date of manufacture).

- PACKAGE
 - Package : 50 Capsules.
 - Average contents per capsule : 0.30g powder and 0.10g (0.085mL) liquid. Minimum net volume of mixed adhesive per capsule : 0.1mL.
 - Option : GC CAPSULE APPLIER (1 piece).

- CAUTION
 - 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
 - 2. In case of contact with eyes, flush immediately with water and seek medical attention.

Last revised : 06/2015

MANUFACTURED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC. 3737 West 127th Street, A1sip, IL 60603 U.S.A. DISTRIBUTED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V. Researchpark Hasbroedre-Luven 1240, Interzonevleaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC. 3737 West 127th Street, A1sip, IL 60603 U.S.A. TEL: +1 708 597 0900

GC ASIA DENTAL PTE. LTD. 11 Tampines Concourse, #03-05 Singapore 528729 TEL: +65 6546 7388

GC AUSTRALASIA DENTAL PTY. LTD. 1753 Botany Rd, Bankmeadow, NSW 2019, Australia TEL: +61 2 9316 4499

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