D. Devices

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
Description	Professional Ultrasonic Cleaner	
Model	D_Ultrasonic Cleaning Bath 6L	
Tank Capacity	6.0l / 1.6 gallon (US)	Max. 5.0 L / 1.3 gallon (US)
		Min. 3.35 L / 0.9 gallon (US)
Tank Size	32.3x20.5x9.9cm/12.7"x8.1"x3.9"	
Longest Item Fits inside Tank	34.0 cm / 13.4"	
Power Supply	260W (AC100~120V 50/60Hz)	
	300 W (AC 220~240V 50/60Hz)	
	250 W (AC 100V 50/60Hz)	
Digital Timer Settings	<input type="text" value="10 min"/> → <input type="text" value="15 min"/> → <input type="text" value="20 min"/> → <input type="text" value="25 min"/> → <input type="text" value="30 min"/> → <input type="text" value="5 min"/>	
Drainage	Drainage valve	
Ultrasonic Frequency	35,000 Hz	
Tank Material	Stainless Steel SUS304	
Housing Material	ABS	
Net Weight	7.5 kg / 16.5 lb	
Gross Weight	8.6 kg / 19.0 lb	
Unit Size	51.5x30.5x24.0cm/20.3"x12.0"x9.4"	
Inner Carton Size	57.5 x 37.0 x 32.5 cm / 22.6" x 14.6" x 12.8"	
Qty per Master Carton	1 pcs /ctn	
Master Carton Size	60.0 x 38.5 x 34.5 cm / 23.6" x 15.2" x 13.6"	

Introduction

Use tap water. Special solutions are not necessary in most cases

Principles of ultrasonic cleaning: Millions of tiny air bubbles are generated within liquid by high frequency vibration. The air bubbles burst when in contact with object and dislodge the debris to achieve the cleaning effect.

- Using tap water is sufficient. Purified water or distilled water has the same cleaning effect as regular tap water for ultrasonic cleaning.
- When cleaning silver or copper items where oxidation has darkened the items, special solutions needs to be added to the water to remove the oxidation.

Main features

- Tank opening 32.5 x 20.5 x 9.9 cm / 12.7"x8.1"x9.9". Longest item that can fit inside the tank is 34.0 cm / 13.4". 6.0 L / 1.6 gallon tank capacity.
- Multiple options for accessories, suitable for various applications.
- 2 industrial grade ultrasonic transducers: 2x(55~60)W=(110~120)W, 2 ceramic heaters: 2x(78~87), W=(156~174)W.
- 3-color LED display, 4-key controls, 6 timer settings, 3 temperature settings, simple and straightforward to operate.
- Control panel with capacitive sensing technology. Durable , reliable, resistant to water and harsh chemicals.
- Electrical circuit features:
 1. Three transducers, each has an independent PCB control. The output power is increased by 40% compared to less expensive control method that uses one PCB to control multiple transducers.
 2. Multiple circuit protectors:
 - a. Thermal switch to prevent running dry. If temperature is too high, the switch cuts off power for approximately 6 minutes then reconnect the power.
 - b. Thermal protector to prevent PCB from overheating due to overload.
 - c. Timer to prevent heater from running continuously.
- Industrial grade IC. Moisture-proofed PCB and Cooling fan.
- Drainage valve.

Read the manual first

The manual should be carefully reviewed before starting to use the device. Warnings should be observed carefully. Please follow the manual for operations.

Table of contents

1. Features.....	20
2. Safety precautions.....	21-22
3. Items Not Suitable for Ultrasonic Cleaning	22
4. Product Structure and Accessories.....	23
5. Product Exploded View	24-25
6. Control Panel and Operations	26
7. Common Cleaning Methods.....	27-29
8. Operation Guidelines.....	29-30
9. Care and Maintenance	31
10. Optional Accessories.....	32-33

1. Features

Large tank capacity	Tank capacity 6.0L / 1.6 gallon; Tank size 32.3x20.5x9.9cm
Control panel with capacitive sensing technology	Durable, reliable, resistant to water and harsh chemicals.
Drainage valve	Convenient to use.
Two industrial grade ultrasonic transducers	Uniform distribution of ultrasonic waves, strong cleaning, high durability.
Two ceramic heaters	Mounted underneath the bottom of the tank for safe and even heating; durable construction.
3-color LED display	Easy to operate, with 4 key control panel, 6 timer settings, 3 temperature settings.
Independent control circuit for each ultrasonic transducer	40% higher efficiency, not susceptible to interference, more durable.
Multiple circuit protectors	When overloaded or improperly used, the protectors shut down the power to certain areas to protect the machine.
Moisture-proofed PCB with industrial IC	PCBs are placed vertically for better moisture-proof performance. Better anti-interference.
Two cooling fans	One for cooling of PCB boards, one for cooling of ultrasonic transducers. Improve heat dissipation and beneficial for continuous operations.
Multiple optional accessories	Suspendible plastic basket, stainless steel rack and tray, indirect cleaning tank, beaker and double beaker holder.
Solid and advanced structural design	4.5 mm thick wall housing made from engineering grade plastic with tongue in groove design for better water-proof and drop-proof performances. 60% quieter than steel housing designs, suitable for indoor uses.

2. Safety precautions

Keep it away from children!

This appliance can be used by children aged from 8 years and above if they have been given supervision or instruction concerning use of the appliance in a safe way and if they understand the hazards involved. Cleaning and user maintenance shall not be made by children unless they are aged from 8 years and above and supervised. Keep the appliance and its cord out of reach of children aged less than 8 years.

- Please store the ultrasonic cleaner where it is not reachable by children.
- Danger to children! Danger for death through suffocation! Keep the packaging material away from children. This appliance shall not be used by children. Keep the appliance and its cord out of reach of children.

To prevent life-threatening electrical shock, please observe the following:

Danger of electrical shock! Do not use while bathing. Never immerse the device or the power cord in water or other liquid.

- Danger of electrical shock! Never touch the power plug with wet hands, especially when inserting or removing the plug.
- Danger of electrical shock! If the unit has fallen into water during operation, do not touch the unit. Remove the power plug from the socket first.
- Danger of electrical shock! Do not spray water or liquid over the device.
- Never operate the device unattended.
- Follow the manual to operate the device.
- Do not use components unapproved by the manufacturer.
- When removing the power cord from the socket, grab the power plug not the cord.
- To protect the power cord from damage, do not cause it to get caught by things such as a cupboard door or a chair leg; do not drag across a hot surface.
- If there is damage to the power plug, cord, housing, or other parts of the device, do not use the device.
- Do not disassemble the device, except by professionals.
- If the unit is damaged, non-operational or has fallen into water, take it to a qualified service provider.
- Remove the power plug from the socket:
 - if malfunction occurs
 - before cleaning the device
 - if the device is not going to be used for prolonged period
 - after each use (recommended).
- The installation of an earth leakage circuit breaker with a rated tripping current of no more than 30 mA provides further protection against an electrical shock. The installation should only be carried out by a trained electrician.

To prevent fire hazards, please observe the following:

- Never block the vents on the device. Keep the vents free from lint, hair and other materials.
- Do not place the device on a soft surface, such as a bed or a couch, where the vents could be blocked.
- Observe the other warnings in the previous section.
- If the supply cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid a hazard.

Other observations:

- Do not operate the product without filling the tank with water. Running dry will damage the unit.
- Do not plug in the power cord before adding water to the tank. Do not fill the tank above the Max line to avoid overflow.
- Do not use solution containing abrasive substances or strong corrosive chemical solution not recommended by the manufacturer or the supplier.
- Place the device on a dry and flat surface for operation.
- When the device is subjected to severe electromagnetic interference, it may malfunction, stop operating or lose control functions. If this happens, unplug the power cord then reinsert it to restart the device.

Important safety instructions


When using electric appliances, basic precautions should always be followed, including the following:

- Read all the instructions before using the appliance.
- To reduce the risk of injury, close supervision is necessary when an appliance is used near children.
- Only use attachments recommended or sold by the manufacturer.
- Do not use outdoors.
- To disconnect, turn all controls to the off (O) position, then remove plug from outlet. Do not unplug by pulling on cord. To unplug, grasp the plug, not the cord. Unplug from outlet when not in use and before servicing or cleaning.
- To reduce the risk of electrical shock, do not put the appliance in water or other liquid. Do not place or store appliance where it can fall or be pulled into a tub or sink.
- All servicing of this product, including transducer replacement, is to be conducted by qualified service personnel.
- Do not operate any appliance with a damaged cord or plug, or after the appliance malfunctions or is dropped or damaged in any manner. Return appliance to the nearest authorized service facility for examination, repair, or electrical or mechanical adjustment.

3. Items Not Suitable for Ultrasonic Cleaning

<p>Soft Jewelry: Pearls, emerald, ivory, coral, agate, sea turtle shells, etc.</p>	<p>These items are not hard, so scratches may occur during cleaning.</p>
<p>Welded, Plated and Glued Items: Welded or plated metal items, glued items</p>	<p>Ultrasonic cleaning may enlarge the gaps inside the welded joints, plated coating or glued items and may cause separation.</p>
<p>Watches: Except diver's watches with depth rating over 50 m (150 ft).</p>	<p>Because the strong penetration capability of the ultrasonic waves, water may get into the watches if they are not truly waterproof. Use the watch stand supplied as a precaution if not sure.</p>
<p>Others: Wood; coated glass, ceramic, camera filters with preexisting cracks.</p>	<p>Ultrasonic cleaning may enlarge the cracks pre-existed in the coating on the glasses, ceramic, and glass. If the items have no pre-existing cracks, then it is okay.</p>


4. Product structure and accessories

Standard accessories:


Power cord


Hose



Optional accessories:


Indirect cleaning tank


Stainless steel rack and tray

Suspendible plastic basket




Beaker and beaker holder

Plastic basket


5. Product exploded view

No.	Item	Material
1	Top handle	ABS
2	Lid-window	AS
3	Lid-top	ABS
4	Lid-bottom	ABS
5	Top ring	ABS
6	Stainless steel handle	SUS 304
7	Basket foldable handle 1	ABS
8	Basket foldable handle 2	ABS
9	Plastic basket	ABS
10	Silicon seal	Silicone rubber
11	Stainless steel tank	SUS 304
12	Left cover	ABS
13	Left handle	ABS
14	Left panel	ABS
15	Inner tank	ABS
16	Front cover	ABS
17	Spout	SUS 304
18	Washer	AL
19	Nut	AL
20	Hose	Silicone rubber
21	Nut	ABS
22	Spout	ABS
23	Ball valve	Standard part
24	Bracket	ABS
25	Elbow	ABS
26	Valve handle	ABS
27	Bottom cover	ABS
28	Power PBC 1	94v0
29	Power PBC 1	94v0
30	Cooling fan	Standard part
31	Right handle	ABS
32	Right panel	ABS
33	Right cover	ABS
34	Transducer PBC	94v0
35	Right panel	ABS
36	Power switch	94v0
37	Power switch	Standard part
38	Power socket	Standard part
39	Cover holder	ABS
40	Back cover	ABS
41	Heater cover	PA66
42	Heater pocket	PA66
43	Thermocouple holder	ABS
44	Thermal cutoff	A0-3A-F
45	LED	Fire retardant ABS
46	Control panel box	ABS
47	Control panel cover	PMMA

6. Control panel and operations

- a. **10:00 – Working Time Display:** it counts down after work is started.
- b. **Multistage Thermometer:**
-  Solid green light indicates that the water temperature is 0-30% of the set temperature.
 -  Solid orange light indicates that the water temperature is 30-60% of the set temperature.
 -  Solid red light indicates that the water temperature is 60-90% of the set temperature.
 -  Flashing red light indicates that the water temperature is close to the set temperature. All lights indicate that the water temperature is at the set temperature (90-100%).
- c. **060 – Set Temperature Display:** press Temp button to select one of the 3 set temperatures.
- 50°C → 55°C → 60°C
→ 122°F → 131°F → 140°F
- d. **CF – °C or °F indicator:**
- °C: set temperature in °C.
 - °F: set temperature in °F.
- e. **▲ – Overheat protection warning light:** if the unit is working improperly or has been running for too long, the red warning light will illuminate, indicating the unit is under overheat protection. The unit can not be operated until idling for about 15 to 20 minutes. The warning light then turns off and the operation can be resumed.
- f. **Heater status:** press Heater button,  flashes, indicating the heater is turned on. Press Heater button again,  dims, indicating the heater is turned off.
- g. **Normal cleaning status:**  illuminated, normal cleaning is underway.
- h. **Heater button:** pressing Heater key twice,  flashes, indicating the heater is turned on. To cancel heating, press Heater key twice again,  dims, indicating the heater is turned off.

 To prevent accidentally turning on or off the heater, the Heater key needs to be touched twice within 6 seconds to turn it on or off  flashes, indicating the heater is on.

When the water reaches the set temperature, the heater will turn off automatically.

The ceramic heater used in this unit has the advantages of having small size, high heat and longevity. To avoid damages to the heater, do not add cold water when the heater is hot otherwise ceramic may break due to large temperature differences. Avoid the following misuses of the heater:

1. Drain dirty water out then add cold clean water when the heater is on.
2. Turn on the heater when the tank is empty: heater can be damaged if it is run for more than 15 seconds without water in the tank. Always turn off the heater before adding cold water. Make sure there is enough water in the tank before turning on the heater.

- i. **Temperature Set Button:** 3 set temperatures (°C/°F) can be selected. Press Temp button to select the temperature in the following sequence:

→ 50°C → 55°C → 60°C
→ 122°F → 131°F → 140°F


- j. **On/Off Button:** After powering on, the LED displays **10:00**, which is the default time (as shown), and the unit is ready with the normal working status. Press On/Off key twice once and the cleaning starts. The cleaning stops when the timer counts down to **00:00**. If the unit needs to be stopped before the timer runs out, press On/Off key twice.

 To prevent accidentally turning on or off the unit, the On/Off key needs to be touched twice within 6 seconds to turn it on or off.  flashes, indicating the unit is running under Normal Cleaning mode.

- k. **Timer Quick Set Button:** press Time button, display shows **10:00**. Each pressing increases the time by 5 minutes.

7. Common cleaning methods

REGULAR CLEANING

Only tap water is needed.


Suitable Locations: Jewelry makers, optical labs, biology / chemical labs, homes, speech therapy clinics, kindergartens.

Cleaning method

1. Put the items to be cleaned inside the basket and put the basket inside the tank.
2. Add water to the cleaning tank to a level between "MIN" and "MAX" lines and above the area to be cleaned.
3. Turn the power switch on. Press Time button to set the timer to 5 -10 minutes. Press On/Off button to start cleaning.

Notes on using the basket:

1. Basket can reduce the friction between the items and the tank. But stainless steel basket absorbs approximately 8 % of ultrasonic energy, plastic basket absorbs about 35 % of ultrasonic energy. This significantly reduces the cleaning effectiveness.
2. During cleaning, dirt will come off like smoke and the water will become murky over time. When "smoke" stops coming, the cleaning is basically done.

INDIRECT CLEANING

Put the items in a separate container. Use ultrasonic waves to penetrate the container to do the cleaning.

Suitable Locations: medical instrument sterilization rooms, precision electrical parts manufacturers, watch and precision parts manufacturers, mobile phone and electronics service shops, large printing shops.

Different industries use different fluids to be used inside the container:

1. Watch service shops - Watch oil (to prevent rust)
2. Medical instrument sterilization – Enzymatic solution (cleaning prior to sterilization)
3. Precision electrical parts – Hexane (to dissolve grease, to evaporate quickly)
4. Mobile phone and electronic service shops - Pure alcohol (to evaporate quickly)
5. Printing shops – Acetone are often used(to dissolve ink)
6. Homes - Rubbing alcohol (to remove odor, to clean and to disinfect at the same time)


Avoid acetone from contacting the plastic housing to prevent corrosion. Use Indirect Cleaning and a glass or metal container for acetone.

Cleaning method:

1. Put the item inside a container. Add proper fluid to submerge the area to be cleaned.
2. Place the container in the tank directly. Add water to a level between MIN and MAX, not over the container.
3. Turn the power switch on. Press Time button to set the timer to 10 minutes. Press On/Off button to start cleaning. Ultrasonic waves will penetrate the container and clean the items.

Selection of container:

1. Plastic containers: plastics are soft and will absorb about 30-40% of ultrasound energy and reduce the cleaning effect.
2. Aluminum containers – Absorbs about 20% of ultrasonic energy.
3. Glass containers– Absorbs about 15% of ultrasonic energy.
4. Stainless steel containers – Absorbs about 8% of ultrasonic energy.

**ENHANCED CLEANING**

Debris accumulated over a long period of time, greasy or heavily soiled.

Suitable Locations: metal processing and jewelry makers, hardware manufacturers, shooting clubs, homes, golf clubs, automotive repair shops.

Cleaning method:

1. Large items can be put in the tank directly. Small items can be put in the basket then put in the tank. Do not put the items on top of each other to avoid rubbing during cleaning.
2. Add water to a level between MIN and MAX and above the area to be cleaned.
3. Add 50 ml (3 tablespoons) of dishwashing liquid.
4. Turn the power switch on. Press Temp button to set the water temperature to 55°C (131°F). Press Heater button to turn on the heater,  flashes. Close the lid. LED  lights up. When water reaches the set temperature,  dims, indicating the heater is turned off.

Warm water and dishwashing liquid can soften grease and improve cleaning efficiency.

5. Press Time button to set the timer to 15 minutes. Press On/Off button to start cleaning. Grease will dissipate and appear to be like black smoke in the water.
6. When it stops, open the drainage valve to release the dirty water. Close the valve afterwards.
7. Clean the tank and add fresh water. Wash the items for another 3 minutes to remove the residual debris and the detergent.


SPECIAL CLEANING

Silver, copper or brass items with oxidation which has darkened the items will need special cleaning.

Suitable Locations: silver, copper or brass product manufacturers, home, shooting clubs.

Cleaning method:

1. Put the items in the basket and then in the tank, add water to a level between MIN and MAX and above the area to be cleaned.
2. Add special solution that can remove silver or copper oxidation according to the ratio recommended (typically 1:10). Use special solutions according to the instructions for brass cartridges.


- Turn the power switch on. Press Temp button to set the water temperature to 50°C (122°F). Press Heater button to turn on the heater, flashes. Close the lid. LED lights up indicating the percentage range of the set temperature being reached. When water reaches the set temperature, dims, indicating the heater is turned off.

The heater has dual protection. It will be turned off automatically after 50 minutes.

- Press Time button to set the timer to 15 minutes.
- When it is done, remove the basket and the items. Rinse the items with tap water or follow steps 6 and 7 in the Enhanced Cleaning section.


During cleaning, debris will come off the items like "smoke". Water will become murky. This method will remove silver or copper oxidation and debris in the decorating patterns and crevices and make the items like new again.

8. Operation guidelines

- Remove the cover and put it on the cover holder vertically. Put items in the basket then put them in the cleaning tank. Add water to a level between MIN and MAX and above the area to be cleaned.

If the unit is turned on without water, ultrasonic energy will not be absorbed. This model also has high power. Once on for over 15 seconds, it may damage the unit or severely reduce the life of the unit.




- Connect the power cord to an outlet and turn the switch on. LED display shows 10:00. This is the most common timer setting. If the timer needs to be adjusted, press Time button. Select desired time setting by pressing Time button:

→ 10 min → 15 min → 20 min → 25 min → 30 min → 5 min

Prolonged cleaning time may result in:

- Loosening of the screws if used.
- If items are plated and have existing cracks, the cracks may become more pronounced.
- Coating with preexisting damages may experience more peeling.

3. Press On/Off button to start cleaning.

During cleaning, buzzing sound can be heard from the tank. Closing the lid will reduce the noise level.

The digital timer will count down to show the remaining cleaning time. When it displays **00:00** the cleaning is done. To stop cleaning at any time, press On/Off button.

The unit is equipped with overheat protector. If it has been running for too long, the red warning light **△** will illuminate, indicating the unit is under overheat protection. The unit can not be operated until idling for about 15 to 20 minutes. The **△** turns off and the operation can be resumed.

4. Select one of the four cleaning methods recommended earlier.

- a. When the heater is needed, press TEMP button to select one of the 3 preset temperatures. Press Heater button, **🔥** flashes, indicating the heater is turned on. When the water reaches the set temperature, **🔥** dims, indicating the heater is turned off.


To ensure safety, the unit is designed with dual protections. The heater will be turned off automatically if it has been running for 50 minutes.

	<p>Thermometer lights up to indicate the percentage of the set temperature reached.</p>
--	---

- b. When debris “smoke” is no longer visible, cleaning is done. If additional cleaning is needed, reset the timer and repeat the steps above.


5. When cleaning is completed, switch off the power, open the lid and retrieve the basket and the items. Alternatively, raise and suspend the basket above water as illustrated in the Optional Accessories section to drain water. Connect the hose to the drainage pipe as shown. Open the drainage valve to drain the dirty water. Clean and wipe the tank dry. Close the drainage valve.

9. Care and maintenance


1. **Do not turn on the unit without water in the tank.** Even though the unit is designed with multiple protections, if it is turned on for over 30 seconds without water in the tank, it may damage the unit or severely reduce the life of the unit.


⚠ A thermal switch is attached to the bottom of the cleaning tank. If the unit is turned on accidentally without water in the tank, the ultrasonic energy will not be transferred out through the liquid and instead will be concentrated on the tank, the temperature of the tank will continue to rise. If the heater is turned on at the same time, the temperature of the tank will rise rapidly. This will trigger the thermal switch to cut off power to stop the ultrasonic transducer and the heater from running, hence to prevent damages to the unit and fire hazard. The switch cuts off power for approximately 6-10 minutes then reconnect the power automatically.


2. **Do not run the unit for extended time or continuously.** The unit is designed with overheat protection. If the unit has been running for 45 minutes, it is recommended to stop the unit for about 20 minutes to prolong the life of the unit.

3. **Do not keep water in the cleaning tank for a long time.** After cleaning is completed, open the drainage valve to drain the dirty water. Clean and wipe the tank dry.


4. **Do not spray water over the housing.** Use a towel to wipe the tank and the housing dry.

5. **Do not expose the unit under direct sunshine for a long time.** Keep the unit in a dry, cool and ventilated area.


10. Optional accessories

SUSPENSIBLE PLASTIC BASKET


The basket has foldable handles. It can be raised and suspended above water to drain water.

STAINLESS STEEL RACK AND TRAY


The rack can be used to support cassettes. 2 layers of cassettes can be used on top of each other. The tray can be used in place of the cassette. The tray is stackable so multiple layers of trays can be used.

INDIRECT CLEANING TANK


Place the indirect cleaning tank across the top of the main cleaning tank as shown. Add water inside the main tank so the water level reaches the indirect cleaning tank. Add solution to the indirect cleaning tank and clean the items inside. Cleaning efficiency can reach 90% inside the Indirect Cleaning Tank due to our technology and quality. It can be used to isolate the cleaning and the solution from the main cleaning action. It also reduces the amount of the solution needed.

BEAKER AND BEAKER HOLDER

Place the beaker holder with beaker(s) across the top of the main cleaning tank as shown. Add water to the main tank so the water level reaches the beaker(s). Add solution inside the beaker(s). It can be used to clean small items such as dental burs and files, nozzles, etc.

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