## **CORE D & CORE D FLOW**

### **Description/Indications**

Dual-curing highly resistant composites for core build-up.

Core D & Core D Flow are available in a tooth-colored shade which is ideal in situations where esthetics of the core is important. Permanent cementation of endodontic fiber-posts.

Core D Flow: permanent cementation of crowns, bridges, inlays, onlays (ceramic, metal, and composite), in conjunction with Healbond Duo SE, self-etching bonding system.

### **Physical data/Characteristics of the material**

Working time: 90 seconds

Intra oral setting time + working time (self-curing mode): 3 minutes 10 seconds

Flexural strength: 120 Mpa

Compressive strength: 300 Mpa

Radiopaque

### **Safety notes**

#### Warnings

Core D & Core D Flow contain polymerizable methacrylate resins which may irritate to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.

Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.

Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

#### Precautions

Use of rubber dam or adequate isolation is recommended.

Wear suitable protective eyewear, mask, clothing and gloves.

Protective eyewear is recommended for patients.

### **Contraindication**

Core D & Core D Flow are contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

### **Instruction for use**

#### Core built up

- Prepare the cavity in the usual manners.
- Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond Duo (with etching) or Healbond Duo SE (self-etching).
- Once the surfaces have been properly treated, they must be kept uncontaminated.
- Remove syringe or cartridge cap, dispense and discard a small amount of material from the dual-barreled syringe (Core D, Core D Flow), or from the cartridge (Core D). Be sure the material is flowing freely from both ports.
- Install a mixing tip + an intra-oral tip on the cartridge (Core D), or on the dual-barreled syringe (Core D, Core D Flow).
- Core D - Insert the dual-barreled syringe, or the cartridge, in the dedicated dispensing gun. Gently depress the trigger of the dispensing gun, to begin the flow of material. Do not use excessive force.
- Core D Flow - Gently depress syringe plungers, to begin the flow of material. Do not use excessive force.
- Apply Core D or Core D Flow in 2mm layer, directly in the prepared cavity, with a banded matrix retainer or clear crown form in place, if necessary (Core D Flow).
- Polymerize each layer for 30 seconds at a light intensity of at least 600 mW/cm<sup>2</sup> or 15 seconds at a light intensity of at least 1100mW/cm<sup>2</sup>.
- Once has set, use rotary instruments to complete the crown preparation (final shaping and finishing) in the same clinical trial.
- After use, leave the mixing tip in place: it will protect the nozzle until the next application.

#### Cementation of endodontic fiber-post

- Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond Duo (with etching) or Healbond Duo SE (self-etching). If core build-up will immediately follow post cementation, treat all available tooth structure beyond post space with adhesive for core build-up bonding.
- Treat surface of post according to adhesive manufacturer's and/or dental laboratory's instructions.
- Bleed material and discard it. Put a mixing tip + an intra-oral tip on the nose of dual-barreled syringe or dual-barreled cartridge.
- Gently depress syringe or cartridge plungers, to begin the flow of material. Do not use excessive force.
- Dispense the same material directly from the syringe or from the cartridge, onto the post surface and into the post hole preparation.
- Immediately place post into final position and stabilize.
- Once stabilized, light-cure, using a visible light-curing unit, all accessible areas of the post for 20 seconds, with a minimum output of 600 mW/cm<sup>2</sup>.

- Proceed with the core build-up.

**Interaction**

Eugenol-containing dental materials should not be used with Core D & Core D Flow because they may interfere with hardening and cause softening of the polymeric components of the materials.

**Storage & shelf-life**

Store at a temperature not exceeding 23°C. When stored under such conditions, Core D & Core D Flow have a shelf-life of 24 months. Refrigerate when not in use (for example, overnight and on weekends). When cold, these materials have a stiffer consistency. For easier handling, remove from refrigerator at least 15 minutes prior to use.

**Presentations & references**

**Core D** - 25ml dual cartridge + 20 yellow mixing tips + 20 intra-oral tips

**CDD-25** (dentin)

**CDB-25** (white)

**Core D** - 10ml dual syringe + 10 brown mixing tips + 10 intra-oral tips

**CDD-10** (dentin)

**CDB-10** (white)

**Core D Flow** - 2 x 2.5ml dual syringe + 10 brown mixing tips + 10 intra-oral tips

**CDF-5** (dentin)

**Important notice**

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

**Disclaimer**

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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