Temporary crown and bridge resin

SELF-CURING

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EN Instructions for use.

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DE

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INSTRUCTIONS FOR USE

EN

SELF-CURING RESIN FOR TEMPORARY CROWNS AND BRIDGES

Definición

Temporary crowns and bridges material available in 50ml cartridge.

Physical data/Characteristics of the material

Intra oral setting time + working time : 2' 25" - 2' 40"

⚠ Safety notes

Warnings

- Temporary crown and bridge resin contains polymerizable methacrylates which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindicación

Temporary crown and bridge resin is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

Instructions for use

- Take a pre-operative impression. Interdental spaces should be rid of silicone material.
- Prepare the abutment teeth.
- Remove the cap from the cartridge and check that the two exits are not obstructed.
- Attach a mixing tip (yellow ring type) on the cartridge.
- Place the cartridge in the mixing gun.
- When using the cartridge for the first time, bleed a small amount of material onto a mixing pad to ensure that both base and catalyst will be extruded equally. Discard the purged amount.
- Blow dry impression and fill the concerned teeth area, with Temporary crown and bridge resin. Always fill from the bottom up, until the impression is full. It is recommended to inject also some resin into the sulcus of the concerned teeth, with a yellow intra-oral tip, to avoid any lack of product, at the limits of the temporary crown.
- Place the impression back in patient's mouth and allow Temporary crown and bridge resin to set undisturbed for about 40/45 seconds.
- Remove the impression tray from mouth. At that stage, the material will be somehow flexible and will pull over any undercut or non-parallel preparations.
- The material will fully cure in 2 additional minutes.
- Remove cured prosthesis from the impression and trim excess flash.
- Remove smear layer on the surface (if any) by whipping with a cloth moistened with alcohol.
- Polish, if desired, using preferably pumice or corundum filled polishing paste.
- Cement the temporary prosthesis with a temporary cement.
- After use, leave the mixing tip in place: it will protect the nozzle until the next application.

Storage & shelf-life

Store at a temperature not exceeding 23°C. When stored under such conditions, Temporary crown and bridge resin has a shelf-life of 24 months. Refrigerate when not in use (for example, overnight and on weekends). When cold, these materials have a stiffer consistency. For easier handling, remove from refrigerator at least 15 minutes prior to use.

Presentations & references

- 50ml dual cartridge + 10 yellow mixing tips
- 69400 (shade A2)
- 69401 (shade A3)

⚠ Important notice

Do not throw away the box and the instruction for use until the product is completely finished. High humidity could damage the packaging on which important instructions and legal information are printed. Storage at a too high temperature could cause to age prematurely and inhibit its polymerization. Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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