PR SCELL GLASS-IONOMER

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PR SCELL GLASS-IONOMER

PR SCELL GLASS-IONOMER is a radiopaque glass ionomer cement recommended for cementing crowns, bridges, inlays or orthodontic appliances securely in place.

Special features include: excellent handling properties, low film thickness, less moisture sensitivity during the initial setting stage and prolonged working time followed by snap-set reaction. It releases fluoride.

1. Powder/liquid proportion

One level scoop of powder with two drops of liquid gives a powder to liquid ratio of 2.0g powder to 1.0g liquid.

2. Dispensing

Use the green powder scoop and dispense one level scoop of powder (using the ledge inside the powder bottle to level off) onto a glass slab or mixing pad. Remove cap from liquid bottle and dispense two drops of bubble-free liquid (remove any air-bubbles trapped in the neck of the bottle by holding it vertically and gently squeezing the bottle, beforehand) near the powder, then replace the caps back onto the bottles.

3. Preparation

Thoroughly clean the preparation with oil-free pumice paste. Rinse with water and air dry. If there is near, or total, exposure of the pulp, line with calcium hydroxide before cementation. Clean the surfaces (of the post, crown, bridge, etc.) that are to be cemented to the tooth structure, in an ultrasonic bath, rinse with water, wipe with plegget moistened in alcohol, and then air dry.

4. Mixing/Mixing times

Mix on a mixing pad or glass slab at room temperature (21-25°C/70-77°F), in a humidity range of 50% +/- 10%. Mix together the dispensed powder and liquid until homogeneous, by applying short, rapid strokes with a plastic spatula. Because of the excellent integration of the powder and liquid, it is not always necessary to divide the powder up into portions for gradual mixing. Total mixing time should be 30 seconds.

5. Placement

Use the spatula to apply a thin lining of cement to the internal surface of the crown. Do not overfill with cement. Seat the interfacing surfaces and have the patient bite down reasonably hard on a cotton ball to remove air and excess mixed cement. Care should be taken to protect the cement from moisture (water or saliva) particularly at the margins.

6. Working and setting times

Working time is about 3.5 minutes from start of mix. The net setting time is 4 minutes from end of mix. For a longer working time, pre-chill the powder and liquid and use chilled glass slab to mix on.

7. Storage

Since this cement is slightly hygroscopic, reseal the powder and liquid bottles immediately after use. Store both in a cool place and do not subject to excessive sunlight or heat.

Note:

- Do not use after expiry date.
- This product is only to be used for recommended indications.

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