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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Parameters determined at time of manufacture according to ISO 11138-1:2006, ISO 11138-3:2006 and IRAM 37102:1999 (Parts 1 and 3). The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier. It also has a purple growth indicator medium contained within the glass ampoule.

Product description

Bionova® BT20 Biological Indicator is specifically designed for monitoring vacuum assisted Steam sterilization processes at 121 - 134 °C and gravity air displacement Steam sterilization cycles at 121 - 134 °C. If the sterilization process was not successful, the indicator medium will turn from purple to yellow after incubation at 60±2 °C, thus indicating the presence of living *Geobacillus stearothermophilus* spores. If the sterilization process was successful, indicator medium will remain purple after incubation. The final readout should be performed at 24 hours of incubation at 60±2 °C.

Warning!

Do not use BT20 Biological Indicators for monitoring Ethylene Oxide, Dry Heat, Chemical Vapor, Radiation or other sterilization processes. Do not re-use biological indicators.

Storage

Store in a dark place at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze.

Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

- Identify the Bionova® BT20 indicator by writing the sterilizer number (in case of having more than one), load number and processing date on the indicator label.
- Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in areas which a priori you consider most inaccessible for the sterilizing agent (high temperature steam). Generally, the center of the load and areas near the door are problematic.
- Sterilize as usual.
- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the biological indicator from the package. **CAUTION!** Wear safety glasses and gloves when removing the BT20 biological indicator from the sterilized package. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
- Allow the biological indicator to cool down to room temperature.
- Check the chemical indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
- Crush the glass ampoule contained in the biological indicator and incubate at 60±2 °C.
- IMPORTANT:** Use a non-sterilized biological indicator as a positive control to ensure that appropriate incubation conditions were met.
- Incubate the processed biological indicator and the indicator used as a positive control for a maximum of 24 hours at 60±2 °C. Readout should be carried out at convenient intervals of 10 hours.

Color change of the growth indicator medium from purple to yellow means that a sterilization process failure has occurred. If at 24 hours there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was efficient). The positive control indicator should show a purple to yellow color change for the results to be valid.

NOTE: To extend the incubation time beyond 24 hours, we recommend using a humidified environment to prevent complete evaporation of the culture medium contained within the indicator.

Record the positive results and discard the biological indicators immediately as described below.

WARNING! Do not use the sterilizer until the biological indicator test results are negative (indicator medium of the processed tube remains purple).

Disposal

Discard biological indicators after use, according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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