**1 Intended use****Intended purpose**

Adhesive cementation of indirect anterior and posterior restorations

**Patient target group**

- Patients with permanent teeth

**Intended users / Special training**

- Dentists
- No special training required.

**Use**

For dental use only.

**Description**

Variolink® Esthetic DC is a colour-stable, adhesive luting composite for the permanent cementation of dental restorations and endodontic posts. The special filler composition gives Variolink Esthetic DC a very high radiopacity (300% Al).

The dual-curing Variolink Esthetic DC is suitable for

- the permanent cementation of glass-ceramic, lithium disilicate glass-ceramic and composite resin restorations (inlays, onlays, partial crowns, crowns, bridges)
- the permanent cementation of opaque restorations made of e.g. oxide ceramics or metal, if they are used in combination with a dual-curing adhesive or an adhesive that is separately light-cured
- the permanent cementation of endodontic posts in combination with a dual-curing adhesive.

The following overview provides recommendations for the selection of the most suitable Variolink Esthetic DC shade:

Shade	Effect	Area of use	
		Inlays / Onlays	Crowns / Bridges
Light +	considerably lightens	—*	✓
Light	slightly lightens	✓	✓
Neutral	no shade effect	✓	✓
Warm	slightly darkens / more yellowish	✓	✓
Warm +	considerably darkens / more yellowish	—*	✓

\* The use of the shades Light+ and Warm+ may lead to visible restoration margins when cementing inlays/onlays.

# Variolink® Esthetic DC

[sv]	<b>Bruksanvisning</b> Dualhärdande dental kompositcement	[no]	<b>Bruksanvisning</b> Dualhärdende odontologisk sementeringskomposit
[da]	<b>Brugsanvisning</b> Dualhärdende resinretentionscement	[nl]	<b>Gebruiksaanwijzing</b> Tandheelkundig, dual uithardend bevestigingsmateriaal op basis van kunsthar
[fi]	<b>Käyttöohjeet</b> Kaksioiskovetäinen yhdistelmämuovisementti		

# Variolink® Esthetic DC

[en]	<b>Instructions for Use</b> Dual-curing resin-based dental luting material	[it]	<b>Istruzioni d'uso</b> Cemento composito dentale ad indurimento duale
[de]	<b>Gebruuchsinformation</b> Dualhärtendes zahnärztliches Befestigungskomposit	[es]	<b>Instrucciones de uso</b> Material dental de cementado en base a resina de polimerización dual
[fr]	<b>Mode d'emploi</b> Composite de collage à prise duale à usage dentaire	[pt]	<b>Instruções de Uso</b> Material para cimentação à base de resina de cura dual

## Technical data

### Classification acc. to ISO 4049

- Type 2 / Class 3
- Dual-curing luting composite that is light-activated intraorally
- Optionally, Variolink Esthetic DC can be cured with light in the wavelength range of 400–500 nm.

### Shades

Variolink Esthetic DC is available in the following five shade gradations:

Shade		Effect
Light +	Translucency approx. 5% (opaque)	Lighter / more whitish
Light	Translucency approx. 10%	✓
Neutral	Translucency $\geq$ 17%	None
Warm	Translucency approx. 12%	✓
Warm +	Translucency approx. 8.5% (opaque)	Darker / more yellowish

The gradation of Variolink Esthetic DC shades is based on the effect a cement exerts on the brightness value of the final restoration. Variolink Esthetic DC Neutral does not affect the brightness value. At the same time, it demonstrates the highest translucency and is thus neutral in shade. Light and Light + make the restorations lighter, while Warm and Warm + create a darker overall shade impression. In order to evaluate the overall effect of the restoration in conjunction with the various Variolink Esthetic DC shades prior to permanent cementation, we recommend using the Variolink Esthetic Try-In Pastes.

### Working time

The working and curing times depend on the ambient temperature. The following times apply as soon as Variolink Esthetic DC has been extruded from the automix syringe:

When used purely self-curing	Room temperature 23 °C $\pm$ 1 °C 73 °F $\pm$ 2 °F	Intraoral 37 °C $\pm$ 1 °C 99 °F $\pm$ 2 °F
Working time	approx. 5 min	approx. 2 min
Curing time (incl. working time)	approx. 8 min	approx. 4 min

### Mixing ratio

Variolink Esthetic DC is always dispensed from the automix syringe in the optimum ratio.

### Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

### Areas of application

- Permanent cementation of glass-ceramic, lithium disilicate glass-ceramic and composite resin restorations (inlays, onlays, partial crowns, crowns, bridges)
- Permanent cementation of opaque restorations made of e.g. oxide ceramics or metal, if they are used in combination with a dual-curing adhesive or an adhesive that is separately light-cured.
- Permanent cementation of endodontic posts in combination with a dual-curing adhesive.

### Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

### Limitations of use

- Do not use if a dry working field cannot be established or the prescribed application technique cannot be applied.
- The product is not intended to be reprocessed or reused.

### Side effects

Systemic side effects are not known to date. In individual cases allergic reactions to individual components have been reported.

### Interactions

- Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, do not disinfect the preparation using oxidative agents.
- Alkaline jet media (e.g. Airflow) on dentin can compromise the effect of self-etching adhesives.
- Hemostatic substances may inhibit polymerization and/or lead to discolouration. Therefore, the directions for use of these substances must be observed.

### Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

### Composition

Si-Zr mixed oxide, ytterbium trifluoride, UDMA, aromatic methacrylate, GDMA, aromatic-aliphatic UDMA, D3MA

The total content of inorganic filler is approx. 38 vol%.

Particle size of inorganic fillers: between 0.15  $\mu$ m and 15.5  $\mu$ m.

## 2 Application

For more detailed information also refer to the instructions for use of the products used in conjunction with Variolink Esthetic DC.

### 2.1 Removal of the temporary restoration and cleaning of the cavity

Remove any possible residue of the temporary luting cement from the cavity or the preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl<sup>®</sup> fluoride-free prophylactic paste). Rinse with water spray. Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

**Note:** Cleaning with alcohol can lead to dehydration of dentin.

## 2.2 Try-in of the restoration and isolation

Seat the restoration using the desired Variolink Esthetic Try-In Paste and check the shade, fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make adjustments with fine diamonds at medium speed and with light pressure and adequate water cooling.

Polish ground surfaces.

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam Plus, is required.

If the dental hard tissue is contaminated with blood or saliva during the try-in procedure, it needs to be cleaned again as described in section 2.1.

## 2.3 Pre-treatment of the restoration

2.3.1 If a restoration which has been pre-treated in the dental laboratory is tried in, the bonding surface of the restoration needs to be prepared as follows after the try-in:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Shake Ivoclean before use and completely coat the bonding surface of the restoration with the product using a microbrush or brush.
- Leave Ivoclean to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.
- Next, prime the bonding surface of the restoration with a suitable bonding agent (e.g. Monobond® Plus). Make sure to observe the Instructions for Use of the bonding agent used.

2.3.2 Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.

2.3.2.1 Subsequently, condition the bonding surface of the restoration as follows:

- a) Glass-ceramic restorations (e.g. IPS Empress®)
  - Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
  - Thoroughly rinse the restoration with water spray and dry with oil-free air.
- b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)
  - Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
  - Rinse the restoration thoroughly with water spray and dry with oil-free air.
- c) Restorations made from zirconium oxide (e.g. IPS e.max ZirCAD) or aluminium oxide ceramics
  - Sandblast the bonding surface (sandblasting parameters according to the instructions of the manufacturer of the restorative material).

Rx ONLY

CE 0123

**Manufacturer:**

Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

**Date information prepared:**

2022-06-29/Rev. 0  
746578/WE

**ivoclar**

**ivoclar**

- If necessary, clean the restoration in an ultrasonic unit for about one minute.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.
- **IMPORTANT!** In order to achieve an optimum bond, do not clean the zirconium oxide surfaces with phosphoric acid.

2.3.2.2 Apply Monobond Plus with a brush or a microbrush to the pre-treated surfaces, let it react for 60 s and then disperse with a strong stream of air.

2.3.3 Alternatively, restorations made of glass-ceramics and lithium disilicate glass-ceramics (e.g. IPS Empress, IPS e.max Press, IPS e.max CAD) may be conditioned with Monobond Etch & Prime® as follows:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil- and moisture-free air.
- Apply Monobond Etch & Prime on the bonding surface using a microbrush and scrub it into the surface for 20 s. Allow to react for another 40 s.
- Thoroughly rinse off Monobond Etch & Prime with water spray until the green colour has disappeared. If any residue remains in areas showing microporosity, the restoration may be cleaned with water in an ultrasonic unit for up to 5 min.
- Dry the restoration with a strong stream of oil- and moisture-free air for approx. 10 s.

2.3.4  Composite and fibre-reinforced composite restorations should always be conditioned according to the instructions of the manufacturer of the restorative material.

With respect to Tetric® CAD restorations, please observe the directions stated in the Tetric CAD Instructions for Use.

## 2.4 Pre-treatment of the preparation and application of the adhesive (e.g. Adhese® Universal)

 Please observe the instructions for use of the adhesive used.

2.4.1 Conditioning with phosphoric acid gel (optional)

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must be conditioned with phosphoric acid gel (e.g. Total Etch). Please observe the instructions for use of the phosphoric acid gel.

a) Selective enamel etching

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15-30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

b) Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then to the dentin. The etchant should be left to react on the enamel for 15-30 s and on the dentin for 10-15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

2.4.2 Application of Adhese Universal

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.  
**Important note:** Avoid pooling, as this may compromise the accuracy of fit of the final restoration.
- Light-cure Adhese Universal for 10 s using a light intensity of  $\geq 500$  mW/cm<sup>2</sup> (e.g. Bluephase® G4).

## 2.5 Seating of the restoration

- For each application place a new mixing tip on the syringe. Dispense Variolink Esthetic DC from the automix syringe and apply the desired quantity directly to the restoration. As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new tip before the next use).
- Variolink Esthetic DC should be processed quickly after extrusion from the automix syringe and the restoration seated in place. Apply Variolink Esthetic DC directly to the internal surface of the restoration.

**Note:** Direct application of Variolink Esthetic DC to a preparation that has been pre-treated with adhesive may - depending on the adhesive used - lead to an acceleration of the curing process and affect the accuracy of fit of the restoration.

- Seat the restoration and **fix/hold it in place during excess removal.**
- Remove all excess luting material.
  - a) Wiping technique  
Remove excess material immediately with a brush, dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).
  - b) Quarter technique (max. 6 bridge abutments) – Light-curing of excess and subsequent removal  
Light-cure excess material with the polymerization light (e.g. Bluephase G4) for 2 s per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of max. 10 mm. Thereafter, excess cement is easy to remove with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).

## 2.6 Polymerization

- As with all composite systems, Variolink Esthetic DC is subject to oxygen inhibition. In other words, the surface layer (approx. 50  $\mu$ m) does not polymerize during the curing process, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/Airblock (e.g. Liquid Strip) immediately after excess removal.

- Polymerize Variolink Esthetic DC in segments, starting with the proximal margins:

Light intensity	Curing time per mm ceramic and segment	
500 – 1,000 mW/cm <sup>2</sup>	20 s	
≥ 1,000 mW/cm <sup>2</sup>	10 s	e.g. Bluephase® G4

- Rinse off glycerine gel/Airblock (e.g. Liquid Strip).

## 2.7 Finishing of the completed restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
- Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraGloss®).
- If necessary, finish the restoration margins with suitable polishers (e.g. OptraGloss).

## Special instructions for the cementation of endodontic posts

- For the cementation of endodontic posts, carefully clean the root canal to remove any residue of root canal filling material. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite).
- Optionally, apply phosphoric etching gel (e.g. Total Etch) to the enamel portion of the preparation first, and then to the dentin portion of the preparation and the root canal. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Thoroughly rinse off the gel from the preparation surfaces. Thoroughly remove any residue of the etching gel from the root canal by rinsing with a syringe and endodontic irrigation cannula. Then dry with oil- and moisture-free air. Remove excess moisture from the root canal using paper points
- A microbrush is used to apply the mixed Adhese Universal DC from the Small/Endo Single Dose unit to the root canal and the surface of the preparation. The adhesive must be scrubbed into the tooth surface to be treated for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate. Disperse Adhese Universal DC with oil- and moisture-free compressed air until a glossy, immobile film layer results. Remove excess from the root canal using paper points.
- The endodontic post that has been prepared according to the instructions of the manufacturer is wetted with the mixed Variolink Esthetic DC.

**IMPORTANT!** Do not apply Variolink Esthetic DC with the lentulo into a root canal that has been wetted with Adhese Universal DC. Otherwise, premature curing might occur, which could possibly prevent the post from being properly positioned.

- Place the endodontic post. Excess cement will be displaced when the post is inserted. Distribute excess cement across the occlusal preparation surface for total coverage. Then light-cure Variolink Esthetic DC for 20 s.
- Core build-up material (e.g. MultiCore®) can then be directly applied onto Variolink Esthetic DC which serves as a bonding agent and cured according to the instructions of the manufacturer (light-curing).

## 3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclar.com](http://www.ivoclar.com)).
- Explanation of symbols: [www.ivoclar.com/eIFU](http://www.ivoclar.com/eIFU)
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.  
Basic UDI-DI: 76152082ACEME001DV

## Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website [www.ivoclar.com](http://www.ivoclar.com)).
- Prevent any contact of uncured Variolink Esthetic DC with the skin/mucous membrane and eyes.
- Uncured Variolink Esthetic DC may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

## Disposal information


Remaining stocks must be disposed of according to the corresponding national legal requirements.

## Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Failure of the adhesive bond

## 4 Shelf life and storage

- Storage temperature 2–28 °C / 36–82 °F 
- Store the Variolink Esthetic DC automix syringe with the mixing tip attached after use.
- Do not use the product after the indicated date of expiration.
- Date of expiration: see information on syringes and packages

Before use, visually inspect the packaging and the product for damage. In case of any doubt, consult Ivoclar Vivadent AG or your local trade partner.

## 5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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