# Reinforcement ribbon

## Instructions for use

- DE** Gebrauchsanweisung.
- EN** Instructions for use.
- ES** Instrucciones de uso.
- FR** Mode d'emploi.
- IT** Istruzioni per l'uso.
- PL** Instrukcja postępowania.
- LT** Instruções de uso.

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## INSTRUCTIONS FOR USE

### FLEXIBLE REINFORCEMENT RIBBON FOR DENTAL USE

#### DESCRIPTION

Proclinic reinforcement ribbon is made of long longitudinal white and high crystallized polyethylene fibers, unidirectional that give high mechanical characteristics, ideal for dental use. The special interlacing enables a neat cut and the complete padding with all the most common resin composites, with unfilled adhesives and with acrylic resin.

The product is available in different heights:

1 - 1 tubular - 2 - 3 - 4.

The length of the single roll is available in 50 cm.

#### DESTINATION OF USE

The most important uses of Proclinic reinforcement ribbon are:

- The carrying out of endodontic posts
- Splinting
- Spacers and retentions
- Orthodontic retainer
- Temporary bridge
- Maryland bridge
- Fixing of bridges and prostheses
- Reinforcement of hand-made products in composite or resin
- Temporary stabilization of artificial teeth.

#### INSTRUCTIONS OF USE

Draw the right quantities of ribbon using clean tweezers. Don't touch it with hands.

Proclinic reinforcement ribbon can be cut using sharp scissors, cutting nippers, lancet, knife for wax that must be brought to 200°C and perform a cut that gives a shape as a head of a recorder.

Proclinic reinforcement ribbon is compatible with a greatest number of flowable cement - composite. Usually it is not necessary to bonding.

See the cement manufacturer's instructions for use. When Proclinic reinforcement ribbon is used with acrylic resins it must be wet with an suitable monomer.

The hand-made product of Proclinic reinforcement ribbon must be covered with a resin - composite. If the resin is worn out, we advise to apply a new layer to avoid exposure of fibres. The product has no expiration date, if it is protected from intense light and pollutants.

#### WARNINGS

Proclinic reinforcement ribbon must be covered with a resin - composite.

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