## Anatomic posts

### CARBON FIBER REFILL

#### Instructions for use

- DE** Gebrauchsanweisung.
- EN** Instructions for use.
- ES** Instrucciones de uso.
- FR** Mode d'emploi.
- IT** Istruzioni per l'uso.
- PL** Instrukcja postępowania.
- PT** Instruções de uso.

Last Update: 05/2016  
More information available at:  
[www.proclinc-products.com](http://www.proclinc-products.com)

BIOLOREN S.R.L.  
Via Volta, 59 Saronno  
21047 Italy

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**DE**

### ANATOMISCHE POSTEN

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**EN**

### ANATOMIC POSTS

Once the endodontic treatment is completed, being very careful not to use eugenol-containing cements, proceed to prepare the canal by using the traditional classic instruments (Gates or Peesos). Use of rubber dam is advisable. • Once the post diameter is chosen, proceed to prepare the canal with a calibrated bur, which will allow a perfect fit. It is advisable to drill the canal at least 6-7 mm. The introduction of the post must take place effortlessly and without oscillations. • Once the right height of the post is calculated, the post must be cut with the diamond separating disk. Burs or nose

pliers must never be used. It is advisable to perform the cut outside the oral cavity. • Once these mechanical operations are completed in the canal, then go on to the chemical phase. It is advisable to disinfect the dentin with the combined action of EDTA and sodium hypochlorite, with an adequate irrigation of the canal. Rinse and dry very delicately by using paper cones and/or blowing light air spurts. Do not dehydrate. Handle the post only with proper instruments (clamps), avoiding contact with your hands. In the last case, wash the post with alcohol or a degreasing agent. Pencil in a dentin-enamel adhesive. Use a BIS-GMA resin-based composite-cement due to its affinity for the post resin. The consistency of the paste must be fluid, in order to plug the dentin tubules better. Introduce the cement into the canal through a lentulo spiral. Light-cure the visible part for 20 seconds. Wait for the bonding to complete (approximately 4 minutes). • The prosthesis' core can be made by using self-curing or light-curing compounds.

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