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Bestdent Putty is a silicone elastomer of the vinyl-polysiloxane type whose reticulation takes place based on a polyaddition reaction. The reticulated product's mechanical properties and reproduction of detail are excellent. Bestdent Putty presents a low development, in addition to demonstrating high resistance to inorganic agents and UV rays, as well as good separation properties from organic resins and gypsum.

CLASSIFICATION: Bestdent Putty meets ISO 4823 Type 0 and ADA Specification No. 19 Type II Very High Viscosity.

TECHNICAL INFORMATION:

| | |
|---|-------------------------------|
| Appearance Part A..... | Blue paste. |
| Appearance Part B..... | Yellow paste. |
| Mix ratio..... | 1:1 |
| Mix Time | 30" |
| Working Time at 23°C/73°F | 1'45" (Including mixing time) |
| Permanence in oral cavity..... | 2' |
| Total Cure Time | 2'45" |
| (Time elapsed from the start of the mixture to full reticulation in oral cavity). | |
| Shore A hardness after 24 hr..... | 70 approx. |
| Volumetric mass g/cc..... | 1,55 g/cc approx. |
| Linear contraction after 24 hr..... | 0,2% max. |
| Elastic recovery | 99,1% |
| Compression deformation | 3% approx. |

USE: Bestdent Putty Hard is used as a tray material in the one-step technique and Two-step print technique together with the Polyaddition Light line.

MIXTURE: Use Part A and Part B in equal portions by using the two appropriate measures. Mix the two parts with your hands until achieving a homogeneous color (30 sec. max.).

TAKING THE IMPRESSION: The mixed material is placed on a tray or directly on the teeth. Next introduce the tray with the material into the mouth before 1', counted from the start of the mixture. Keep it in that position for at least 2'. The time periods mentioned are valid for a temperature of +23°C/+73°F and 50% Relative Humidity. At higher temperatures time is shortened and vice versa. If the first print is placed separately for more than 30' before taking the second print or performing a correction, it is preferable to apply a very fine layer of adhesive. The impression may be cast 30' after removal from the oral cavity. Avoid extreme storage conditions such as proximity to heaters or exposure to sunlight.

WARNING: The use of latex gloves may disturb the addition reaction, so the use of polyethylene gloves is advisable. Containers must be closed immediately after use. Do not exchange container caps. This prevents the products kept in different containers from contacting each other. Silicone adheres to tissue and it is difficult to remove even with chemical agents.

STORAGE: For the product to keep its optimum characteristics, it is advisable to store it at 23°C/73°F, and use it within a period of 15 months.

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[EN] Bestdent Light Body Fast Set is a silicone elastomer which reticulates as a result of a polyaddition reaction. The hydrophilic properties of Bestdent Light Body Fast Set enable the user to obtain a higher degree of accuracy and a lower risk of bubble formation due to the presence of saliva. The mechanical properties and reproduction of detail of the reticulated product are excellent. Bestdent Light Body Fast Set has shown high resistance to inorganic agents and UV rays, as well as good separation properties from organic resins and gypsum.

CLASSIFICATION: Bestdent Light meets ISO 4823 Type 0 and ADA Specification No. 19 Type II Very Low Viscosity.

TECHNICAL INFORMATION:

| | |
|--------------------------------|----------------|
| Appearance Part A..... | White fluid. |
| Appearance Part B..... | Colored fluid. |
| Mix Ratio..... | 1:1 |
| Working Time at 23°C/73°F..... | 1'30" |
| Permanence in Oral Cavity..... | 2' |
| Total Cure Time..... | 2' |

(Time elapsed from the start of the mixture to full reticulation in oral cavity).

| | |
|-------------------------------------|-------------------|
| Shore A Hardness after 24 hr..... | 53 approx. |
| Volumetric Mass g/cc..... | 1,37 g/cc approx. |
| Linear Contraction after 24 hr..... | 0,2% max. |
| Elastic Recovery..... | 99,6% |
| Compression Deformation..... | 3% aprox. |

USE: Bestdent Light Body Fast Set is used as a wash material in the one-step technique and the Two-step technique together with the Polyaddition Putty line.

MIXTURE: Loading the gun: • Keep the blocking stem in the upward position. • Insert the plunger into the gun while making sure that blocking stem is pointing upward. • Push the plunger all the way to the bottom and quickly release the blocking stem. • Place the cartridge on the support. • Pressure the control lever so the plunger will touch the cartridge. • Remove the cartridge hood by giving it one-quarter turn.

USE: • Pressure the control lever to mix the two parts. • Remove the first millimeter of material. • To stop the flow, release the control lever. • Leave the mixer in place until the next use. • Replace the mixer before reusing the material. • Remove all the hardened material residues before replacing the mixer. • To remove the empty cartridge, keep the blocking stem in an upward position and pull the plunger out to release the cartridge. • Before introducing the new mixer, insert the mixing drum. • The gun may be sterilized in an aseptic solution.

TAKING THE IMPRESSION: • The mixed material is placed on a tray or directly on the teeth. • Next introduce the tray with the material into the mouth before 30", counted from the start of the mixture. • Keep it in that position for at least 2'30". • The time periods mentioned are valid for a temperature of +23°C/+73°F and 50% Relative Humidity. • At higher temperatures time is shortened and vice versa.

STORAGE: For the product to keep its optimum characteristics, it is advisable to store it at 23°C/73°F, and use it within a period of 18 months.

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