SELF-ETCH ADHESIVE

DE

SELF-ETCH ADHESIVE

EN

SELF-ETCH ADHESIVE

ES

SELF-ETCH ADHESIVE

FR

SELF-ETCH ADHESIVE

IT

PRODUKTBESCHREIBUNG

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PRODUCT DESCRIPTION

Self-etch dental adhesive is a simple to use self-etching light curing one component adhesive. It ensures a strong bonding of composites to enamel and dentin.

Self-etch dental adhesive is compatible with all current brands of visible light cure composite restorative materials.

INDICATIONS/INTENDED USE

- Adhesive for:
- Restorations with light cure composites

PERFORMANCE FEATURES

The performance features of the product meet the requirements of the intended use.

CONTRAINDICATIONS

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

PATIENT TARGET GROUP

Persons who are treated during a dental procedure.

INTENDED USERS

This medical device should only be used by a professionally trained dental practitioner.

INCOMPATIBILITY WITH OTHER MATERIALS

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of Self-etch dental adhesive. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

APPLICATION

1. Isolation

Rubber dam is the recommended method of isolation.

2. Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5 - 1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

3. Pulp Protection

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

4. Application of Self-etch dental adhesive

Apply **Self-etch dental adhesive** with a brush onto the slightly moist enamel and dentin surfaces for at least 15 seconds with agitation.

The material should build a homogeneous layer. Air thin gently for 5-10 seconds to remove the volatile components and to disperse the adhesive. Then light cure with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2 for 20 seconds before placement of a composite.

5. Restorative Placement

Apply the restorative material according to the instructions of the manufacturer. Best results are obtained with application of a thin layer of a light cure flowable composite followed by the application of a moldable composite. Light cure each composite layer separately according to the corresponding user instructions.

STORAGE

Do not store above 25 °C (77 °F)! Protect from direct sunlight. Do not use after expiry date.

ADDITIONAL NOTES/WARNINGS

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

COMPOSITION

Hydrophilic methacrylates, MDP, ethanol, photo initiators.

DISPOSAL

Disposal of the product according to local authority regulations.

REPORTING OBLIGATION

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

NOTE

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

WARRANTY

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior Dental Materials GmbH product.

LIMITATION OF LIABILITY

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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