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

MTA PROCLINIC EXPERT		
REFERENCES	78551: white colour	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class IIb, Directive 93/42/EEC. Notified Body: Kiwa Cermet Italia S.p.A., No. 0476	
DESCRIPTION	Bioceramic endodontic cement composed of several mineral oxides. It is constituted by thin hydrophilic particles. Non sterile product.	
COMPOSITION	Tricalcium silicate, Dicalcium silicate, Tricalcium aluminate, Tetracalcium aluminoferrite, Calcium oxide, Bismuth oxide. Distilled water.	
PACKAGING	1 unit MTA, 1 unit distilled water, 1 measuring spoon	
SHELF LIFE	3 years	
INTENDED USE	Indicated in cases of root canal lateral and furcation perforations, internal resorption, reverse root filling, pulp capping, pulpotomy, apexification, and apexogenesis.	
INSTRUCTIONS FOR USE	<ol style="list-style-type: none"> 1. Sterilize a glass slab, a metal spatula and all instruments for the insertion of Proclinic Expert MTA. 2. Mix for 30 seconds the content of 1 sachet of Proclinic Expert MTA (or 1 spoon of Proclinic Expert MTA) with 1 drop of distilled water. The mixture should be homogeneous and with a consistency similar to wet sand; 3. Place the cement on the selected site with a sterilized amalgam carrier or other appropriate instrument; 4. Condense the cement with instruments such as amalgam condensers, a number 1 spatula or absorbent paper points moistened with distilled water. <p><i>(*) Important: If Proclinic Expert MTA is not used immediately after mixing, its dehydration can be prevented and the working time increased by covering the mix on the glass slab with a wet gauze. If the mixed MTA dehydrates, it must be discarded.</i></p>	

<p>WARNINGS & PRECAUTIONS</p>	<ul style="list-style-type: none"> • Use eye protection, mask and gloves when handling Proclinic Expert MTA. In case of eye or skin contact, wash abundantly with water; • Only open the sachet or bottle immediately prior to use. The Proclinic Expert MTA powder is very sensitive to humidity. Close the bottle tightly after each use; • Do not use Proclinic Expert MTA to fill a root canal. Its viscosity is inadequate for this procedure and very difficult to remove in case of a reentry; • Do not use Proclinic Expert MTA in areas of the tooth in contact with the gingival sulcus or it will be completely dissolved; • Do not use MTA in areas with exposure to light in order to avoid darkening of the dental structure; • Only use Proclinic Expert MTA after remission of acute signs and symptoms of the endodontic disease. The acidic pH of endodontically compromised sites (lesions) prevents its setting reaction; • Apply Proclinic Expert MTA carefully. Similar to other endodontic cements, it is resorbed if extruded. However, the excess of any cement may impair the healing process; • Do not use MTA on patients with reported sensitivity to any of its components. <p>For exclusive use of the dentist.</p>
<p>HANDLING AND STORAGE CONDITIONS</p>	<p>Store in a cool dry place, away from direct sunlight. Protect from moisture and excessive heat.</p>