**FICHA TÉCNICA**

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## TECHNICAL DATA SHEET

<b>TEMPORARY CEMENT PROCLINIC EXPERT</b>		
<b>REFERENCES</b>	<b>78576:</b> Temporary Cement Implant.	
<b>MANUFACTURER</b>	<b>B.J.M. Laboratories Ltd.</b> 12 Hassadna St. Industrial Park, Or-Yehuda 6022011 - ISRAEL Tel. 972-74-7000111 Fax. 972-3-7353020	
<b>CLASSIFICATION OF THE PRODUCT</b>	Medical Device, Class IIa, Directive 93/42/EEC.	
<b>DESCRIPTION</b>	Temporary cement for implant retained crowns	
<b>COMPOSITION</b>	Multifunctional methacrylates, Urethane dimethacrylate, Catalysts	
<b>PACKAGING</b>	1 automix syringe 6 g 5 mixing tips	
<b>SHELF LIFE</b>	2 years from the date of manufacture	
<b>INTENDED USE</b>	Specially formulated cement for luting of provisional crowns and bridges on implant abutments.	
<b>INSTRUCTIONS FOR USE</b>	<ol style="list-style-type: none"> <li>1) Finish and polish the restoration. Check occlusion and proximal contacts and complete final adjustments.</li> <li>2) Dry and isolate the area being treated. Absolute dryness of the restoration is not necessary, but contamination should be avoided.</li> <li>3) Make sure that the restoration is clean and dry. It is advisable to apply a thin layer of water-based lubricant to the abutment surface to facilitate future removal of implant restorations.</li> <li>4) Bleed the syringe before attaching the automix tip. Bleed the syringe again after placing Proclinic Expert Temporary Cement's automix tip to ensure correct flow of base and catalyst components. Place a thin layer of cement in the restoration. Remember that the working time for this cement is between 45-60 seconds.</li> <li>5) Keep the crown in the desired position for 2-2.5 minutes by asking the patient to close on a cotton ball to ensure accurate, secure placement. Excess cement can safely be removed after 2.5 minutes.</li> </ol>	

	6) Total setting time of the cement is between 4 - 5 minutes, depending on conditions.
--	--

<b>SAFETY INFORMATION</b>	Though uncommon, some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use, remove from the patient's mouth and refer patient to a doctor.
---------------------------	---

<b>HANDLING AND STORAGE CONDITIONS</b>	Store below 24°C (68°F). Refrigerate for maximum shelf life. It is not recommended for use beyond the expiration date marked on the packaging.
--	--