LuxaCore

Z-Dual

Inhalt | Contents

Gebrauchsinformation – DEUTSCH	4
Instructions for use – ENGLISH	14
Mode d'emploi – FRANÇAIS	24
Istruzioni d'uso – ITALIANO	34
Instrucciones de empleo – ESPAÑOL	44
Instruções de uso – PORTUGUÊS	54
Gebruiksaanwijzing – NEDERLANDS	64
Brugsanvisning – DANSK	74
Användarinstruktioner – SVENSKA	84
Instrukcja użycia – POLSKI	94
Информация по применению – РУССКИЙ	105
使用说明-简体中文	116
사용 설명서-독일어	125
使用説明書-日本語	134

Produktbeschreibung

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Product Description

LuxaCore Z-Dual is an automatic mixing, dual-hardening composite that has been specially developed for all types of stump build-ups and construction fillings. The curing time can be self-determined due to the option of light curing. It is possible to apply LuxaCore Z-Dual directly using the Intraoral tip and the Endo tip. LuxaCore Z-Dual is radiopaque.

Indications

- All types of stump build-up
- Cementation of root posts

Contraindications

- Do not use the material in cases of allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.
- Do not use the material with a single-bottle bond because an optimum adhesive bond is generally not assured with chemical or dual-hardening materials.
- Do not use the material if a dry working environment is not possible.

Basic safety instructions

- Keep out of reach of children!
- For dental use only!
- In order to avoid any pulp reactions, use suitable pulp/dentin protection!
- Avoid contact with skin and eyes. In case of accidental contact, irrigate immediately with

plenty of water and consult a physician, if necessary!

Side effects

There have been no systemic side effects to date.

Interactions

- Materials containing eugenol may inhibit the polymerization of LuxaCore Z-Dual and lead to discoloration.
- Air containing water and oil can inhibit the polymerization of LuxaCore Z-Dual at the contact surfaces.

Application system

- Automix cartridge: See “Product Description” on page 14.
- Smartmix injector: See “Using the Smartmix Injector” on page 20.

Timing

Working time*	≈ 1:30 min
Chemical curing	intraoral 5:00 min
Light curing	20 s (layer ≤ 2 mm) 40 s (layer ≤ 4 mm)

***Please note:** The times apply at a room temperature of 23 °C and a normal relative air humidity of 50%. Higher temperatures shorten, lower temperatures increase this time.

Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.

Recommended use

Carry out stump build-up

1. Isolate the tooth prior to preparation using a rubber dam.
2. Remove existing fillings and caries.
3. If the bonding agent used requires an additional etching step then etch using the etching gel per the manufacturer's instructions.

Please note: DMG recommends the use of the dual-hardening adhesive system LuxaBond-Total Etch. In this case it is essential to etch the affected

enamel and dentin surfaces. Observe the manufacturer's instructions.

Caution! Unprotected pulp.

- ▶ Prior to bonding in deep cavities in areas close to pulp, provide protection by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement for example.
4. Apply bonding agent in accordance with the manufacturer's instructions.
 5. In order to simplify the application of LuxaCore Z-Dual, place a matrix around the prepared tooth and apply LuxaCore Z-Dual directly into the cavity. If a pre-prepared stump build-up form is used, apply LuxaCore Z-Dual directly into the stump build-up form and place on the prepared tooth.

Please note: LuxaCore Z-Dual can be formed within the working time of 1:30 minutes using a standard composite instrument, e.g. a Heide-mann spatula. In order to enhance processing,

coat the Heidemann spatula with a bond as a separating medium.

6. Leave LuxaCore Z-Dual to bond intraorally for approx. 5 minutes or expose to a suitable light curing unit for a minimum 40 seconds (maximum layer thickness ≤ 4 mm). In the case of thinner layers (≤ 2 mm) exposure to light for 20 seconds is sufficient.

Please note: With chemical curing, low temperatures (e.g. when working in the area of the front teeth) lead to longer bonding times.

7. If used, remove a matrix only once the material has hardened fully.
8. Produce the preparation for the desired restoration form using standard preparation instruments.
9. If a plastic temporary restoration is to be applied to the prepared tooth, isolate the preparation using a suitable separation medium (e.g. Vaseline).

Cementation of root posts

1. Carry out any preparation procedures necessary for root posts in accordance with the manufacturer's instructions.
2. Prepare the root canal in accordance with the requirements of the respective root post.
3. If the bonding agent used requires an additional etching step then etch inside the canal using the etching gel per the manufacturer's instructions.

Please note: It is advisable to use the dual-hardening adhesive system LuxaBond-Total Etch. In this case it is essential to etch the affected enamel and dentin surfaces. Observe the manufacturer's instructions.

4. Apply bonding agent in accordance with the manufacturer's instructions. If an adhesive system with a light-curing component is used, use a paper tip to remove any excess from the canal prior to the application of light. Otherwise the lumen may be blocked.


5. Apply LuxaCore Z-Dual to the canal lumen. When doing so keep the end of the tip constantly immersed in the material. If desired, it is additionally possible to coat the root post with LuxaCore Z-Dual. The LuxaCore Z-Dual working time is 1:30 minutes.

Please note: When using LuxaBond-Total Etch, the working time of LuxaCore Z-Dual in the root canal is reduced to approx. 30 seconds due to the contact with the pre-bond. It is also possible that the working time may be reduced with other bonding agents. Observe the manufacturer's instructions.

6. Set the root post.
7. Using a suitable lamp, apply light to the material for 20 seconds.
8. Leave the material to harden for approx. 5:00 minutes. This guarantees complete curing, also in the areas that cannot be reached by light.

Please note: With LuxaCore Z-Dual it is possible to work on the coronal portion of the tooth as a stump build-up. Due to the easy-flowing consistency, it is recommended that work be carried out using a matrix (see “Carry out stump build-up” on page 16).

Using the Automix cartridge

Insert cartridge

1. Push up the lever [A] at the back of the Automix dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C], insert the cartridge and lock it using the plastic lever [C].

Please note: Make sure that the notches on the cartridge line up with those on the Automix dispenser.

Attaching a mixing cannula

1. Remove the cartridge cap [D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Please note: In order to achieve an optimum mixing result, DMG recommends the use of mixing tips available from DMG. All mixing tips are also available as refill packs.

Please note: Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it.
4. Attach an Intraoral tip [F] or Endo tip.

Applying the material

- ▶ When using a cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- ▶ The material is mixed in the mixing tip by pulling the trigger [G] of the Automix dispenser and it is then applied directly.


Please note: After the application, leave the used mixing tip on the cartridge to cap unused material!

Remove cartridge


1. Push up the lever [A] at the back of the Automix dispenser and pull back the slide [B] completely.

2. Lift up the plastic lever [C] and remove the cartridge.

Using the Smartmix Injector


1. Before putting on the mixing cannula, remove the cap [A] or used mixing tip after twisting it 90° counterclockwise, and then discard it.


2. Attach a new mixing tip [B].

Please note: In order to achieve an optimum mixing result, DMG recommends the use of mixing tips available from DMG. All mixing tips are also available as refill packs.

Please note: Make sure that the notches on the mixing tip line up with those on the Smartmix injector [C].

3. Turn the mixing tip clockwise 90° to lock it.

4. Attach an Intraoral tip or Endo tip.
5. The material is mixed in the cannula when it is expressed and it can be applied directly.

Please note: When using an injector for the first time, squeeze out a pea-sized amount of material and discard it!

Leave the used mixing cannula on the Smartmix injector after use. It serves as a cap!

Disinfection of the Automix dispenser

The Automix dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix dispenser is also autoclavable.

Technical specifications

Compression strength	380 MPa
Flexural strength	150 MPa
Diametral tensile strength	63 MPa
Water absorption	20 $\mu\text{g}/\text{mm}^3$

Ingredients

Barium glass, pyrogenic silicic acid, nanofiller and zirconium dioxide in a Bis-GMA based matrix from dental resins. Filler content: 71 % by weight = 50 % by vol. (0.02 – 2.4 μm)

Storage and shelf life

- Store in a dry place at 2 - 25 °C/36 - 77 °F!
- If possible, store unopened material in a refrigerator.
- Do not use after the expiration date.

Classification

ISO 4049:2000, type 2, class 3

Available product sizes

Refill packs

1 cartridge @ 48 g paste, 30 Automix tips, 20 Intraoral tips, 10 Endo tips	
Color A3	REF 213330
Color blue	REF 213331
Color light opaque	REF 213332
2 tips @ 9 g paste, 20 Smartmix tips, 10 Intraoral tips, 10 Endo tips*	
Color A3	REF 213334
Color blue	REF 213335

2 tips @ 9 g paste, 20 Smartmix tips 10 Intraoral tips, 10 Endo tips*	
Color light opaque	REF 213336

Accessories

1 Automix dispenser type 25 1:1	REF 110253
50 Automix tips yellow 1:1	REF 909200
50 Intraoral tips yellow 1:1	REF 909201
50 Endo tips yellow	REF 212280
50 Smartmix tips combi	REF 212042
50 Smartmix tips Endo	REF 212041

* In the USA and Canada, 10 Smartmix tips and 10 Intraoral tips are included.

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DMG Chemisch-Pharmazeutische Fabrik GmbH