Prior to use, carefully read the instructions for use.

GB

GC RELINE BONDING AGENT

GC RELINE BONDING AGENT is a specially designed methacrylate base bonding agent that enhances the bond between GC Reline Material and denture base materials.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS:

For use as a bonding agent when processing GC RELINE Hard Denture Chairside Reline Material.

CONTRAINDICATIONS:

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE:

Note: Consult the instructions for use for the GC RELINE Hard Denture Chairside Reline Material for complete processing information.

- Roughen the area of the denture to be relined. Clean and dry the denture thoroughly.
- Coat the artificial teeth, labial and buccal surfaces of the denture with COE LUBRICANT or Vaseline using a cotton swab being careful to stay at least 3 mm (1/8 inch) from the peripheral border. *Do not contaminate the tissue side of the denture with lubricant.*
- Using the brush provided, coat the area to be relined with a uniform coating of the GC RELINE BONDING AGENT. Air-dry the bonding agent for approximately 10 seconds to allow the solvent to evaporate.

Note: GC RELINE BONDING AGENT is highly volatile; therefore be sure to replace and tightly close the cap after each use.
- Continue processing the GC RELINE Hard Denture Chairside Reline Material following the recommended procedures.

CAUTIONS:

- The solvent in GC RELINE BONDING AGENT is highly volatile and flammable. Avoid sources of ignition and keep tight capped when not in use.
- In case of contact with skin, wash the affected area thoroughly.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Wash the brush with GC RELINE liquid or ethanol after using the GC RELINE BONDING AGENT.

STORAGE:

Store at normal temperature 15-25 °C (60 - 77 °F)
(Shelf life: 3 years from date of manufacture)

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Alsip, Illinois 60803 U.S.A.
PRINTED IN U.S.A.