





ProphyMax - ProphyMax Lux


User's manual • Notice d'utilisation • Manual de uso
Benutzerhandbuch • Manuale d'uso • Handleiding voor de gebruiker

Illustrations

Illustrations

1. Introduction

1. Introduction	1
2. Installation	5
3. Starting and stopping the device	6
4. Using ultrasound mode.....	7
5. Using polishing mode	8
6. Maintenance	9
7. Troubleshooting	10
8. Technical specifications	13

The Prophy Max™ is a multi-clinical piezoelectric ultrasound generator combined with an air polisher. It enables treatment to be carried out in optimal sterile conditions.

Prophy Max™ is equipped with a Suprasson® Newtron handpiece (ultrasound) and a ProphyPen handpiece (polishing).

Prophy Max™ Lux is equipped with a Suprasson® Lux light handpiece (ultrasound) and a ProphyPen handpiece (polishing).

General characteristics

- Four ultrasound power ranges, each adjustable with 10 levels.
- Color display for visual and ergonomic representation of the power.
- Boost function for instant selection of maximum power in the selected mode.
- Memory for each of the four power modes for easy recall of a chosen power level.
- Polishing function via the powder drum.
- Irrigation via the water supply with flow regulation on the handpiece connector.
- Main functions on the multifunction footswitch.
- Integrated heating system for the spray for optimal patient comfort.

The Prophy Max is a multifunction device for scaling, endodontic and periodontic treatment in sterile conditions, as well as for polishing.

Irrigation kit

Available as an option

- Irrigation kit comprising a peristaltic pump, an irrigation line and two reservoirs for irrigation liquids.
- Independent adjustment of irrigation flow.
- Syringe effect.
- Easy manual selection of the active reservoir or complete shut-off.
- Support for bottles or pouches of sterile solution.

Parts included

- A Prophy Max or Prophy Max Lux device (depending on the option)
- A mains cord
- A single or multifunction pedal footswitch (depending on the option)
- Suprasson® Newtron or Suprasson® Lux handpiece depending on the version
- ProphyPen handpiece with titanium nozzle
- Tip kit depending on the sales version
- ProphyPen kit depending on the version
- Powder Max introductory powder kit
- Pack three plastic nozzle

General view

Refer to figure 1 on the fold-out part of the cover.

- | | |
|---------------------------------|---|
| 1. Rear face | 5. Irrigation kit and removable cover |
| 2. Powder drum | 6. Button to control pump flow |
| 3. Suprasson® Newtron handpiece | 7. Display area and sensitive control panel |
| 4. Multifunction footswitch | 8. ProphyPen handpiece |

1. Introduction

Front panel - Display

9. Boost indicator **BOOST**

Lights when the BOOST function is activated.

10. Standby indicator **⏻**

Lights when the device is in standby mode (energy saving).

11. Display

The display background color indicates:

- **power on**, the color indicator appears bluish white. Any other color signals a malfunction (see chapter troubleshooting).
- **in use**, the color displayed (green, yellow, blue and orange) corresponds to the color code of the ultrasound power range selected.

The figures displayed indicate:

- **power on**, the number 88 indicates the correct operation of the segments of the display. Any other value signals a malfunction of the indicator (see chapter troubleshooting).
- **in use**, the ultrasound power level (1 to 10) for the range selected, controlled by the **- +** keys (Ref. 22).
- When programming:
 - the length of time of inactivity before the device enters standby mode;
 - for Prophy Max™ Lux only, the length of time before the light is switched off on the Suprasson® Lux handpiece, when ultrasound is no longer in use;

the display blanks out (black background and no figures) when the device is in air polisher mode.

The following figure shows an example of the display, from left to right, power on, in use (setting the ultrasound power level) and when programming (light time limit).

12. Water supply indicator **~**

Lights when the **~** key is pressed (ref. 19), and indicates irrigation via the drinking water supply.

13. Irrigation kit indicator **⊂**

Lights when the **⊂** key is pressed (ref. 20), and indicates irrigation via one of the two reservoirs.

14. Spray or dry mode indicator **⚡**

Lights when the **⚡** key is pressed (ref. 21).

15. Air polisher mode indicator **↘**

Indicates activation of air polishing mode when the **↘** key is pressed (ref. 23).

Sensitive controls panel

16. Ultrasound power ranges keys

Each yellow, green, blue and orange key corresponds to an ultrasound power range, from the weakest to the strongest as follows:

- Green background: low power (ref. 16a).
 - Yellow background: intermediate power (ref. 16b).
 - Blue background: high power (ref. 16c).
 - Orange background: very high power (ref. 16d).
- The range selected:
- a reminder is given by the color of the display backlighting;
 - is programmed in the factory at level 5 except for the orange key which is programmed at level 1 and uses irrigation from the water supply.

The range may also be selected with the multifunction footswitch (see chapter on the multifunction footswitch pedal).

17. Boost **BOOST**

Applies the power (value 10) of the selected range; the **BOOST** lights up. A second press on the key exits the mode. This mode can also be accessed from the multifunction footswitch (cannot be memorized and use must be occasional).

18. Memory key **⏻**

Memorizes a power level for each of the four power ranges, with the exception of switching on irrigation and Boost.

19. Water supply key **~**

Pressing this key selects drinking water (water supply) as the irrigation solution. The indicator lights up. Irrigation will only take place after the spray or dry **⚡** key (ref. 21) is pressed; the rate of flow may be set via the handpiece connector. The two keys 19 and 20 are both exclusive and may be memorized.

20. Irrigation kit key **⊂**

Pressing this key selects the irrigation kit (optional). Irrigation will only take place after the spray or dry **⚡** key (ref. 21) is pressed.

The rate of flow may be adjusted via the control knob on the pump. The two keys 19 and 20 are both exclusive and may be memorized.

21. Spray or dry key **⚡**

Switches on irrigation after keys 19 or 20 have been selected. The indicator (ref. 14) lights up. A second press turns irrigation off.

22. The **- +** keys

Reduce (**-** key) or increase (**+** key) ultrasound power in the selected range. The values go from 1 to 10.

1. Introduction

23. Polishing mode key

Activates polishing mode. The indicator lights up and the screen goes out. This mode can also be activated by the right pedal on the footswitch.

24. Purge key

A short press purges the handpiece powder circuit for a predetermined length of time of a few seconds.

25. Button to enter or leave standby mode

This button is mainly used for leaving standby mode. The initialization sequence is started. *Only this button can be used to leave standby mode.* This button can also be used to put the device into standby mode, but the device goes into standby automatically after a set time (see § *Setting the standby delay of the device*).

Rear panel

26. Team Up connector (optional)

Connection to the electroanalgesia device used during prophylactic treatment.

27. Compressed air supply filter

Before use the compressed air first passes through a filter (ref. 2a) with a purge button (ref. 2b). A rapid disconnection system with an automatic valve enables it to be connected to the mains.

28. Footswitch connector

Connection for the single or multifunction footswitch pedal.

29. Irrigation kit connector

Refer to kit instructions for details.

30. On/off switch (O/I)

31. Air pressure adjustment

Air pressure is set in the factory. For any adjustment of the pressure, contact an approved installer of dental equipment.

32. Drinking water supply connector

Its rapid disconnection system with an automatic valve enables it to be connected to the supply. A removable filter is located on the water inlet pipe.

33. Purge of compressed air filter

Press this purge button to eliminate air remaining in the filter

34. Fan

Do not obstruct the fan at the risk of causing severe damage to the device.

35. Power socket

Connection to mains electricity is by a standard connector conforming to IEC 60-320.

36. Fuse housing

Contains the device's fuses.

Suprasson® and Suprasson® Lux handpieces

Refer to the relevant instructions.

ProphyPen handpiece

Refer to the relevant instructions.

Multifunction footswitch

37. Left pedal

Activates irrigation without ultrasound. When the irrigation kit is available, pressing the pedal initiates irrigation (syringe effect - with a flow rate of about 120 ml/min).

38. Connector

Connection of the device.

39. Range of functions

Each press selects one of the four power modes. The ultrasound power setting used will be the one memorized, or by default the factory setting, (the mid position value of 5, except for the top range which is set to 1). Check the display to view the function selected.

40. Boost

See notes relating to the key 17.

41. Irrigation

See notes relating to the key 21.

42. Right pedal

- Selected ultrasound function: activates ultrasound and associated functions that have been previously selected.
- Selected polishing function: activates polishing.

1. Introduction

Irrigation kit

The pump (ref. 48) sends the contents of the selected reservoir to the Suprasson® handpiece fitted with the Sterijoint (54).

It is possible to regulate the flow by turning the limiter located above the pump motor (ref. 46).

The selector (ref. 47) is used as follows:

- selector to the left : left reservoir.
- selector to the right: right reservoir.
- selector in the center: both reservoirs closed.

The two reservoirs (ref. 44) are independent and fixed on either side of the pump. They can be sterilized in an autoclave.

The connector (ref. 43) links the irrigation kit to the device.

2. Installation

Recommendations

For your safety, connections (water, air and electricity) should be carried out by an approved installer of dental equipment. The device requires the following connections (see technical specifications chapter).

Connect to the drinking water supply using an automatic cutoff valve of the *Aquastop* type, or to the dental chair water supply.

Caution: always turn off the water supply when you leave the dental practice.

Setting up the equipment

- Position the device in the place where it is going to be used.
- Connect the footswitch (ref. 38) to the connector on the device(ref. 28).
- Connect the optional equipment (*Team Up* ref. 26, irrigation kit ref. 29).
- Connect the water supply pipe with its filter to the connector located at the rear of the device (ref. 32).
- Connect the air supply pipe to the air filter located at the rear of the device (ref. 27).
- Plug in the power socket (ref. 35)
- Open the supply taps (water-air).


3. Starting and stopping the device

This device should only be used by professionals (dental practitioners and hygienists). Before being used for the first time it is:

- imperative to clean, disinfect using alcohol or disinfectant dental wipes and/or sterilize the accessories (casing of the device, handpiece, tips, wrenches, etc.);
- recommended that the accessories be plugged in and the spray operated for a few minutes to rinse the device's internal water circuits.

Important: It is imperative that the water supply is turned off when leaving the dental practice.

Switching on the power

1. **Set the power switch to the on position (I) (rep. 30).**
The device automatically goes into standby mode: the green indicator  (ref. 10) on the front lights up. The device is in energy saving mode.

To operate

1. **Press the  (ref. 25) key on the front of the device**
The initialization phase lasting a few seconds is launched. The information displayed is as follows (any other values indicate malfunction):
 - Display: the figure 88 is displayed.
 - Background color of the display: bluish-white.
 - Indicators: they all light up.At the end of the initialization phase, the display will show one of the following:
 - by default, low power mode is at level 5 on a green background, no irrigation, water supply selected, light out in three seconds, (Prophy Max Lux), enters standby mode automatically after nine minutes;
 - or the last configuration used.The device is ready to use.

Entering standby mode

1. **Press the  (rep 25) key**
 - The front face goes out.
 - The device is in energy saving mode.
 - The  (rep 10) indicator lights up.Reminder: the device enters standby mode on its own after a predetermined time. See § *Setting the device's standby delay*.

Leaving standby mode

1. **Press the  key (rep 25)**
The  (rep 10) indicator goes out and the device returns to the configuration being used before entering standby.

Changing the delays

Delays refer to the time the light in the Suprasson® Lux handpiece stays on after ultrasound has been used, and the time taken before the device goes into standby mode after a period of inactivity.

Delay before the light goes out


This function only applies to the Prophy Max Lux. It specifies the length of time before the light goes out after ultrasound has been used (three seconds by default). To change this setting:

1. **Put the device into standby ( key).**
2. **Press the yellow and boost keys together, then the  key while keeping the other keys depressed.**
The display will indicate «L» (Light) followed by a figure giving the current time delay in seconds.
3. **Release the keys.**
4. **Set the time with the + or - keys.**
The length of time can be set from 0 (L-) to nine seconds (L9). The new value comes into effect when the device is restarted.

Therefore a light out time delay of 7 seconds will be displayed at "L7".

Delay before the device enters standby mode

This setting specifies the length of time before the device enters standby mode. To change this setting:

1. **Put the device into standby ( key).**
2. **Press the yellow and blue keys together, then the  key, while keeping the other keys pressed.**
The display will indicate "A" (Stop) followed by a figure giving the current time delay in minutes.
3. **Release the keys.**
4. **Set the time with the + or - keys.**
The length of time can be set from no standby to 15 minutes. The new value comes into effect when the device is restarted.

"AC" indicates that standby mode will be activated after 12 minutes. The table below shows the value displayed according to the desired time delay in minutes.

Display	min	Display	min
A-	0 (*)	A8	8
A1	1	A9	9
A2	2	AA	10
A3	3	AB	11
A4	4	AC	12
A5	5	AD	13
A6	6	AE	14
A7	7	AF	15

(*) 0 = no standby

Switching off the device

1. **Set the on/off switch to the 0 position.**

4. Using ultrasound mode

Choosing the ultrasound power range

1. Select a power range by pressing one of the four keys (green, yellow, blue, orange) (ref. 16) on the keyboard.
By default the display indicates the color of the range selected at power setting 5 (1 for high power) and irrigation from the water supply.
Reminder: the ultrasound power range may be selected with the multifunction footswitch.

Setting ultrasound power

Incremental setting

1. Press the **+** key to increase ultrasound power or the **-** key to decrease it from 1 to 10, by steps of 1.

Setting the power to maximum (BOOST)

1. Press the **BOOST** key on the front of the unit (ref. 17) or the multifunction footswitch (ref. 40) to set the power to maximum.
The indicator lights up (ref. 9) and the display indicates the value of 10 in the selected range.
To exit *Boost* mode, press the key on the front of the unit or press the multifunction footswitch once more, or change the power range.

Spray or dry

1. Press the **⚡** key on the face of the unit. (ref. 21) or the multifunction footswitch (ref. 41) to activate irrigation.
 - In water supply mode, controlling the water flow is via the connector ring on the Suprasson[®] handpiece (ref. 55)
 - In irrigation kit mode, controlling the water flow is via the flow limiter located on the pump (ref. 46).

Memorizing the ultrasound power level

1. Select the desired power range.
2. Change the power level with the **- +** keys (ref. 22)
3. Press the memory key **⏏**. (ref. 18)
The power display shows "--" for one second to confirm memorization.

Recalling a level from memory

1. Press the key for the desired power range (ref. 16)
The memorized value is displayed.
Refer to the Tipbook for power settings of Satelec tips.

5. Using polishing mode

Filling the powder drum

1. Put the device into standby mode to depressurize the drum.
2. Unscrew the cap (ref. 50).
3. Fill the drum with *Satelec PowderMax* powder up to the level marker (ref. 51).
4. Press the  key (ref. 23). The display blanks out and the indicator (ref. 15) lights up.

Testing the effectiveness of the powder jet

1. Point the handpiece at the test card and press the pedal.
2. Set the water flow on the handpiece connector.
 - If there is not enough water: powder accumulates on the card.
 - If there is too much water: the polishing effect is less effective.
 - If the setting is correct: the powder is immediately removed by the water.
 - Where the powder jet is not powerful enough, clean the nozzle (see chapter "Cleaning the polishing system") or have the air pressure adjusted by an approved installer of dental equipment.

Polishing

1. Use the PropyPen polishing handpiece.
2. Press the pedal to start polishing. It is important to observe a jet/surface angle of from 45° to 60° and to keep to a spray/tooth distance of 4 to 6 mm.

Polishing mode can also be selected with the multifunction footswitch (ref. 39).

Purging the polishing circuit

1. Release the pedal to stop polishing.
2. Point the PropyPen handpiece towards the cuspidor.
3. Press the  key (ref. 24) to purge the circuit.

6. Maintenance

Cleaning and disinfecting the device

The casing of the Prophy Max as well as the cords should be cleaned and disinfected every day, using alcohol or disinfectant dental wipes.

Cleaning, disinfection and sterilization of ultrasound systems

In the rear of the handpiece an irrigation circuit waterproof o'ring (1.15x1) is accessible and visible thanks to a groove of the rear socket. This allows for easy replacement of the o'ring in the case of damage or wearing. You can replace it yourself or have it done. After a number of disconnections/connections and/or sterilizations, it is recommended to lubricate this o'ring with dental instrument quality lubricant to maintain its efficiency and avoid internal leakage.

Without an irrigation kit

1. Separate the Suprasson® handpiece (ref. 53) from the cord (ref. 55) by pulling, without turning.
2. Unscrew the tip (ref. 52) from the Suprasson® handpiece.
3. Clean and sterilize the Suprasson® handpiece and tip after each use.

With an irrigation kit

It is imperative that all parts in contact with the solution should be carefully rinsed after treating each patient as follows:

1. Disconnect the irrigation line from the kit.
2. Immerse it in clean water.
3. Press the pedal for one minute to rinse the Sterijoint (ref. 54), the Suprasson® handpiece (ref. 53) and the tip (ref. 52).
4. Disconnect the irrigation line from the Sterijoint, and discard.

5. Disconnect the Sterijoint unit (ref. 46), Suprasson® handpiece (ref. 47) and tip (ref. 44) from the cord.
6. Take out and empty the reservoirs.
7. Disinfect them.
8. Sterilize the reservoirs and the Sterijoint in an autoclave.
9. Clean the motor and peristaltic pump with alcohol every day.

Cleaning the polishing system

It is imperative that the powder circuit is purged before any cleaning takes place.

The powder drum

Carefully empty the drum of powder if the device is not to be used for several hours (overnight for example). Residual dampness of the air alters the powder's properties.

To clean the powder drum, unscrew the cap, empty the powder into a waste bin or suck out the residue using the dental chair's suction system.

The ProphyPen handpiece

After each use, unscrew the ProphyPen handpiece nozzle (ref. 20) and clean in an ultrasound tank for 10 minutes (fragments the powder), dry well and sterilize the nozzle and ProphyPen handpiece in an autoclave.

Cleaning the filters

The water filter

Cleaning the water filter must be done in relation to the hardness index of the water. Proceed as follows:

1. Shut off the water at the tap on the water supply side and purge the circuit.
2. Set the on/off button to the "O" position and unplug the mains lead.
3. Unscrew the two parts of the filter (ref. 57 and 59) using two 10mm open-ended spanners.

4. Extract and replace the cartridge (ref. 58).
5. Reassemble by reversing the operation, retightening the two parts of the filter (ref. 57 and 59).

The air filter

When water is present, purge by turning the black screw (ref. 33) clockwise.

To stop the purge, turn the screw in the opposite direction.

Replacing the fuses

1. Set the on/off button to the "O" position (ref. 60) and unplug the mains lead (ref. 61).
2. Insert the end of a screwdriver (ref. 62) in the slot to open the protective cover.
3. Replace the fuses (ref. 63) with fuses of the same type and value.
4. Push the fuse holder back into its housing.

Fan

1. Check regularly that the fan is free from dust.
2. Clean if necessary by gentle blowing after having disconnected the mains supply (ref. 34).

7. Troubleshooting

Fault	Possible cause	Solution
Does not operate. (general)	Unplugged or defective mains cord.	Check and/or replace the mains cord.
	Fuses blown.	Replace the fuses.
	Defective wall power outlet.	Contact your electrician.
	Device switched off (0).	Set the switch to the on position (I).
	Footswitch disconnected and/or defective.	Check and/or replace the footswitch.
No water but ultrasound working.	No water in the water supply.	Contact your approved installer.
	Water inlet pipe not connected.	Connect the water inlet.
	Incorrect spray setting.	Adjust the setting.
	Handpiece and/or tip blocked.	Disconnect the handpiece and check if water is coming out of the cord: <ul style="list-style-type: none"> - if not, return the device to Satelec After Sales Service; - if it is, unscrew the tip from the handpiece, using the syringe blow air through the tip and the handpiece. If the problem persists, return to Satelec After Sales Service.
Water not atomized by the ultrasound.	Unsuitable treatment or tip, or selected ultrasound power too low.	Change the treatment or tip or adjust the ultrasound.
	Tip not tightened correctly.	Re-tighten the tip.
	Bent or worn tip.	Replace the tip.
	Incorrect use of the handpiece and/or the tip.	Refer to the instructions.
	Faulty contact between the handpiece and its connector.	Push the handpiece into the connector more firmly, without turning it.
	Water present in the handpiece connector.	Dry the connector. Where there are black marks between the pins of the connector and/or the contacts of the handpiece, return to Satelec After Sales Service.
Ultrasound too strong or too weak.	Power not set correctly in relation to treatment.	Adjust the power (see instructions).
	Tip not tightened correctly.	Re-tighten the tip.
	Bent or worn tip.	Replace the tip.
	Incorrect use of the handpiece and/or tip.	Refer to the instructions.

7. Troubleshooting

Fault	Possible Cause	Solution
No ultrasound	Tip not tightened correctly.	Re-tighten the tip.
	Defective handpiece.	Replace the handpiece.
	Oxidization and/or black marks on the pins of the connector.	Clean the contacts or replace the cord or return to Satelec After Sales Service.
	Cord cut.	Replace the cord or return to Satelec After Sales Service.
Water leaking from the handpiece connector.	Worn handpiece o'ring.	Replace the o'ring.
Leaking from the end of the tip, pedal released.	Dirt in the water circuit (electric valve blocked open)	Clean the circuit and electric valve (check the electric valve is working properly) or return to Satelec After Sales Service.
No light with the Lux handpiece.	Blown bulb.	Replace the bulb.
	Incorrect use or setting of the mode.	Refer to the instructions.
	Electrical problem.	Return to Satelec After Sales Service.
Overheating of the Lux handpiece.	No air in the supply.	Contact your approved installer.
	Air inlet pipe not connected.	Connect the air inlet.
No powder, air or water.	No air or water coming from the pipes connected to the mains supply.	Contact your approved installer.
	Handpiece and / or nozzle blocked.	Disconnect the handpiece and check that air and powder are coming out of the cord: <ul style="list-style-type: none"> - If not, return the device to Satelec After Sales Service. - If they are, unscrew the handpiece nozzle, unblock the nozzle with a file and place the handpiece in an ultrasound tank for 10 minutes.
No powder or water but air present.	Powder: drum empty or too full.	Fill the drum or remove the excess.
	Water: water inlet not connected or incorrect water setting.	Check the water inlet pipe or adjust the water flow on the handpiece cord.
Energy saving mode not working properly.	Incorrect use or setting of the mode.	Refer to the instructions.
	Exiting from energy saving mode not possible.	Return to Satelec After Sales Service.
	Adjustment impossible.	Return to Satelec After Sales Service.

8. Technical Specifications

Electrical

Class 1 type BF device

Intermittent duty 10 min / 5 min

Information display: LCD digital display with colored backlighting.

Option selection: press buttons.

Halogen bulb: 3.5 V / 750 mA / 15 000 lux (Prophy Max Lux only)

Light out delay: can be set between 0 and 9 seconds (Lux version only)

Mains power supply:

- 110 to 150 V AC, 50/60 Hz
- 230 to 260 V AC, 50/60 Hz

Fuse protection: dimensions 5x20 mm, 0.5 A (230 V) or 1 A (115 V).

Power consumption: 75 VA

Insulation: 4kV (1 ms)

Factory configuration:

- Power 5 in all ranges
- Light out delay: after three mins (Prophy Max Lux only)
- Standby: after nine mins

Manufactured in accordance with current legislation

CEM: conforms to the general IEC 601-1 and CEM CEI 60601-1-2 standards

Retention of the memorized configuration: 10 years at +25 °C

Mechanical

Water pressure at inlet: 1 to 5 bar (14 to 72 PSI)

Air pressure at inlet: 4.5 to 6 bar (65 to 87 PSI)

Operational

Water flow rate at exit: 0 to 120 ml/min

Frequency of vibrations: 27 to 33 kHz

Rate of flow of powder: 1 to 3 g/mn

Rate of flow of powder: 1 to 3g/min

Specifications of the powder: white to slightly yellow powder, impalpable, slightly salty taste and mint or lemon flavor.

Powder reservoir capacity : 50 g

Peristaltic pump: adjustable from 5 to 40 ml/min (optional)

Syringe function: rate of flow from about 100 ml/min (with the peristaltic pump)

Dimensions and use

Overall dimensions:

- Width: 270 mm
- Height : 110 mm
- Depth: 300 mm

Weight: 3 kg

Fixing: mount on a table


Operating temperature: +10 to +40 °C / 30 to 75% RH.

Storage/transport temperature: -40 to +70 °C / 10 to 100% RH

Symbols and markings

 Alternating current

 Type BF device

 Warning, please refer to the accompanying documentation.

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