# Cavitron® JET Plus Ultrasonic Scaler & Air Polishing Prophylaxis System with Tap-On™ Technology

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Please read carefully and completely before  
operating unit.

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

# Quick Start Installation Instructions


**Cavitron®**


- 1** Set unit to **OFF** position.  
Symbol: ○ means off.


- 2** Connect **POWER** supply.


- 3** Connect **WATER** supply by pushing **BLUE** water hose into water port.


- 4** Connect **AIR** supply (if equipped) by pushing **BLACK** air hose into air port.

- 5** Install **BATTERIES** in wireless foot pedal unit.


- 6** Connect **HANDPIECE** to cable.  
**TIP:** Handpiece is detachable and autoclavable. Make sure the electrical connections are dry before connecting.


- 7** Set unit to **ON** position.  
Symbol: | means on.  
The indicator light will be green when unit is on.


- 8** Hold upright and depress foot pedal until handpiece **FILLS WITH WATER**.  
**TIP:** adjust lavage control to higher flow rate to increase water flow.  
Tap-On™ Technology does not function unless an insert is in the handpiece.

- 9** Place **ULTRASONIC INSERT** into handpiece.

- 10** Adjust **POWER LEVEL**.

- 11** Adjust **LAVAGE CONTROL** on handpiece to desired setting.


- 12** Consult **Directions For Use** for daily maintenance, usage techniques and more.


# Quick Start User Guide

Cavitron®


## Tap-On™ Technology

To **activate** Tap-On™ Technology:  
Tap wireless foot pedal one time quickly.

Relax foot as you scale or air polish.


To **deactivate** Tap-On™ Technology:  
Tap wireless foot pedal one time.

## User controls

**Turbo:** When pressed, increases power delivered to the system by up to 25%; turbo power remains on until the button is pressed again.


**Boost:** Temporary, hands-free activation by pressing the foot pedal all the way to the floor; enables quick removal of tenacious calculus.  
**Tip:** Quick release of pedal from Boost will not deactivate Tap-On™ Technology.

**Purge:** Remove insert from Handpiece and press Purge button; water will purge through system for two minutes or until you press the foot pedal or press the Purge button again.

**Finer Lavage Water Control:**  
Conveniently adjust lavage water control directly on the handpiece;  
1 = lowest water flow  
6 = highest water flow  
> 6 = flush or purge

**Prophy Mode Auto Cycles**  
(available on Cavitron® JET Plus Systems): Automatically cycles between air polish and rinse without pumping the foot pedal. Choose short, medium or long prophy cycle times. Tap-On™ Technology will activate prophy mode auto cycles.

# Quick Reference Guide Diagnostic Display

## ON/OFF

Illuminates when the Main Power On/Off Power switch is in the “ON” (I) position.

## TURBO

Offers the ability to increase power to the system by up to 25% with the push of a button. Purple arrows illuminate when in use.

## BOOST

Illuminates when the Boost Mode is activated by the Tap-On™ Wireless Foot Pedal. To activate, fully depress Tap-On™ Foot Pedal to the second position (all the way to the floor). To deactivate, release Tap-On™ Foot Pedal to the first position.

## PURGE BUTTON

Illuminates when Purge function is activated. To activate Purge, remove insert from handpiece, turn the Handpiece Lavage Control to maximum water flow, press the Purge Button on the Diagnostic Display. Water will purge through system lines for two minutes. To deactivate during two minute cycle, press Purge button again or press Tap-On™ Foot Pedal.

## SERVICE

Illuminates when the system is not functioning properly. This display has three distinct modes:


- Slow blink (1 blink per second) means the system is not operating within factory specifications.
- Fast blink (3 blinks per second) indicates an improper set-up.
- Steady light indicates the system is overheating.

## LOW BATTERY

Illuminates when the Tap-On™ Foot Pedal battery power is approaching end of life. Replace batteries as instructed in section 7.9.

# Power Control

## Power Level Control

Turn knob to select ultrasonic power level for operation. Turning the knob clockwise increases the distance the insert tip moves (the stroke) without changing frequency; turning knob counter-clockwise decreases the distance the insert tip moves (the stroke) without changing the frequency.

## RINSE

Rinse mode is used during an ultrasonic scaling procedure when lavage is required to flush the procedural area. To activate, turn Power Level Control Knob fully counter-clockwise until a “click” is heard.

## BLUE ZONE

Provides an extended low-power range for improved patient comfort when subgingivally scaling.

# TABLE OF CONTENTS

## QUICK START INSTALLATION

<b>INSTRUCTIONS</b> . . . . .	2
<b>QUICK START USER GUIDE</b> . . . . .	3
<b>QUICK REFERENCE GUIDE:</b>	
<b>DIAGNOSTIC DISPLAY</b> . . . . .	4
<b>INTRODUCTION</b> . . . . .	6
<b>PRODUCT OVERVIEW</b> . . . . .	6
<b>TECHNICAL SUPPORT</b> . . . . .	6
<b>SUPPLIES &amp; REPLACEMENT PARTS</b> . . . . .	6
<b>INDICATIONS FOR USE</b> . . . . .	7
1.1 Ultrasonic Procedures . . . . .	7
1.2 Air Polishing Procedures . . . . .	7
<b>CONTRAINDICATIONS</b> . . . . .	7
<b>WARNINGS</b> . . . . .	7
<b>PRECAUTIONS</b>	
4.1 System Precautions . . . . .	8
4.2 Procedural Precautions . . . . .	8
Ultrasonics . . . . .	8
Air Polishing . . . . .	8
<b>ADVERSE REACTIONS</b> . . . . .	9
<b>INFECTION CONTROL</b>	
6.1 General Infection Control . . . . .	9
6.2 Water Supply Recommendations . . . . .	9
<b>INSTALLATION INSTRUCTIONS</b>	
7.1 Water Line Requirements . . . . .	9
7.2 Air Line Requirements &: Recommendations . . . . .	9
7.3 Electrical Requirements . . . . .	10
7.4 Unpacking the System. . . . .	10
7.5 System Installation . . . . .	10
7.6 Power Cord Connection. . . . .	10
7.7 Water Supply Line Connection . . . . .	10
7.8 Air Supply Line Connection . . . . .	11
7.9 Tap-On™ Foot Pedal Battery . . . . .	
Installation/ Replacement . . . . .	11
7.10 Tap-On™ Foot Pedal Synchronization . . . . .	11
<b>CAVITRON® JET Plus COMBINATION SYSTEM DESCRIPTION</b>	
8.1 System Controls . . . . .	13
8.2 Diagnostic Display Indicators and Controls . . . . .	14
8.3 Handpiece/Cable . . . . .	15
8.4 Cavitron® 30K™ Ultrasonic Inserts . . . . .	15
8.5 Cavitron JET Air Polishing Inserts	16

8.6 Tap-On™ Wireless Foot Pedal Operation . . . . .	16
8.7 Accessories and User Replaceable Parts . . . . .	16
8.7.1 Accessories . . . . .	16
8.7.2 User Replaceable Part Kits . . . . .	16

## SYSTEM SETUP, OPERATION AND TECHNIQUES FOR USE

9.1 Handpiece Setup . . . . .	17
9.2 Turbo Mode. . . . .	17
9.3 Boost Mode . . . . .	17
9.4 Patient Positioning . . . . .	17
9.5 Performing Ultrasonic Scaling Procedures . . . . .	17
9.6 Patient Comfort Considerations . . . . .	18
9.7 Air Polishing Powder Bowl . . . . .	18
9.8 Performing Air Polishing Procedures . . . . .	18
9.9 Proper Angulation of the Air Polishing Insert . . . . .	19

## SYSTEM CARE

10.1 Daily Maintenance. . . . .	20
Start-up procedures at the beginning of the day . . . . .	20
Between patients . . . . .	20
Shut-down procedures at the end of the day . . . . .	20
10.2 Weekly Maintenance . . . . .	21
10.3 Monthly Maintenance . . . . .	21
Water Line Filter Maintenance . . . . .	21
10.4 Air Supply Line Filter Maintenance . . . . .	21
10.5 Powder Bowl Maintenance . . . . .	21

## TROUBLESHOOTING

11.1 Troubleshooting Guide . . . . .	22
11.2 Technical Support and Repairs. . . . .	23

## WARRANTY PERIOD. . . . . 23

## SPECIFICATIONS. . . . . 23

## SYMBOL IDENTIFICATION . . . . . 24

## CLASSIFICATIONS . . . . . 24

## DISPOSAL OF UNIT . . . . . 24

## ELECTRO-MAGNETIC COMPATIBILITY

## PRECAUTIONS . . . . . 25-27

## QUICK REFERENCE GUIDE:

## TROUBLESHOOTING . . . . . 28

# INTRODUCTION

Congratulations!

Your decision to add the Cavitron® JET Plus Ultrasonic Scaler and Air Polishing Prophylaxis Combination System with Tap-On™ Technology to your practice represents a wise investment in good dentistry.

For over four decades, dental professionals have preferred the clinical benefits and labor-saving advantages inherent in Cavitron ultrasonic scalers. Clinical studies and independent research have confirmed the speed, efficiency and versatility of ultrasonic scaling.\*

With the addition of air polishing capabilities in the Cavitron JET Plus Combination System, your Cavitron JET Plus system becomes a compact prophylaxis center that optimizes the time spent performing scaling and polishing procedures and minimizes the need for strenuous calculus and stain removal with hand instruments. Clinical studies have proven that air polishing is far superior to traditional cup and pumice for stain and plaque removal.\* With proper technique and simple daily maintenance, your Cavitron® JET Plus™ Combination System will immediately become an indispensable component in your practice of modern preventive dentistry.

DENTSPLY Professional is an ISO 13485 registered company. All DENTSPLY Professional medical devices sold in Europe are CE marked in conformance with Council Directive 93/42/EEC.

Website: [www.professional.dentsply.com](http://www.professional.dentsply.com)

**CAUTION: United States Federal Law restricts this device to sale by or on the order of, a licensed dental professional.**

**For dental use only.**

## PRODUCT OVERVIEW

The Cavitron® JET Plus Combination System is a precision engineered and manufactured instrument. It contains controls and components for ultrasonic scaling and air polishing modes. In the scaling mode, the system produces 30,000 strokes per second at the ultrasonic insert's working tip that when combined with the cavitation effect of the coolant lavage creates a synergistic action that is designed to "power away" even the heaviest calculus deposits while maintaining operator and patient comfort. In the air polishing mode, the system delivers a precise air/water/powder mixture at the JET air polishing insert tip that polishes the tooth enamel without contact so there is less abrasion to enamel and no physical pressure or heat build-up that could cause discomfort in sensitive patients.

The Cavitron JET Plus Combination System is equipped with the Sustained Performance System™ (SPS Technology), which offers a constant balance between scaling efficiency and patient comfort by maintaining power when the insert tip encounters tenacious deposits, allowing the clinician to effectively scale even at a decreased/lower power setting. The Cavitron JET Plus System has extended the SPS technology by expanding the Blue Zone range, providing finer resolution to the power settings.

Advanced features that make the Cavitron JET Plus a wise investment include a Tap-On™ Wireless Foot Pedal with Tap-On™ Technology, Turbo Mode, Prophy Mode Auto Cycles, illuminated diagnostic display, rinse setting, automated purge function, JET-Mate™ detachable sterilizable handpiece, and 330° swivel handpiece cable with more precise lavage water control. These features combined with a low power range (Blue Zone™) and hands-free Boost Mode are designed to deliver a positive ultrasonic scaling and air polishing experience for your patients while providing your practice with the quality and reliability you've come to expect from Cavitron brand ultrasonic systems.

The Cavitron JET Plus Combination System is UL/ULC certified and approved. The Cavitron JET Plus Combination System is classified by Underwriters Laboratories Inc. with respect to electric shock, fire, mechanical hazards in accordance with the IEC 60601 Standard. The Cavitron JET Plus Combination System complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: 1) this device may not cause harmful interference, and 2) this device must accept any interference received, including interference that may cause undesired operation. Cavitron JET Plus base FCC certification/registration number: FCC ID: TF3-DPD81842; IC: 4681B-81842. Cavitron JET Plus Tap-On™ Foot Pedal FCC certification/registration number: FCC ID: TF3-DPD81861; IC: 4681B-81861. The term IC before the certification/registration number signifies that the Industry Canada technical specifications were met.

## TECHNICAL SUPPORT

For technical support and repair assistance in the U.S., call the DENTSPLY Professional Cavitron Care<sup>SM</sup> Factory Certified Service at 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For other areas, contact your local DENTSPLY® Professional representative.

## SUPPLIES & REPLACEMENT PARTS

To order supplies or replacement parts in the U.S., contact your local DENTSPLY Professional Distributor or call 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For other areas, contact your local DENTSPLY Professional Representative.

\* data on file

# SECTION 1: Indications For Use

The Cavitron® JET Plus System is designed for use in general prophylaxis and periodontal treatments and other areas of operative dentistry. It can also be used for:

## Ultrasonic Procedures

- All general supra and subgingival scaling applications
- Periodontal debridement for all types of periodontal diseases
- Endodontic procedures

## Air Polishing Procedures

- Removal of a variety of extrinsic stains, including tobacco, coffee, tea, and chlorhexadine
- Prophylaxis of orthodontic patients
- Preparing tooth surfaces prior to bonding and sealant procedures

# SECTION 2: Contraindications

- Ultrasonic Systems should not be used for restorative dental procedures involving the condensation of amalgam.
- Cavitron® PROPHY-JET Prophy Powder is a water-soluble Sodium Bicarbonate powder. Therefore, this powder is not recommended for patients on a sodium restricted diet. Cavitron® JET-Fresh Prophy Powder is a sodium free powder and can be used on patients who are on sodium restricted diets.
- Patients who have severe respiratory illness should consult their physician before undergoing air polishing prophylaxis procedures. Use is not recommended for patients with known allergies to the components of Cavitron Brand Prophy Powders.

# SECTION 3: Warnings

- The use of High Volume Saliva Evacuation to reduce the quantity of aerosols released during treatment is highly recommended.
- Prior to beginning treatment, patients should rinse with an antimicrobial such as Chlorhexidine Gluconate 0.12%. Rinsing with an antimicrobial reduces the chance of infection and reduces the number of microorganisms released in the form of aerosols during treatment.
- It is the responsibility of the Dental Healthcare Professional to determine the appropriate uses of this product and to understand:
  - the health of each patient,

- the dental procedures being undertaken,
- and applicable industry and governmental agency recommendations for infection control in dental healthcare settings,
- requirements, and regulations for safe practice of dentistry; and
- these Directions for Use in their entirety, including Section 4 Precautions, Section 6 Infection Control, and Section 10 System Care.

- Where asepsis is required or deemed appropriate in the best professional judgment of the Dental Healthcare Professional, this product should not be used.
- Handle Cavitron insert with care. Improper handling of insert, specifically the insert tip, may result in injury and/or cross contamination.
- Failure to follow properly validated sterilization processes and approved aseptic techniques for Cavitron inserts may result in cross contamination.
- Persons fitted with cardiac pacemakers, defibrillators and other active implanted medical devices, have been cautioned that some types of electronic equipment might interfere with the operation of the device. Although no instance of interference has ever been reported to DENTSPLY, we recommend that the handpiece and cables be kept 6 to 9 inches (15 to 23 cm) away from any device and their leads during use. There are a variety of pacemakers and other medically implanted devices on the market. Clinicians should contact the device manufacturer or the patient's physician for specific recommendations. This unit complies with IEC 60601 Medical Device Standards.
- Do not direct the air polishing stream at soft tissue or into the sulcus. Tissue emphysema has been reportedly caused when the air/water/powder stream was directed at the soft tissue or into the sulcus. If tissue emphysema occurs, consult a healthcare Professional.
- Insufficient water flow could result in elevated water and tip temperature. When operated at the input water temperature specified in the Water Line Requirements Section 7.1 and with sufficient water flow, the water and tip temperature should not exceed 50° C (122° F). Failure to follow recommendations for environmental operating conditions, including input water temperature, could result in injury to patients or users. If temperature is elevated, increase water flow. If temperature remains elevated, discontinue use.
- Empty the air polishing powder bowl at the end of the day to prevent "caking" of the powder and clogging of the lines and JET Air Polishing insert nozzle. Residual prophy powder in threads of the bowl and cap can result in excessive wear and disengagement of the cap during unit operation. Be sure to clean threads regularly as per Section 10 System Care. Check o-ring and threads on powder bowl cap and replace if worn.
- During boil-water advisories, this product should not be operated as an open water system (e.g. connected

to a public water system). A Dental Healthcare Professional should disconnect the system from the central water source. The Cavitron DualSelect system can be attached to this unit and operated as a closed system until the advisory is cancelled. When the advisory is cancelled, flush all incoming waterlines from the public water system (e.g. faucets, waterlines and dental equipment) in accordance with the manufacturer's instructions for a minimum of 5 minutes.

- Per FCC Part 15.21, changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.
- This portable transmitter with its antenna complies with FCC/IC RF exposure limits for general population / uncontrolled exposure.
- This Device complies with Industry Canada License-exempt RSS standard(s). Operation is subject to the following two conditions: 1) this device may not cause interference, and 2) this device must accept any interference, including interference that may cause undesired operation of the device.

## SECTION 4: Precautions

### 4.1 System Precautions

- Cavitron® Prophy Powders are specially formulated for use in Cavitron® Air Polishing Systems. Do not use any other materials in the air polishing powder reservoir. Use of other materials may damage the system and void the warranty.
- Do not modify powder with additives. Additives may clog your air polishing system.
- Close manual shut-off valve on the dental office water supply every night before leaving the office.
- Do not place the system on or next to a radiator or other heat source. Excessive heat may damage the system's electronics. Place the system where air is free to circulate on all sides and beneath it.
- The system is portable, but must be handled with care when moving.
- Equipment flushing and dental water supply system maintenance are strongly recommended. See Section 10: System Care.

### 4.2 Procedural Precautions

#### General

- As with all dental procedures, use universal precautions (i.e., wear face mask, eyewear, or face shield, gloves and protective gown).

## Ultrasonics

- The Cavitron JET Plus unit works with Cavitron inserts as a system, and was designed and tested to deliver maximum performance for all currently available Cavitron brand ultrasonic inserts. Companies that manufacture, repair or modify inserts carry the sole responsibility for proving the efficacy and performance of their products when used as a part of this system. Users are cautioned to understand the operating limits of their insert before using in a clinical setting.
- Like bristles of a toothbrush, ultrasonic inserts "wear" with use. Inserts with just 2 mm of wear lose about 50% of their scaling efficiency. In general it is recommended that ultrasonic inserts be discarded and replaced after one year of use to maintain optimal efficiency and avoid breakage. A DENTSPLY Professional Insert Efficiency Indicator is enclosed for your use.
- If excessive wear is noted, or the insert has been bent, reshaped or otherwise damaged, discard the insert immediately.
- Ultrasonic insert tips that have been bent, damaged, or reshaped are susceptible to in-use breakage and should be discarded and replaced immediately.
- Retract the lips, cheeks and tongue to prevent contact with the insert tip whenever it is placed in the patient's mouth.

## Air Polishing

- Patients should wear safety glasses or eye protection during air polishing treatment.
- Patients wearing contact lenses should remove them prior to air polishing treatment.
- Patients who have severe respiratory illness should consult their physician before undergoing air polishing prophylaxis procedures.
- Avoid use on cementum or dentin.
- Direct contact of prophy powder with surfaces and marginal areas of dental restorations should be avoided.
- Set the air polishing powder flow control to the maximum (H) position only when it is necessary to remove particularly difficult stains. Return the powder flow control to the medium position upon the completion of the procedure.
- JET Air Polishing Insert nozzles that have been bent, damaged or re-shaped, are susceptible to in-use breakage and should be discarded and replaced immediately.

## SECTION 5: Adverse Reactions

- Tissue emphysema may occur if the air/water/powder stream was directed at the soft tissue or into the sulcus.
- Allergic reactions to components of Cavitron Brand Prophy Powders have been known to occur. If an allergic reaction occurs, seek medical attention immediately.

## SECTION 6: Infection Control

### 6.1 General Infection Control

- For operator and patient safety, carefully practice the infection control procedures detailed in the Infection Control Information Booklet accompanying your System. Additional booklets can be obtained by calling Customer Service at 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For areas outside the U.S., contact your local DENTSPLY Professional representative.
- As with high speed handpieces and other dental devices, the combination of water and ultrasonic vibration from the Cavitron JET Plus Combination System will create aerosols. Following the procedural guidelines in Section 9 of this manual can effectively control and minimize aerosol dispersion.

### 6.2 Water Supply Recommendations

- It is highly recommended that all dental water supply systems conform to applicable CDC (Centers for Disease Control and Prevention) and ADA (American Dental Association) standards, and that all recommendations be followed in terms of flushing, chemical flushing, and general infection control procedures. See Sections 7.1 and 10.
- As a medical device, this product must be installed in accordance with applicable local, regional, and national regulations, including guidelines for water quality (e.g. drinking water). As an open water system, such regulation may require this device to be connected to a centralized water control device. The Cavitron® DualSelect™ Dispensing System may be installed to allow this unit to operate as a closed water system.

## SECTION 7: Installation Instructions

Anyone installing a Cavitron JET Plus System should observe the following requirements and recommendations.

### 7.1 Water Line Requirements

- A water supply line with user-replaceable filter is supplied with your system. See Section 10 System Care for replacement instructions.
- Incoming water supply line pressure to the system must be 20 psi (138 kPa) to 40 psi (275 kPa). If your dental water system's supply line pressure is above 40 psi, install a water pressure regulator on the water supply line to your Cavitron JET Plus Combination System.
- A manual shut-off valve on the dental water system supply line should be used so that the water can be completely shut-off when the office is unoccupied.
- In addition to the water filter supplied, it is recommended that a filter in the dental water system supply line be installed so that any particulates in the water supply will be trapped before reaching the Cavitron system.
- After the above installations are completed on the dental water supply system, the dental office water line should be thoroughly flushed prior to connection to the Cavitron system.
- Incoming water temperature to the Cavitron System should not exceed 25°C (77°F). If needed a device should be installed to maintain a temperature within this specification, or a Cavitron DualSelect Dispensing System attached to allow this system to be operated as a closed water system.

### 7.2 Air Line Requirements & Recommendations

- An air supply line with a user-replaceable filter assembly is supplied with your Cavitron JET Plus Combination System. Refer to Section 7.8 Air Supply Line Connection.
- Incoming air supply line pressure to the system must be 65 psig (448 kPa) to 100 psig (690 kPa). If your office air line pressure is above 100 psig (690 kPa), install an air pressure regulator on the supply line to your Cavitron JET Plus Combination System.
- A manual shut-off valve on the office air supply line should be used so that the air line can be completely shut-off, and the line pressure relieved when the office is unoccupied.
- The Cavitron System must be supplied with clean, dry air to help prevent water condensation from forming in the air supply line which may cause it to malfunction. In addition to the air filter supplied with your System, it is strongly recommended that an air dryer be used on the compressor line supplying the Cavitron System to prevent moisture from reaching the System, which can cause "caking" of the air polishing powder and clogging of the lines and air polishing insert nozzle.

## 7.3 Electrical Requirements

- Incoming power to the system must be 100 volts AC to 240 volts AC, single phase 50/60 Hz capable of supplying 1.0 amps.
- The system power should be supplied through the AC power cord provided with your system.
- **WARNING:** To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

## 7.4 Unpacking the System

Carefully unpack your Cavitron JET Plus Combination System and verify that all components and accessories are included:

1. Cavitron® JET Plus™ Combination System with Handpiece Cable Assembly with swivel
2. Air Line Assembly (Black) with Filter and Quick Disconnect
3. Water Line Assembly (Blue) with Filter and Quick Disconnect
4. Additional Water Line Filter
5. Detachable AC Power Cord (not shown)
6. Cavitron® Tap-On™ Wireless Foot Pedal
7. "AA" Batteries (4-Pack)
8. Auxiliary Cable for Tap-On™ Foot Pedal
9. Cavitron® JET Air Polishing Insert with cleaning tool
10. JET-Mate Detachable Sterilizable Handpiece
11. Prophy Handpiece Cleaning Wire (not shown)
12. Cavitron® Ultrasonic Inserts (quantity optional)
13. Efficiency Indicator for Cavitron Inserts
14. Literature Packet
15. PROPHY-JET® Sodium Bicarbonate Prophy Powder
16. JET-Fresh® Aluminum Trihydroxide Prophy Powder (may not be included in all kits)
17. Powder Removal Container

## 7.5 System Installation

- The Cavitron JET Plus Combination System is designed to rest on a level surface. Be sure unit is stable and resting on four feet.
- The Cavitron JET Plus Combination System should not be positioned such that access to the power input and AC power cord are limited.
- Placing unit in direct sunlight may discolor plastic housing.

- The system has been equipped with a Cavitron® Tap-On™ Wireless Foot Pedal which was factory synchronized to operate with the system's base unit. If your office has more than one Cavitron system with Tap-On Technology, it is recommended that you mark the Tap-On Foot Pedal and base unit for easy reference as to which Tap-On Foot Pedal operates with which base unit. Should resynchronization be necessary, follow the instructions in section 7.10.

## 7.6 Power Cord Connection

- Verify the Main Power ON/OFF switch, located at the center front underside of the System, is set to the OFF (O) position before proceeding.

- Insert the AC power cord into the power input on the back of the System.
- Insert the pronged plug into an AC wall outlet.


## 7.7 Water Supply Line Connection

- Grasp the Water Supply Line (blue hose) by the end opposite the quick-disconnect and insert it into the water inlet connector until fully seated.

- Connect the quick-disconnect to the dental office water supply or a Cavitron DualSelect Dispensing System.
- Inspect all connections to make certain there are no leaks.
- To remove the water line from the Cavitron JET Plus Combination System, turn off the dental office water supply. Disconnect the water supply line from the dental office water supply. If a quick-disconnect connector is attached to the end of the hose, relieve the water

pressure by pressing the tip of the connector in an appropriate container and allow water to drain. To remove the hose from the system, push on the outer ring of the system's water inlet and gently pull out the water line.

Press ring to release water supply tube.

## 7.8 Air Supply Line Connection

- Grasp the Air Supply Line (black hose) by the end opposite the quick-disconnect and insert it into the air inlet connector until fully seated.

- Connect the quick-disconnect to the dental office air supply or a Cavitron DualSelect Dispensing System.
- Inspect all connections to make certain there are no leaks.
- A filter mounting bracket is included for hanging the air filter. Mount the bracket to a suitable vertical surface and slide the filter onto the bracket. The clear bowl should hang downward allowing for moisture separation and drainage of water from the air filter. See Section 10 System Care for replacement instructions.
- To remove the air supply line from the Cavitron JET Plus Combination System, turn off the dental office air supply. Disconnect the air supply line from the dental office air supply, then push on the outer ring of the system's air inlet and gently pull out the air line. If a quick-disconnect connector is attached to the end of the hose, relieve the air pressure by pressing the tip of the connector and allowing the air to escape.

Press ring to release air supply tube.

## 7.9 Tap-On™ Foot Pedal Battery Installation/Replacement

- Turn Tap-On™ Wireless Foot Pedal over and using a Philips screwdriver carefully remove battery cover screw and battery cover. If applicable, remove used batteries and install two new “AA” batteries as shown. Do not depress Tap-On™ Foot Pedal while installing batteries.

Look for blinking communications light.

- The communication light will blink for approximately two seconds to indicate the Tap-On™ Foot Pedal's ability to communicate with the unit. If the light does not blink, check the batteries. If the batteries are good and the light doesn't blink, a communications error may exist. To re-establish communication with Tap-On™ Foot Pedal review Synchronization procedure, section 7.10.
- The remote frequency communication can be bypassed using the auxiliary Tap-On™ Foot Pedal cable. Refer to Section 11.2 Technical Support and Repair for further action.
- Replace the battery cover and screw and hand tighten cover with Philips screwdriver.
- Remove batteries if Tap-On™ Foot Pedal is to be stored for an extended period of time.

## 7.10 Tap-On™ Foot Pedal Synchronization

The Tap-On™ Wireless Foot Pedal supplied with your system has been factory synchronized with the base unit. Should a replacement Tap-On™ Foot Pedal be necessary, synchronization will be required prior to system operation. Perform the following steps to synchronize the Tap-On™ Foot Pedal with the base unit.


- Turn the MainPower switch located at the center front underside of the system to the OFF (O) position.

2. Install a new set of “AA” batteries into the foot control (see section 7.9). Leave the battery cover of the Tap-On™ Foot Pedal open so the red push button is accessible.
3. Maintain a distance of no more than 10 feet (3 meters) between the base unit and Tap-On™ Foot Pedal during the synchronization process.
4. Remove any inserts from the handpiece and adjust the Power Level Control out of Rinse Mode. Turn the Main Power switch to the ON (I) position and wait for the Diagnostic Display graphics to light (refer to Section 8.2).
5. While all graphics are lit, press the Purge button, located on the Diagnostic Display.

The graphics will begin to blink in a sequential pattern, representing the synchronization mode. This mode will last 5 to 6 seconds.

6. During this mode, press the red button located in the battery compartment of the Tap-On™ Foot Pedal. This will complete the synchronization process.

7. Synchronization is successful when all graphic lights blink at the same time.
8. To verify proper communication, press the foot control to the Boost position (Tap-On™ Foot Pedal fully depressed – 2<sup>nd</sup> position) and ensure the Boost graphic on base unit illuminates.
9. Attach battery cover and tighten the screw.
10. In the event communication cannot be established, temporarily use the supplied Auxillary Tap-On™ Foot Pedal Cable to connect the Tap-On™ Foot Pedal directly to the unit.

# SECTION 8: Cavitron JET Plus Combination System Description

## 8.1 System Controls

### Ultrasonic Power Level Control

Turn knob to select the ultrasonic power level for operation. Turning the knob clockwise increases the distance the insert tip moves (stroke) without changing the frequency; turning the knob counterclockwise decreases the distance the insert tip moves (stroke) without changing the frequency.

**The Blue Zone** is an extended low-power range for improved patient comfort when subgingivally scaling.

### Rinse

Turn the ultrasonic power level control knob fully counterclockwise until a “click” is heard. Rinse mode is for use during an ultrasonic scaling procedure when lavage is desired with minimal cavitation.

### Powder Flow Control

Rotate clear Powder Flow Control to adjust powder flow rates. For minimum powder flow turn control clockwise to “L”. For maximum flow turn to “H”.

### Diagnostic Display

See section 8.2 for more details.

### Handpiece Holder

Securely holds the system’s handpiece, or cable connector when handpiece is not installed.

### JET-Mate™ Handpiece

Accepts all Cavitron® 30K™ Ultrasonic inserts and Cavitron JET Air Polishing Inserts. Automatically selects air polishing or scaling mode. See section 8.3 for more details.

### Main Power ON/OFF Switch


ON/OFF Switch located at the center front underside of the system.

### Prophy Mode Auto Cycles

Automatically cycle between air polish and rinse without pumping foot pedal. Choose short, medium or long prophy cycle times. See section 9.8 for more details.

### Tap-On™ Wireless Foot Pedal

Eliminates the need to hold down or pump foot pedal. See section 8.6 for more details.

## 8.2 Diagnostic Display Indicators and Control

### Service Indicator

Illuminates when the system is not functioning properly. This display has three distinct modes.

- A fast blink (3 blinks per second) indicates an improper set-up.
- A slow blink (1 blink per second) means the system is operating out of factory specifications.
- A steady light indicates the system is overheating.

Refer to section 11.1 for Troubleshooting guidelines.

### Power Indicator


Illuminates (3 sec. delay) when the Main Power ON/OFF Control Switch is ON ("I" position).

### Low Battery Indicator

Illuminates when the Tap-On™ Foot Pedal battery power is approaching end of life. Replace batteries as instructed in section 7.9.

### Boost Indicator

Illuminates when the Boost Mode has been activated with the Tap-On™ Foot Pedal.

### Turbo Mode

When pressed, extra power is delivered to the system, up to 25%. This will increase the stroke of the tip. Turbo power remains on until the button is pressed again or unit is turned off. (Purple arrows illuminate when Turbo Mode is on).

### Purge Control

Illuminates when the Purge function is activated. To activate Purge, remove insert from the handpiece, turn the handpiece lavage control to maximum water flow and press the Purge button. Water will purge through system for 2 minutes. To deactivate mode during the 2 minute cycle, press Purge button again or press Tap-On™ Foot Pedal.

The Purge Control is also used during the Tap-On™ Foot Pedal Synchronization process. See section 7.10 for more details.

## 8.3 Handpiece / Cable


### Lavage Control

Turn the Lavage Control to select flow rate during system operation. Flow rate is based on a scale from 1 to 6. Turn clockwise toward 6 to increase flow at insert tip. Turn counter-clockwise toward 1 to decrease flow. The flow rate through the handpiece also determines the temperature of the lavage. Lower flow rates produce warmer lavage. Higher water flow rates produce cooler lavage.

If the handpiece becomes warm, increase the flow rate. With experience the Dental Healthcare Professional will be able to determine the best flow rate setting for optimum operating efficiency and patient comfort.

### Swivel Feature

Reduces cable drag as handpiece rotates during procedures.

### Soft Nozzle Grip

Ergonomically designed to provide for a comfortable grasp of the handpiece. The grip is a replaceable wear component. Prior to use, verify that the grip is flush with the hard plastic of the insert port.

### Powder Delivery Port


Creates an airtight seal between the air polishing insert and the handpiece. Replace when wear is noticed or powder is leaking at nozzle interface.

### Insert Port

The Cavitron® JET-Mate Sterilizable Handpiece accepts all Cavitron® 30K Ultrasonic Inserts and JET Air Polish inserts.

## 8.4 Cavitron 30K Ultrasonic Inserts

The many styles of Cavitron and Cavitron Bellissima 30K Ultrasonic Inserts are easily interchangeable for various procedures and applications. See enclosed literature for specific information.

## 8.5 Cavitron JET Air Polishing Insert

**Insert Heater Rod:** Heats the delivered water for patient comfort.

**Prophy Powder Delivery Tube:** Directs air/powder flow to insert tip.

**Insert Marking:** Manufacturer, Date (YDDD=Single digit year and triple digit day of the year).

**Air Polishing Insert Nozzle:** Tube-in-a-tube design delivers precise air/water/powder mixture at point of delivery.

**O-Ring:** Seals water when Insert is fully seated in Handpiece. O-ring should be replaced when worn.

## 8.6 Tap-On™ Technology Wireless Foot Pedal Operation

### Using Foot Pedal in Tap-On™ Mode

For scaling procedures, Tap-On™ Technology eliminates the need to hold the pedal down. Tapping the foot pedal once activates ultrasonic power or rinse mode for approximately 4 minutes. Tapping the foot pedal while in Tap-On™ mode disables the ultrasonic power and water flow. Boost is still available while scaling in Tap-On™ mode. To use boost, simply depress the foot pedal to the second position (all the way to the floor) to activate and hold as long as boost is desired. Release foot pedal to return to Tap-On™ mode.

For prophy procedures, Tap-On™ Technology and Prophy Mode Auto Cycles eliminate the need to pump the pedal by automatically cycling between rinse and polish. Tapping the foot pedal once will enable an automatic air polishing/rinse cycle that lasts for approximately 1 minute. Tapping the pedal a second time disables the automatic air polishing/rinse cycle. See Prophy Mode Auto Cycles in section 9.8 for cycle details.

### TIPS:

-  Tap-On™ Technology will not run water unless an insert or air polishing insert is in the handpiece
-  A sensor in the handpiece holder will prevent Tap-On™ Technology from operating when the handpiece is in the holder.
-  If the foot pedal is not tapped quickly, it will function in a conventional manner.

### How to Disable and Enable Tap-On™ Technology

Tap-On™ Technology feature can be disabled by simultaneously holding the Purge and Turbo buttons for a period of approximately 5 seconds. The two buttons will blink approximately 6 times. When the buttons are released, they will blink an additional 6 times to confirm Tap-On™ has been disabled.

Tap-On™ Technology feature can be enabled by simultaneously holding the Purge and Turbo buttons for a

period of approximately 5 seconds. The two buttons will blink approximately 6 times to confirm Tap-On™ has been enabled.

### Using Foot Pedal without Tap-On™ Mode

For scaling operation, the first position activates both the ultrasonic energy and lavage at the insert tip. The second position activates the Boost Mode. The Boost Mode (fully depressed Tap-On™ Foot Pedal) increases the ultrasonic power level for quick removal of tenacious deposits without adjusting the power level knob. To deactivate Boost Mode, release Tap-On™ Foot Pedal to first position.

For prophy operation, the first position activates Rinse Mode. The second position activates Air Polishing Mode. (Boost Indicator will not illuminate.)

## 8.7 Accessories and User Replaceable Parts

### 8.7.1 Accessories

1. AC Power Cord
2. Tap-On™ Technology Wireless Foot Pedal
3. Auxiliary Tap-On™ Foot Pedal Power Cable
4. Cavitron JET-Mate Sterilizable Handpiece
5. Prophy Handpiece Cleaning Wire
6. Cavitron 30K Ultrasonic Inserts
7. Cavitron DualSelect Dispensing system
8. Cavitron JET Air Polishing Insert
9. Cavitron JET Nozzle Cleaning Tool
10. Cavitron Prophy-Jet Sodium Bicarbonate Prophy Powder
11. Cavitron Jet Fresh Aluminum Trihydroxide Prophy Powder

### 8.7.2 User Replaceable Part Kits

1. Powder Bowl Cap O-Ring, Part Number 628052001
2. Powder Bowl Cap, Part Number 81728
3. Cavitron Insert Replacement O-ring Kits, 12/packs  
Part Number 62351 (black) for plastic and soft grips  
Part Number 62605 (green) for metal grips and air polishing insert
4. Handpiece Cable O-Ring, Part Number 79357
5. JET-Mate Handpiece Nozzle Grip, 81717
6. Lavage (Water) Filter, 10/Pack, Part Number 90158

For detailed information, contact your local DENTSPLY Professional Representative or authorized DENTSPLY Professional Distributor.

# SECTION 9: System Setup, Operation and Techniques for Use

## 9.1 Handpiece Setup

- Follow Precautions listed in the general and ultrasonic sections of 4.2 Procedural precautions. This handpiece is sterilizable. Refer to Infection Control Information Booklet for sterilization instructions prior to using handpiece.
- Connect the Handpiece to the Cable Assembly by aligning the electrical connections. If Cable Assembly does not seat into the handpiece, gently rotate the handpiece until contacts align, then fully insert handpiece.
- Hold empty handpiece in a semi-upright position over a sink or drain. Activate the Tap-On™ Foot Pedal until water exits to release any air bubbles that might be trapped inside the handpiece. Avoid letting water into the Powder Delivery Port as clogging may result. NOTE: Tap-On™ Technology only operates when an insert or an air polishing insert is in the handpiece.
- Lubricate the O-ring on the insert with water before placing it into the handpiece. Fully seat insert with a gentle push-twist motion. DO NOT FORCE. If using the air polishing insert, align the powder delivery tube with the powder delivery port and gently push into the handpiece until fully seated. DO NOT FORCE.
- Turn the Lavage Control to select flow rate during system operation. Flow rate is placed on a scale from 1 to 6. Turn control clockwise toward 6 to increase flow at insert tip. Turn control counter-clockwise toward 1 to decrease flow. The flow rate through the handpiece also determines the temperature of the lavage. Lower

water flow rates produce warmer lavage. Higher flow rates produce cooler lavage. If the handpiece becomes warm, increase the flow rate. With experience the Dental Healthcare Professional will be able to determine the best flow rate setting for optimum operating efficiency and patient comfort.

## 9.2 Turbo Mode

Pressing the “Turbo” button on the scaling unit increases the unit’s ultrasonic power up to 25%. When you need more power for an extended period of time, simply press the Turbo button on the display panel (arrows will illuminate to show you are in “Turbo” mode).

DENTSPLY recommends that the clinician familiarize themselves with the available power levels throughout the power knob rotation in both normal and turbo modes. To do this, simply hold your favorite insert over the sink and adjust the power knob while observing the inserts spray pattern and toggling between both normal and turbo modes.

## 9.3 Boost Mode

Boost provides a temporary increase in ultrasonic scaling power for quick removal of tenacious calculus without touching the unit. Boost is activated by fully depressing the Tap-On™ Foot Pedal to second position (all the way to the floor). When boost is activated, the boost icon will illuminate on the display panel. Boost remains on as long as the clinician has the foot pedal pressed all the way down. In order to deactivate boost, release the Tap-On™ Foot Pedal to first position.

## 9.4 Patient Positioning

For optimal access to both the upper and lower arches, the backrest of the chair should be adjusted as for other dental procedures. This assures patient comfort and clinician visibility. Have the patient turn his/her head to the right or left. Also position chin up or down depending upon the quadrant and surface being treated. Evacuate irrigant using either a saliva ejector or High Volume Evacuator (HVE).

## 9.5 Performing Ultrasonic Scaling Procedures

**Note: Refer to the Infection Control Information booklet supplied with your system and Section 10 of this manual for general procedures to be followed at the beginning of each day and between patients.**

- Follow precautions listed in the General and Ultrasonics sections of 4.2 Procedural Precautions.
- The edges of Cavitron Ultrasonic Inserts are intentionally rounded so there is minimal danger of tissue laceration with proper ultrasonic scaling technique. Whenever the insert tip is placed in the patient’s mouth, the lips,

cheek and tongue should be retracted to prevent accidental (prolonged) contact with the activated tip.

- Turn Power Level Control to select ultrasonic power level for operation. Clockwise increases system power. Power level will increase throughout the full range of the control. Hold the handpiece over a sink or drain. While in Tap-On™ mode, simply tap the Tap-On™ Foot Pedal to activate the system. (If Tap-On™ mode is turned off, press and hold the Tap-On™ Foot Pedal down to activate the system.) Check water spray to verify fluid is reaching the working end of the insert tip. Adjust the water lavage control until the water (lavage) flows with a rapid drip or small spray. Higher water flow settings provide cooler irrigation.
- It may be necessary to adjust lavage with the system in “Boost” mode (Tap-On™ Foot Pedal fully depressed) so adequate fluid will be available to cool tip to tooth interface.
- In general, it is suggested that a “feather-light-touch” be used for ultrasonic scaling. The motion of the activated tip and acoustic effects of the irrigating fluid, in most cases, are adequate to remove even the most tenacious calculus.
- Periodically check the Cavitron Ultrasonic Insert for wear with the Cavitron Insert Efficiency Indicator.
- The use of a saliva ejector or High Volume Evacuator (HVE) is recommended during all procedures.
- Set the system’s Power Level Control to the lowest efficient power setting for the application and the selected insert.
- Keep the foot pedal near your foot to make it convenient to access.

## 9.6 Patient Comfort Considerations

### **Reasons for sensitivity**

- Incorrect tip placement. The point should never be directed toward tooth root surfaces.
- Not keeping tip in motion on tooth. Do not allow the insert to remain in a static position on any one area of the tooth. Change the insert’s path of motion.
- Applying excessive pressure. Use a very light grasp and pressure, with a soft tissue fulcrum whenever possible, especially on exposed cementum.

If sensitivity persists, decrease power setting and/or move from the sensitive tooth to another and then return.

## 9.7 Air Polishing Powder Bowl

- Use only Cavitron® Prophy Powders in your Cavitron JET Plus Combination System. Any other substance or additives may clog the system and will void the warranty. For your convenience, the prophy powders are supplied in bottles. Keep stored in a location that does not exceed 95° F.

- A special container is provided with your System for use in emptying the powder bowl.
- It is strongly recommended that the powder bowl be emptied at the end of each day. This will reduce moisture absorption and minimize clogging.

### **To fill, or refill, the powder bowl:**

- Turn the System OFF.
- Unscrew the Powder Bowl Cap.
- With the cap of the powder bottle closed, shake the powder bottle vigorously to break up any lumps that may have formed from settling. Carefully pour powder into the bowl until the level reaches the top of the center tube.
- Using a soft dry cloth, remove powder adhering to the cap and bowl threads. Secure the cap on the powder bowl.
- Turn the System ON.

**NOTE: Use only Cavitron® Prophy Powders in the system. Powder should be kept dry and stored between 32°F/0°C and 95°F/ 35°C.**


### **To adjust the flow of powder:**

- Adjust the powder flow rate by positioning the control pointer on the cap at H (12 o’clock), M (9 o’clock) or L (6 o’clock).
- For heavy stain removal, set the control to H.
- For light stain removal, set the control to L.
- The control can be set at any position between H and L.
- The view window at the center of the pointer lets you observe the powder flow (small white circle of powder) during operation. If no flow is seen, check for clogging or add prophy powder.

## 9.8 Performing Air Polishing Procedures

- Follow precautions listed in the General and Air Polishing sections of 4.2 Procedural Precautions.
- Place a 2 x 2 gauze on lip.
- Select the proper amount of powder and water to create the slurry needed using the Powder Flow Control on the powder bowl cap and the Lavage Control on the handpiece cable. Ensure that the system’s Power Level Control is in the “Prophy Mode” range. Use more powder for heavy stains and less powder for light stains. With experience the Dental Healthcare Professional will be able to determine the best flow rate settings for optimum efficiency and patient comfort. Never operate the system with powder only.
- Flush the patient’s tongue with water to help reduce the saline taste.

- The recommended normal procedure is to clean 1-3 teeth with the air polishing spray and then rinse the area with water in order to inspect the work site before proceeding to the next 1-3 teeth. When using Tap-On™ Technology and Prophy Mode Auto Cycle, the system will automatically cycle between rinse and air polish. When Tap-On™ Technology is turned off, the air polishing spray is activated when Tap-On™ Foot Pedal is depressed to the second position (all the way to the floor) and rinse is activated when the Tap-On™ Foot Pedal is in the first position. If desired, the bleed air passing through the air polishing insert tip can be used to dry the work site during inspection (bleed air occurs when the Tap-On™ Foot Pedal is released).

Prophy Mode	Air Polish	Rinse
MANUAL *	Press foot pedal to second position	Press foot pedal to first position

Prophy Mode Auto Cycles	Approximate Cycle Times	
	Air Polish	Rinse
SHORT **	0.75 sec	1.25 sec
MEDIUM **	2.0 sec	1.0 sec
LONG **	3.0 sec	2.0 sec

- Use your free hand and the patient's cheeks or lips to form a cup to contain aerosols. Tilt the patient's head toward you to help prevent puddling in the cupped lip and minimize aerosol dispersion. Rinse excessive slurry from the patient's mouth thoroughly and often.

- Maintain a 2 to 4 mm tip-to-tooth operating distance. Keep the tip in constant circular motion and maintain a sweeping motion from interproximal to interproximal. When air polishing the anteriors, center the spray on the middle third of the tooth. The edge of the spray will clean the teeth to the gingiva. Refer to section 9.9 Proper Angulations for all tooth surfaces.

- Use adequate evacuation. Use of a high-speed suction (High-Volume Evacuator) with the aid of a dental assistant is recommended. When performing air polishing without the aid of a dental assistant, the use of a saliva ejector and/or aerosol-reduction device is recommended.

- Do not aim directly at the soft tissue.

- Avoid use on surfaces and marginal areas of dental restorations.


- Prophy Mode Auto Cycles: When performing air polishing procedures, it is recommended to clean 1-3 teeth with air polishing spray and then rinse the area with water in order to inspect the work site before proceeding to the next 1-3 teeth. Prophy Mode Auto Cycles allow for automatic cycling between air polishing and rinse while the foot pedal is in Tap-On™ mode. A small burst of air is released before the Prophy Mode Auto Cycle begins to alert the clinician that the cycle is beginning.

- If at any time continuous air polishing is needed, simply depress the foot pedal all the way to the floor. Release of the foot pedal will disable Prophy Mode Auto Cycle.


\*The 'MANUAL' selection does not automatically cycle between air polish and rinse, it allows user to air polish conventionally. Press foot pedal to first position for rinse and then press foot pedal to second position for air polishing. Continue to repeat as needed.

\*\*'SHORT', 'MEDIUM', and 'LONG' selections will automatically cycle between air polish and rinse at approximately the times listed in the table above.


## 9.9 Proper Angulation of the Air Polishing Insert

Recommended angulation on the anterior teeth is 60° with the tip aimed at the middle third of the tooth surface.

Recommended angulation on the buccal and lingual surfaces of posterior teeth is 80° with the tip aimed slightly distally.

Recommended angulation to occlusal surfaces is 90°

# SECTION 10: System Care

It is recommended that you perform the following maintenance procedures.

## 10.1 Daily Maintenance

### START-UP PROCEDURES AT THE BEGINNING OF THE DAY:

1. Open the manual shut-off valve on the dental office water supply system.
2. With the Cavitron JET Plus Combination System OFF, unscrew the powder bowl cap. Verify the powder bowl is empty. Turn the system ON for 15 seconds to eliminate residual moisture in the lines. Turn the system OFF.
3. Shake the powder bottle well to create an even consistency of powder mixture.
4. Pour enough powder into the bowl for the procedure to be performed. With experience the Dental Healthcare Professional will be able to determine the amount of powder required. Do not fill above the top of the center tube.
5. Secure the cap on the powder bowl.
6. Install a sterilized JET-Mate Handpiece onto the handpiece cable.
7. Set the Power Level Control to minimum and the Lavage Control to maximum.
8. Turn the system ON.
9. If powder fluffing is observed when the Tap-On™ Foot Pedal is not in use, this would indicate an air leak. To correct, turn the System OFF, remove the Powder Cap, clean any residual powder from the O-ring seal and threads, replace the Powder Cap, tighten, and turn the System back ON.
10. Hold the sterilized handpiece (without an insert or nozzle insert installed) over a sink or drain. Activate the Purge Control button.
  - The Purge button will light for two minutes indicating proper activation of the purge function.
  - If the Purge button is activated with an insert present in the handpiece, the button will blink for 3 seconds and disable. Remove the insert from the handpiece and press the Purge button again.
  - The Purge function can be interrupted at any time during the two minute cycle by pressing the Purge button again or by pressing the Tap-On™ Foot Pedal.
11. After completing the purge cycle, place a sterilized 30kHz Cavitron® Ultrasonic Insert into the handpiece and set the Power Level Control and Lavage Control to your preferred operating position for ultrasonic scaling. For air polishing, place a sterilized JET Air Polishing Insert into the handpiece and adjust the Power Level Control to Prophy Mode, and the Powder

Flow and Lavage Controls to your preferred operating positions.

### BETWEEN PATIENTS:

1. Remove the used Cavitron® Ultrasonic Insert or JET Air Polishing Insert. Clean and sterilize following the Infection Control Procedures that were enclosed with your insert.
2. Hold the handpiece over a sink or drain and activate Purge function as described in Step 10 of the Start-Up procedure.
3. After the purge cycle is complete, turn the System to the OFF (0) position.
4. Remove the JET-Mate handpiece, clean and sterilize following the procedures outlined in the Cavitron Systems Infection Control Procedures booklet that was enclosed with your system.
5. Disinfect the surfaces of the cabinet, Power Cord, Handpiece Cable, Tap-On™ Foot Pedal and Cable assembly (if applicable), Water Supply and Air Supply lines by applying an approved non-immersion type disinfectant solution\* carefully following the instructions provided by the disinfectant solution manufacturer. To clean the system, generously spray disinfectant solution on a clean towel and wipe all surfaces. Discard used towel. Dry with a clean cloth. To disinfect the system, generously spray disinfectant on a clean towel and wipe all surfaces. Allow disinfectant solution to air dry. Never spray disinfectant solution directly on the system.
6. Inspect the handpiece cable for any breaks or tears.
7. If using a closed water supply or DualSelect Dispensing System, check for adequate fluid volumes for the next patient.
8. Check the powder bowl for sufficient powder for the next procedure.
9. When ready to use, place a sterilized JET-Mate Handpiece onto the handpiece cable assembly and insert a sterilized ultrasonic insert or air polishing insert into the handpiece and adjust system controls as preferred.

**\*NOTE: Broad spectrum hospital grade water-based disinfection solutions are preferred, such as Lysol IC. Some alcohol-based disinfectant solutions may be harmful and may discolor plastic materials.**

### SHUT-DOWN PROCEDURES AT THE END OF THE DAY:

1. Follow the "Between Patients" maintenance procedures, steps 1 through 6. In addition, it is recommended to close the manual shut-off valve on the dental water supply system.

2. Unscrew the powder bowl cap.
3. Remove the powder bowl from the unit and discard the unused powder.
4. Holding the open end of the powder bowl away from you, activate the system for 15 seconds to clear the bowl. A high volume evacuator can be used to remove any residual powder.
5. Remove the o-ring seal from the powder bowl cap and using a soft dry cloth, wipe residual powder from the cap, the o-ring and the powder bowl threads. Be careful not to scratch or otherwise damage the cap.
6. Return the o-ring to the cap and secure the cap on the powder bowl.

## 10.2 Weekly Maintenance

- Remove residual prophylaxis powder from the cap and bowl threads using a soft brush (toothbrush). If not removed, caked powder in threads can result in thread wear and powder bowl cap disengagement.
- It is strongly recommended that this system be disinfected by chemically flushing the waterlines with a 1:10 Sodium Hypochlorite solution (NaOCl) at the end of each week. This can be accomplished by connecting this device to the Cavitron DualSelect Dispensing System or a number of other devices available from your local distributors. Where this device is connected to the Cavitron DualSelect Dispensing System, please follow the DualSelect system's Directions for Use manual. If connected to another device, please follow those directions for use, keeping in mind that a chemical flush should be performed at maximum water flow for at least 30 seconds. The system should be left undisturbed for 10 minutes but no more than 30 minutes to allow the sodium hypochlorite solution to soak in the lines. As a suggestion, it is recommended that a sign be placed on the system stating that the SYSTEM IS BEING DISINFECTED WITH A STRONG DISINFECTANT AND SHOULD NOT BE USED. When ready, flush system with clean water for at least 30 seconds or until sodium hypochlorite odor disappears. ALL CHEMICALS MUST BE FLUSHED FROM THE SYSTEM BEFORE IT IS READY FOR PATIENT USE.

## 10.3 Monthly Maintenance

### WATER LINE FILTER MAINTENANCE:

When the water line filter becomes discolored, the filter should be replaced to prevent reduced water flow to the Cavitron JET Plus system. A 10-pack of replacement filters is available by ordering Part Number 90158 from your local DENTSPLY Professional distributor.

1. Verify that the system is turned OFF.
2. Disconnect the water supply hose from the dental office water supply. If a quick disconnect connector is

attached to the end of the hose, relieve the water pressure by pressing the tip of the connector in an appropriate container and drain the water.

3. Grasp the fittings on either side of the filter disk and twist counterclockwise. Remove the filter section from either side of the water hose.
4. Install the replacement filter onto the water hose fittings. The filter should be positioned to match up with the correct hose fitting.
5. Hand tighten one hose fitting onto filter in a clockwise direction. Tighten second hose onto filter in clockwise direction. Reconnect the water supply line, operate unit to bleed the air and test for leaks.

## 10.4 Air Supply Line Filter Maintenance

Water build up in the air supply line filter should be drained. This can be accomplished by turning the knob on the bottom of the filter counter-clockwise to open. After draining, turn the knob fully clockwise to close (some filters drain by depressing the Schrader valve stem at the bottom of the filter). If the inner filter element becomes discolored or dirty, a new filter assembly should be installed. Replacement filter assemblies are available by ordering Part Number 90088 from your local DENTSPLY Professional distributor.

1. Verify that the system is turned OFF.
2. Disconnect the air supply line from the dental office air source.
3. Using a 7/16" wrench, loosen the nuts on the side fittings of the filter. Unscrew the nuts and slide them down the hose. Disconnect hoses from the filter and discard the used filter.
4. Insert the short hose into the input port of the filter and the long hose into the output port of the filter. Slide the nuts up the hoses and screw onto the fittings. Tighten using pliers or a wrench.
5. Turn the System ON, operate the system and check for leaks.

## 10.5 Powder Bowl Maintenance

1. Turn the System OFF.
2. Allow the powder bowl to depressurize and unscrew the Powder Cap.
3. Empty powder from the bowl and use the high suction to remove any residual powder in the bowl.
4. Turn the System ON and check for strong air flow from the center tube of the powder bowl.
5. If no or low air flow is present, turn the System OFF.
6. Unscrew the knurled ring at the bottom of the bowl assembly and remove the fitting assembly.
7. Using the JET Air Polishing Insert nozzle cleaning wire tool, clean clogged powder from the fitting assembly. Turn the System ON and check for strong air flow.

- Turn the system OFF.
- Check that the o-ring is properly positioned in the groove of the fitting assembly and reassemble the fitting assembly to the bowl. Tighten knurled ring. Place powder bowl into System.
  - Fill the powder bowl with fresh prophylaxis powder and test for flow and leaks.
  - Remove residual powder from thread on the cap and bowl with a soft, dry cloth.

## SECTION 11:

### Troubleshooting

Although service and repair of the Cavitron JET Plus Combination System should be performed by DENTSPLY personnel, the following are some basic troubleshooting procedures that will help avoid unnecessary service calls. Generally, check all lines and connections to and from the System. A loose plug or connection will often create problems. Check the settings on the System's controls.

#### 11.1 Troubleshooting Guide

##### Symptom:

##### System operates: Tap-On™ Technology is not working

- Tap-On™ Technology might be disabled. Refer to section 8.6.
- Check to see if handpiece is in holder. Tap-On™ Technology is disabled when handpiece is in holder.
- Check to see if the insert is secured inside the handpiece. Tap-On™ Technology is disabled when there is no insert in the handpiece.

##### Symptom:

##### System will not operate: No Power ON indicator

- Check that the Main Power Switch is in the ON (I) position, and that the detachable Power Cord is fully seated in the receptacle on back of System.
- Check that the system's power cord plug is fully seated in an approved AC wall outlet.
- Check that the wall outlet is functional.

##### Symptom:

##### System will not operate: Power ON Indicator is illuminated

- If the office has more than one Tap-On™ Foot Pedal, test each to ensure that the proper Tap-On™ Foot Pedal is being used. With a handpiece and insert installed, depress the Tap-On™ Foot Pedal to the first position. The system should dispense water. If none of the Tap-On™ Foot Pedals operate the system, continue to the next step.
- Resynchronize one Tap-On™ Foot Pedal to the system (see Section 7.10 Tap-On™ Foot Pedal Synchronization).

##### Symptom:

##### System operates: No water flow to insert tip or

#### handpiece overheats

- Assure that handpiece lavage control is properly adjusted.
- Check for clogged insert. Replace insert if necessary.
- Check that dental office water supply valves are open.
- If the system is connected to DualSelect Dispensing System, check that fluid level in the selected bottle is sufficient. Make sure valves are open when using external water source.
- Check that the water line filter is clean. Replace filter if needed.

##### Symptom:

##### System operates: No insert cavitation

- Check that the Power Level Control is not in Rinse Mode.
- Check the insert for damage and that it is properly installed in the handpiece.
- Check that the handpiece is properly installed to the cable assembly.
- Verify that the soft nozzle grip is flush with the hard plastic of the insert port.
- Turn the system's Main Power Switch to the OFF (O) position. Wait 5 seconds and turn the system back ON.
- If problem still exists, replace both "AA" batteries in Tap-On™ Foot Pedal with new "AA" batteries (Refer to section 7.9) or connect the Auxiliary Tap-On™ Foot Pedal Cable.

##### Symptom:

##### System operates: Purge Mode will not function – icon flashing

- Check that there is no insert in the handpiece.
- Check that handpiece is properly installed to the cable assembly.

##### Symptom:

##### System operates: Service Indicator blinking

- Fast Blinking (3 blinks per second) – Indicates improper set-up
- If insert is in the handpiece, remove. Verify the handpiece is properly seated and hold the foot control for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step.
  - Attach a NEW handpiece and hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. Discard the old handpiece or return if within warranty. If blinking remains, continue to the next step.
  - Install and fully seat an insert into handpiece. Hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step.
  - Install and fully seat a NEW insert in handpiece and hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, system is ready for use. Discard old insert or return if within warranty. If blinking remains, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

- Slow Blinking (1 blink per second)
  - The system is not operating within factory specifications.
    1. Remove insert.
    2. Turn Main Power Switch OFF, (O) position. Wait five seconds. Turn unit ON, (I) position.
    3. Operate Purge function.
    4. If service indicator still blinks, refer to section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

### Symptom:

#### System operates: Service Indicator illuminated

1. Ensure that the base unit has adequate ventilation and is not near a heat source (i.e. radiator, heat lamp, sunlight or other heat producing operatory equipment).
2. Turn Main Power Switch to the OFF (O) position. Allow system to cool for 10 minutes and turn system to the ON (I) position. Verify light is not illuminated.
3. If light is still illuminated, refer to section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

### Symptom:

#### System operates: Air Polishing Insert nozzle blocks repeatedly

1. Powder is contaminated (lumpy). Discard powder.
2. Air Supply Line Air Filter is contaminated. Refer to Section 10.4 Air Supply Line Air Filter Maintenance.
3. Dental office air source should be serviced to eliminate the source of the contamination.

### Symptom:

#### System operates: No bleed air

1. Blocked JET Air Polishing insert nozzle. Clean nozzle using supplied tool.
2. Blocked air bleed “duckbill” air filter. Refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

### Symptom:

#### System operates: No or poor cleaning action

1. Very low powder level or empty powder bowl. Fill Powder bowl.
2. Blocked JET Air Polishing insert nozzle. Clean nozzle using supplied tool.
3. Powder Cap loose. Turn Main Power switch to the OFF (O) position. Tighten Powder Cap to powder bowl and turn System ON. If cap does not fit tightly, check for thread wear and replace cap, o-ring seal, or bowl assembly.
4. Clogged fitting assembly on powder bowl. Refer to Section 10.5 Powder Bowl Maintenance.
5. Dental office air source should be serviced to eliminate the source of the contamination.

### Symptom:

#### System operates: Continuous powder agitation

1. Powder Cap not securely sealed. Turn Main Power switch to the OFF (O) position and remove Powder Cap.

2. Remove the o-ring seal from the Powder Cap and clean residual powder from the cap. Be careful not to scratch or otherwise damage the plastic cap.
3. Wipe off the o-ring and place it in the Powder Cap. Tighten Powder Cap to Powder Bowl and turn system ON. Worn caps and o-rings should be replaced when wear is noted.

## 11.2 Technical Support and Repairs

For technical support and repair assistance call DENTSPLY Professional Cavitron Care<sup>SM</sup> Factory Certified Service at 1-800-989-8826 Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For areas outside the U.S., contact your local DENTSPLY Professional representative.

## SECTION 12: Warranty Period

The Cavitron JET Plus Combination Ultrasonic Scaler and Air Polishing System is warranted for TWO YEARS from date of purchase. The JET-Mate Handpiece enclosed with your system is warranted for SIX MONTHS from date of purchase. Refer to the Warranty Statement Sheet furnished with your system for full Warranty Statement and Terms.

## SECTION 13: Specifications

Electrical Voltage	Continuous (100-240 VAC)
Current	1.0 Amperes, Maximum
Phase	Single
Frequency	50/60 Hertz
Water Pressure	20 to 40 psig (138 to 275 kPa)
Water Temperature	< 25°C (77°F)
Air Pressure	65 to 100 psig (448 to 600 kPa)
Water Flow Rate	Minimum Setting (CCW) < 15 ml/min Maximum Setting (CW) > 55 ml/min
Weight	4.4 lbs (2 Kg)
Dimensions	Height: 6 in. (15,24 cm) Width: 9.5 in. (24,13 cm) Depth: 8 in. (20,32 cm) Handpiece Cable length: 6.5 ft. (2.0 M) Auxillary Footswitch Cable length: 8 ft. (2.4 M) Water Supply Line length: 8 ft. (2.4 M) Air Supply Line length: 10 ft. (3.04 M)

# SECTION 13: Specifications, continued

Footswitch	Protection Class IPX1. Not for operating theatres.
Remote Communication	Frequency: 2405 to 2480 MHz
	Power: < 1mW
	Channels: 16
Operating Environment	Temperature: 15 to 40 Deg. Celsius (59 to 104 Deg. Fahrenheit)
	Relative Humidity: 30% to 75% (non-condensing)
Transport and Storage Conditions	Temperature: -40 to 70 Deg. Celsius (-40 to 158 Deg. Fahrenheit)
	Relative Humidity: 10% to 100% (non-condensing)
	Atmospheric Pressure: 500 to 1060 hPa

# SECTION 14: Symbol Identification

-  AC POWER
-  TYPE B APPLIED PART EQUIPMENT
-  PROTECTIVE EARTH (GROUND)
- IPX1** Footswitch not for operating theatres  
Protection Class- IPX1  
IPX1 Classification of ingress of water
-  Consult Instructions for Use
- O/I** AC Power Switch (O = Off, I = On)
-  Footswitch

This symbol is a mandatory marking for devices entering the European market to indicate conformity with the essential health and safety requirements set out in European Directives. The symbol may be accompanied by a four-digit identification number of the notified body.

MEDICAL EQUIPMENT WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA-C22.2 No. 601.1, ANSI/AAMI ES60601-1 (2005, 3rd ed.) CAN/CSA-C22.2 No. 60601-1 (2008), 13VA

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:  
1) this device may not cause harmful interference, and  
2) this device must accept any interference received, including interference that may cause undesired operation.

Foot pedal: FCC ID:TF3-DPD81861  
IC: 4681B-81861  
Base: FCC ID:TF3-DPD81842  
IC: 4681B-81842

Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union

This is a wireless device.

# SECTION 15: Classifications

- Type of protection against electric shock: Class 1
- Degree of protection against electric shock: Type B
- Degree of protection against the harmful ingress of water: Ordinary
- Mode of operation: Continuous
- Degree of safety of application in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide: Equipment not suitable for use in the presence of flammable anaesthetic or oxygen.
- According to medical device directive: IIA (rule 9)

# SECTION 16: Disposal of Unit

U.S. - Dispose of the system components in accordance with state and local laws.  
EU - Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union.

# SECTION 17: Electro-Magnetic Compatibility Precautions

## Guidance And Manufacturer's Declaration - Electromagnetic Emissions

The Ultrasonic Scaler model G137 is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrasonic Scaler should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ultrasonic Scaler uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Ultrasonic Scaler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	The Ultrasonic Scaler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

## Guidance And Manufacturer's Declaration – Electromagnetic Immunity

The Model G137 is intended for use in the electromagnet environment specified below. The customer or the end user of the Model should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+6 kV contact +8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## Guidance And Manufacturer's Declaration - Electromagnetic Emissions

The Model G137 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model G137 should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance Level	Electromagnetic environment - guidance
Radiated RF  IEC 61000-4-3	3 V/m  26 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Model G137, including cables, than the Recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><math>d = 1.7 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup> should be less than the compliance level in each frequency range<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: right;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model G137 is used exceeds the applicable RF compliance level above, the Model G137 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Model G137.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distance between Portable and mobile RF communications equipment and the model @ 3Vrms

The model G137 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model G137 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model G137 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{v_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	-	0.12	0.23
0.1	-	0.34	0.74
1	-	1.7	2.3
10	-	3.7	7.4
100	-	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Recommended separation distance between Portable and mobile RF communications equipment and the model @ 10Vrms

The model G137 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model G137 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model G137 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{v_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	-	0.035	0.07
0.1	-	0.11	0.22
1	-	0.35	0.7
10	-	1.1	2.2
100	-	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# QUICK REFERENCE GUIDE: TROUBLESHOOTING

SYMPTOM	ACTION TAKEN
<b>System will not operate: No Power ON Indicator</b>	<ol style="list-style-type: none"> <li>1. Check that the Main Power Switch is in the ON (I) position, and that the detachable power cord is fully seated in the receptacle on back of system.</li> <li>2. Check that the system's power cord plug is fully seated in an appropriate AC wall outlet.</li> <li>3. Check that the wall outlet is functional.</li> </ol>
<b>System will not operate: Power ON Indicator is illuminated</b>	<ol style="list-style-type: none"> <li>1. If the office has more than one Tap-On™ Foot Pedal, test each to ensure that the proper Tap-On™ Foot Pedal is being used. With a handpiece and insert installed, depress the Tap-On™ Foot Pedal to the first position. The system should dispense water. If none of the Tap-On™ Foot Pedals operate the system, continue to the next step.</li> <li>2. Resynchronize one Tap-On™ Foot Pedal to the system (see Section 7.10 Tap-On™ Foot Pedal Synchronization).</li> </ol>
<b>System operates: No water flow to insert tip or handpiece overheats</b>	<ol style="list-style-type: none"> <li>1. Assure that handpiece lavage control is properly adjusted.</li> <li>2. Check for clogged insert. Replace insert if necessary.</li> <li>3. Check that dental office water supply valves are open.</li> <li>4. If the system is connected to DualSelect Dispensing System, check that fluid level in the selected bottle is sufficient. Make sure valves are open when using external water source.</li> <li>5. Check that the water line filter is clean. Replace filter, if needed.</li> </ol>
<b>System operates: No insert cavitation</b>	<ol style="list-style-type: none"> <li>1. Check that the Power Control is not in Rinse Mode.</li> <li>2. Check the insert for damage and that it is properly installed in the Handpiece.</li> <li>3. Check that the handpiece is properly installed to the cable assembly.</li> <li>4. Verify that the soft nozzle grip is flush with the hard plastic of the insert port.</li> <li>5. Turn the system's Main Power Switch OFF, (O) position. Wait 5 seconds and turn the system back ON.</li> <li>6. If problem still exists, replace both "AA" batteries in Tap-On™ Foot Pedal with new "AA" batteries (Refer to Section 7.9) or connect Auxiliary Tap-On™ Foot Pedal Cable.</li> </ol>
<b>System operates: Service indicator blinking</b>	<ol style="list-style-type: none"> <li>1. Fast blinking (3 blinks per second) – Indicates improper set-up <ol style="list-style-type: none"> <li>A. If insert is in the handpiece, remove. Verify the handpiece is properly seated and hold the Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step.</li> <li>B. Attach a NEW handpiece and depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. Discard the old handpiece or return if within warranty. If blinking remains, continue to the next step.</li> <li>C. Install and fully seat an insert into handpiece. Depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for operation. If blinking remains, continue to the next step.</li> <li>D. Install and fully seat a NEW insert in handpiece and depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, system is ready for use. Discard old insert or return if within warranty. If blinking remains, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.</li> </ol> </li> <li>2. Slow blinking (1 blink per second) - The system is not operating within factory specifications. <ol style="list-style-type: none"> <li>A. Remove insert.</li> <li>B. Turn Main Power Switch OFF, (O) position. Wait five seconds. Turn unit ON, (I) position.</li> <li>C. Operate Purge function.</li> <li>D. If service indicator still blinks, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.</li> </ol> </li> </ol>
<b>System operates: Service indicator illuminated</b>	<ol style="list-style-type: none"> <li>1. Ensure that the base unit has adequate ventilation and is not near a heat source (i.e. radiator, heat lamp, sunlight or other heat producing operatory equipment).</li> <li>2. Turn Main Power Switch OFF (O) position. Allow system to cool for 10 minutes and turn system to the ON (I) position. Verify light is not illuminated.</li> <li>3. <u>If light is still illuminated, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.</u></li> </ol>
<b>System operates: Purge mode will not function - icon flashing</b>	<ol style="list-style-type: none"> <li>1. Check that there is no insert in the handpiece.</li> <li>2. Check that Handpiece is properly installed to the cable assembly.</li> </ol>
<b>System operates: Air Polishing Insert Nozzle blocks repeatedly</b>	<ol style="list-style-type: none"> <li>1. Powder is contaminated (lumpy). Discard powder.</li> <li>2. Air Supply Line Air Filter is contaminated. Refer to Section 10.4 Air Supply Line Air Filter Maintenance.</li> <li>3. Dental office air source should be serviced to eliminate the source of the contamination.</li> </ol>
<b>System operates: No bleed air</b>	<ol style="list-style-type: none"> <li>1. Blocked JET Air Polishing insert nozzle. Clean nozzle using supplied tool.</li> <li>2. Blocked air bleed "duckbill" air filter. Refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.</li> </ol>
<b>System operates: No or poor cleaning action</b>	<ol style="list-style-type: none"> <li>1. Very low powder level or empty powder bowl. Fill Powder Bowl.</li> <li>2. Blocked JET Air Polishing insert nozzle. Clean nozzle using supplied tool.</li> <li>3. Powder Cap loose. Turn Main Power switch to the OFF (O) position. Tighten Powder Cap to powder bowl and turn System ON. If cap does not fit tightly, check for thread wear and replace cap, o-ring seal, or bowl assembly.</li> <li>4. Clogged fitting assembly on powder bowl. Refer to Section 10.5 Powder Bowl Maintenance.</li> </ol>
<b>System operates: Continuous powder agitation</b>	<ol style="list-style-type: none"> <li>1. Powder Cap not securely sealed. Turn Main Power switch to the OFF (O) position and remove Powder Cap.</li> <li>2. Remove the o-ring seal from the Powder Cap and clean residual powder from the cap. Be careful not to scratch or otherwise damage the plastic cap.</li> <li>3. Wipe off the o-ring and place it in the Powder Cap. Tighten Powder Cap to powder bowl and turn system ON. Worn caps and o-rings should be replaced when wear is detected.</li> </ol>
<b>System operates: Continuous powder agitation</b>	<ol style="list-style-type: none"> <li>1. Tap and Scale Mode might be disabled. Refer to Section 8.6.</li> <li>2. Check to see if handpiece is in holder. Tap and Scale is disabled when handpiece is in holder</li> <li>3. Check to see if the insert is secured inside the handpiece. Tap and Scale is disabled when there is no insert in the handpiece.</li> </ol>

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**United States of America**

DENTSPLY Professional  
Technical Service and Repair Department  
1301 Smile Way  
York, PA 17404-1785  
Phone: (800) 989-8826 or (717) 767-8502

**Deutschland**

DENTSPLY DeTrey GmbH  
De-Trey-Strasse 1  
78467 Konstanz  
Germany  
Phone: 7531 583 0

**France**

DENTSPLY DeTrey  
17 Michael FARADAY  
78380 Montigny Le Bretonneux  
France  
Phone: (1) 30 14 77 77

**Australia**

DENTSPLY (Australia) Pty. Ltd  
11-21 Gilby Road  
Mount Waverley, Victoria 3149  
Australia  
Phone: (61) 3-9538-8280

**United Kingdom**

DENTSPLY Ltd.  
Building 1  
Aviator Park  
Station Road  
Addlestone  
Surrey  
KT1 52PG  
United Kingdom  
Phone: (0) 800 072 3313

**Italia**

DENTSPLY DeTrey Italia s.r.l.  
Via A. Cavaglieri, 26  
I-00173 Roma  
Italia  
Phone: (06) 723 3626

**Canada**

DENTSPLY Canada  
161 Vinyl Court  
Woodbridge, Ontario  
L4L 4A3 Canada  
Phone: (905) 851-6060

**DENTSPLY**  
**PROFESSIONAL**

Manufactured by:  
DENTSPLY Professional  
DENTSPLY International  
1301 Smile Way  
York, PA 17404-1785 USA

DENTSPLY DeTrey GmbH  
De-Trey-Str. 1  
78467 Konstanz  
Germany

Imported and  
Distributed by:  
DENTSPLY Canada  
Woodbridge, Ontario  
L4L 4A3