

Product : **PATHFILE**
Ref: **A0015****Ref doc. :** **F19 02 105.X/ 03 / 2008 updated 03/2011**

		# Page
GB	DIRECTION FOR USE	3
FR	PROTOCOLE D'UTILISATION	7
ALL	GEBRAUCHSANWEISUNG	11
ESP	INSTRUCCIONES DE USO	15
IT	ISTRUZIONI PER L'USO	19
SWE	BRUKSANVISNING	23
DK	BRUGSANVISNING	27
FIN	KÄYTTÖOHJEET	31
POR	PROTOCOLO DE UTILIZAÇÃO	35
NL	GEBRUIKSVOORSCHRIFT	39
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΓΙΑ	43
POL	INSTRUKCJA OBSŁUGI	48

	GB Symbols	F Symboles	D Symbole	E Simbolos	I Simboli
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	See directions for use	Voir mode d'emploi et instructions d'usage	Siehe Gebrauchsanweisung	Ver modo de empleo e instrucciones de uso	Seguire le istruzioni d'uso
	Opened packages are not replaced	Emballage ouvert, non remplacé	Geöffnete Verpackung, nicht ersetzbar	Los embalajes abiertos no se cambian	Le confezioni aperte non sono sostituibili
	Batch number	N° de lot	LOT - Nummer	Número de lote	Numero di lotto
	Cannot be sold separately	Vente séparée, non autorisée	Wird nicht einzeln verkauft	No puede ser vendido separadamente	Non può essere venduto separatamente
	Assortment	Assortiment	Sortiment	Surtido	Assortimento
	Nickel titanium	NiTi	Nickel-Titan	Niquel Titanio	Nickel-titanio
	Green dot	Point vert	Der Grüne Punkt	Punto verde - recicable	Punto verde

	P Simbolos	S Symboler	FIN Symbolit	DK Symbol	NL Symbolen	GR ΣΥΜΒΟΛΑ	POL Symboler
	Fabricante	Tillverkaren	tekijä	Producent	Fabrikant	Κατασκευαστής	producent
	Ver as instruções de uso	Se bruksanvisningarna	Katsokaa käyttöohjeet	Se brugsanvisningerne	Zie gebruiksaanwijzing	Δείτε οδηγίες χρήσεως	Zobacz instrukcję obsługi
	Embalagens abertas não são substituídas	Öppnade förpackningar kan inte returneras	Avattuja pakkauksia ei 'Sivu: 2 vaihdeta	Åbnede pakninger ombyttes ikke	Geopende verpakkingen worden niet vervangen	Ανοιγμένες συσκευασίες δεν αντικαθίστώνται	Otwarte opakowania nie będą wymieniane
	N° da lote	Tillverkningsnummer	Erän numero	Batchnummer	Lotnummer	Αριθμός συσκευασίας	Numer serii
	Não pode ser vendido separadamente	Säljes inte separat	Ei voida myydä erikseen	Sælges ikke separat	Wordt niet apart verkocht	Δεν μπορεί να πωληθεί μεμονωμένα	Nie można sprzedawać osobno
	Sortido	Sortiment	Lajitelma	Lajitelma	Assortiment	Ασσορτί	Asortyment
	Niquel-titânio	Nickel Titan	Nikkeli-titaani	Nikkel-titan	Nikkel-titanium	Νικέλιο τιτανίου	Niklowo-tytanowy
	ponto verde	grön prick	Kierrätettävä	grøn prik	groene punt	πράσινο σημείο	Zielona kropka

For dental use only

DIRECTION FOR USE - PATHFILE™

0. Composition

The cutting part of these instruments is made of a nickel-titanium alloy.

1. Indications for use

Applications: catheterization of the root canal

These instruments are to be used only in a clinical or hospital environment, by qualified users. Users are responsible for checking whether the product can be used for the purpose intended.

2. Contra-indications

None known to date.

3. Warnings

None known.

4. Precautions

- Multiple use and re-processing cycles may lead to increased risk of file separation
- These instruments should not continuously be immersed in a sodium hypochlorite solution.
- Instrument decontamination: strictly follow decontamination instructions from the manufacturer
- Irrigate abundantly and frequently
- Before using the PathFile™ instruments, establish a glide path using hand files, to at least an ISO size 010.
- Use at constant rotation at a speed of 300 rpm.
- Clean flutes frequently and check for signs of distortion or wear.
- Use the Files with light pressure.
- Use the files to follow the canal to the working length, and then withdraw.

5. Adverse reactions

In the present technical state, no adverse reaction has been reported to date.

6. Step-by-step instructions for use

- Establish a manual Glide Path with standard stainless steel K-File # 008, # 010
- Identify the working length with the # 010 K-File in combination with an Apex Locator
- Irrigate
- Use the PathFile™ # 013 to working length
- Irrigate
- Use the PathFile™ # 016 to working length
- Irrigate
- Use the PathFile™ # 019 to working length
- Irrigate
- Before starting the canal shaping with NiTi ProTaper Universal files (or another system), confirm your working length with a K-File # 015 combined with an Apex Locator.

7. Disinfection, cleaning and sterilization

Reprocessing procedure for dental instruments and implantable radicular devices

Foreword

For hygiene and sanitary safety purposes, all instruments not marked “sterile” must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

Area of application

Disinfection and sterilisation before first usage and reprocessing procedures concerning:

A1. Instruments:

Cutting instruments, (hand and engine driven) such as:

Endodontic instruments (files, broaches, reamers, enlargers, endodontic burs, ultrasonic inserts);

Rotary cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs).

Root canal filling instruments (Pluggers, spreaders, compactors);

Supports, kits and instrument organisers

Hand instruments and clamps.

A2. Implantable devices:

Dentinal and radicular posts made of steel, titanium and glass fibers.

Supports, kits and organiser systems for posts.

A3. Contra angle

B. Filling material: Only chemical disinfection (no sterilisation)

Gutta percha, Thermafil obturation devices.

Exclusion

- Equipment such as Motors, Apex locators and other devices with reprocessing procedures included in the individual Direction for Use.
- MTA, Glyde, TopSeal.

General recommendation

1 - Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents

2 - For your own safety, please wear personal protective equipment (gloves, glasses, mask).

3 - The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.

4 - Limitations and restrictions on reprocessing:

The individual DFU indicates if the useful life of a device might be reduced by the number of reprocessing cycles.

Furthermore, the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of colour coding or marking, are indications that the devices are not able to fulfil the intended use with the required safety level.

5 - Single use marked instruments are not approved for re-use.

6- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.

7 - Tungsten carbide burs, plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H₂O₂) solution.

8 - NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.

9 - Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.

10- The washer-disinfector is not recommended for instruments made of aluminium, tungsten carbide or carbon steel.

Step-by-step procedure

A. Devices

			A3. Contra Angle			
			A2. Implantable devices			
			A1. Instruments		Following uses	
					First use	
	Operation	Operating mode	Warning			
1.	Disassembling	- Disassemble the device, if required	Silicone stops have to be removed.			X
2.	Pre-Disinfection	- Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible.	- Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor. - Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments (See general recommendations). - For visible impurities observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.			X
3.	Rinsing	- Abundant rinsing (at least 1 min)	- Use quality water in accordance with local regulations. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the cleaning.			X
4a.	Automated Cleaning with washer-disinfector	- Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (Ao value > 3000 or, at least 5 min at 90 °C).	- Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer-disinfector use kits, supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.		X	X X X X
OR						
4b.	Manual Cleaning or assisted by an ultrasonic device	- Place the devices in a kit, support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.		X	X X X
5.	Rinsing	- Abundant rinsing (at least 1 min)	- Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.		X	X X X
6.	Inspection	- Inspect devices and sort out those with defects. - Assemble the devices (stops)	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging. - For Contra Angle : lubricate the device with an adequate spray before packaging		X	X X X X

		A3. Contra Angle			
		A2. Implantable devices			
		A1. Instruments		Following uses	
				First use	
	Operation	Operating mode	Warning		
7.	Packaging	- Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in "Sterilisation pouches".	- Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.	X	X
8.	Sterilization	- Steam sterilisation at: 134 °C / 273°F-during 18 min.	- The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use only autoclaves that are matching the requirements of EN 13060, EN 285. - Use a validated sterilisation procedure according ISO 17665 - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure. - Control the efficiency (packaging integrity, no humidity, colour change of sterilisation indicators, physico-chemical integrators, digital records of cycles parameters). - Traceability of procedure records	X	X
9.	Storage	- Keep devices in sterilization packaging in a dry and clean environment	- Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).	X	X

B. Filling material

1.	Operation	Operating mode	Warning
	Disinfection	- Immerse the obturation devices in NaOCl (2,5 % at least) during 5 mn at ambient temperature.	- Do not use disinfecting solutions containing Phenol or any products which are not compatible with the treated filling material (See general recommendation).

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