 Stick Tech Ltd

Lemminkäisenkatu 46  
P.O. Box 114  
FI-20521 TURKU, FINLAND  
Tel: + 358 2 4808 2500  
Fax: +358 2 241 0032  
www.sticktech.com

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
## everStickORTHO

Fibre reinforcement for aesthetic orthodontic retainers

Fibre type: Silanated E-glass fibre impregnated with bis-GMA and PMMA  
Form: Unidirectional fibre bundle  
Diameter: ~ 0.75 mm

Caution: Federal law restricts this device to sale by or on the order of a licenced dental professional.

**everStick®**  
GIANT OF FIBRES


# everStick®ORTHO

Fibre reinforcement for aesthetic orthodontic retainers

PA10266 updated 2011-06

everStick®  
**ORTHO** Orthodontic retainer

---

---

---

---

---

---

---

---

---

---

<b>EN</b>	ENGLISH .....	3
<b>FI</b>	SUOMI .....	7
<b>ES</b>	ESPAÑOL.....	11
<b>NL</b>	NEDERLANDS .....	15

<b>IT</b>	ITALIANO.....	21
<b>NO</b>	NORSK.....	27
<b>PL</b>	POLSKI .....	31
<b>FR</b>	FRANÇAIS.....	37

<b>SV</b>	SVENSKA .....	41
<b>DE</b>	DEUTSCH .....	45
<b>DA</b>	DANSK.....	51
<b>RU</b>	РУССКИЙ.....	55

# everStick® **ORTHO**

everStick®ORTHO

– glass fibre reinforcement for orthodontics:  
lingual retainer on the lower jaw

## **WHAT IS everStick®ORTHO?**

everStickORTHO fibre reinforcement is made of glass fibres and a porous polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which

facilitates handling of the fibre bundle. The fibre bundle is flexible and sticky, which allows it to easily bond tightly to teeth.

As everStickORTHO fibre does not possess any memory, when positioned and polymerised it can be fixed passively on to the teeth. everStickORTHO fibre reinforcement's indication is the retention phase after active orthodontic treatment.

## **INSTRUCTIONS FOR USE:**

**CAUTION:** everStickORTHO fibre bundles must be coated with a thin (0.5 mm) layer of composite, including the interproximal spaces, when fixing it on to the teeth. Spot bonding with composite, which is

a commonly used method when using archwires in orthodontic therapy, does not provide sufficient bonding between the fibre retainer and surface of the tooth. The fibre retainer must be positioned as incisally as possible to minimise the forces that might otherwise loosen it. Also check that the retainer is not in occlusal contact if you place it, for example, on the palatal side of the maxilla.

### **1. Measuring and cutting the fibre**

Using, for example, a piece of dental floss, measure the length of the fibre required for making the everStickORTHO retainer. Talc-free gloves are recommended when handling everStickORTHO fibres. Open the foil pouch and draw the silicone wrapped fibre bundle partly out. Using sharp scissors

cut the desired length. Shield the cut piece of fibre from light by placing it under a cover during preparation of the tooth surfaces. Close the foil pouch with its sticker. Store the pouch in a refrigerator (at a temperature of +2 ... +8°C, +35 ... +46°F) between uses.

## **2. Cleaning the tooth surfaces**

The entire length of the retainer must be attached to the tooth surfaces. Clean the tooth surfaces with pumice and water, and air-dry the area. Alternatively, you can sandblast the surfaces with a microetcher for about 5 seconds per tooth. Sandblasting increases the bond strength. After sandblasting, the surfaces must be rinsed with water and air-dried. Place wedges in the interproximal spaces, if possible, so that they do not

become filled with composite. If you are working without wedges, note that the interproximal spaces must remain free of composite afterwards – see item 5.

## **3. Etching the tooth surfaces**

In the area where the retainer and composite is to be placed, etch the tooth surfaces and interproximal spaces thoroughly with ortho-phosphoric acid. The correct etching time is about 45-60 seconds. Rinse and dry the tooth surfaces carefully after etching. Ensure that the bonding surface is dry before applying any resin.

## **4. Treating the tooth surfaces with resin**

Apply a thin layer of bonding agent to the tooth surfaces in the area of the

retainer. Air dry the bonding agent to a thin layer. Light cure the bonding agent as described by the manufacturer.

## **5. Bonding**

Apply a thin layer of flowable composite (for example, StickFlow or a similar orthodontic composite) on to the tooth surfaces in the area of the retainer, two teeth at a time. Carefully cover the bonding area with a thin layer (about 0.5 mm) of composite, including the interproximal spaces. Leave enough room for cleaning the interdental spaces. Do not cure the composite during this phase.

## **6. Positioning and light-curing the fibre**

Remove the white cover paper and use tweezers to pick the fibre bundle up from the

silicone groove. Remove any silicone granules from the fibre bundle. Place the fibre bundle within the flowable composite on the tooth. Note that you can coat the fibre bundle with an additional thin layer (0.5 mm) of composite if necessary. This can be applied with a brush. Then the need for finishing will be minimal.

Attach the retainer as incisally as possible. Check that it will not be in occlusion. You can first place one end of the fibre bundle in position by pressing it into the flowable composite with an instrument (e.g. StickStepper). Alternatively, you can position the whole retainer at one go using the special silicone instrument (RefixD). Hold the fibre in place, one tooth at a time, using a curing light for about 5-10 seconds. Shield the uncured

fibre from the light with a wide StickStepper instrument. It is recommended to point the light-curer away from the uncured fibre bundle. Press the fibre into the interproximal spaces as well. Keep the fibre bundle as round as possible, particularly in the interproximal spaces, so that the cleaning areas will not be covered with fibres and composite.

#### **7.-8. Coating and finishing the retainer**

After pre-curing, coat the whole retainer with a thin layer of composite. Then light-cure the whole retainer for 40 seconds per tooth. Check and adjust the occlusion. Do not cut the fibre when finishing and polishing the retainer.

**STORING:** everStick products should always be stored in a refrigerator

(+2 ... +8°C, +35 ... +46°F). In addition, the products should be protected from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

**NOTE:** Stick Tech's products should be used

clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing irritation-causing fibers.

The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours.

StickStepper, StickCarrier hand instruments and RefixD, RefixL silicone instruments must be sterilized before use.

**WARNING:** Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous

membrane, or eyes. Unpolymerized everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick before waste disposal.

**WARRANTY:** Stick Tech Ltd replaces defective products. Stick Tech does not assume liability for any damage or loss, whether direct or consequential, caused by improper use of the product or use contrary to the accompanying instructions. Before starting to use the product, the user is responsible for assessing the suitability of the product for the intended purpose. The user bears all risk and liability for any consequences due to use contrary to the instructions

or in inappropriate circumstances.

Federal law restricts this device to sale by or on the order of a licensed dental professional.

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Symbol	PRODUCTS	RECYCLOR DISPOSAL OF PACKING MATERIALS
	StickCARRIER / StickSTEPPER plastic cover	Can be recycled or burned in waste incinerator.
	everStickPOST INTRO plastic tray	Can be recycled but not suitable for incineration.
	Plastic bags, everStick STARTER KIT foam, StickFLOW syringe and StickRESIN bottle	Can be recycled or burned at home or in waste incinerator.
	Stick, StickNET and everStickNET primary packages	Can be recycled or burned in waste incinerator.
	everStickSTARTER KIT acrylic box	Can be recycled or burned in waste incinerator.
	Cardboard cases or printed user's instruction	Can be recycled or burned at home or in waste incinerator.
	Silicones and aluminium foil bags	Can be disposed with mixed waste according to official regulation.

The national, local and European regulation of package and package waste directive 94/62/EEC must be taken into consideration during disposal.

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---