## Essix C+® Plastic

(Instructions for use on Essix® Vacuum Machines only)

### Intended Use:

Essix C+® Plastic is thermoformed to fabricate intra-oral appliances such as retainers, aligners for minor tooth movement, and bite planes. Essix C+® Plastic is especially suited for use with bruxers/grinders.

### Indications for Use:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

**Contraindications:** None known.

### ⚠ Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

### Precautions:

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only.

**Adverse Reactions:** Allergic reactions.

### Step by Step Instructions:

1. All models should be trimmed to 3/4" high. Full arch models should have the palate trimmed away.
2. **IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
3. Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
4. Place the plastic into the frame of the vacuum machine and the model on the base plate.
5. After pre-heating, swing the heating element over the plastic.
6. Observe the plastic as it heats.
7. During heating, Essix C+® Plastic rises toward the heat source, then becomes wavy as it lowers from the raised position.
8. The waves will start to flatten out into its original flat position – this is the optimal thermoforming time. Do not let Essix C+® Plastic sag.
9. Move the heating element away from the plastic.
10. Turn on the vacuum, then lower the plastic over the model.
11. Push into the plastic at each undercut with the Essix Accentuator® Tool.
12. Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it "shrinks" the material around the cast, realizing a tighter fit.
13. Trim and finish appliance as desired.

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic "folds" when thermoformed, reduce heating time until material does not fold.


Thickness	Vacuum Heat Time*	Biostar™ Code*	Essix® SelectVac®	Druformat
1.0mm (.040")	60 seconds	163	2:30	1:30

\*after preheating element

Manufacturer:

**DENTSPLY International**  
**DENTSPLY Raintree Essix**  
7290 26th Court East  
Sarasota, FL 34243 USA  
Tel: 800.883.8733  
www.essix.com

Emergo Europe  
Molenstraat 15,  
2513 BH, The Hague  
The Netherlands  
Tel: (+31) 70 345 8570  
Fax: (+31) 70 346 7299

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