KaVo DIAGNOcam

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### Estudio 1

Am Dent Assoc. 2008 Oct;139(10):1374-81  
The correlation of DIFOTI to clinical and radiographic images in Class II carious lesions.  
Bin-Shuwaish M, Yaman P, Dennison J, Neiva G.  
Department of Restorative Dental Sciences, College of Dentistry, King Saud University, Riyadh, Saudi Arabia.

#### Abstract

**BACKGROUND:** The authors conducted a study to evaluate the correlation between digital imaging fiber-optic transillumination (DIFOTI) (KaVo Dental, Lake Zurich, Ill.) and clinical and radiographic images in estimating the true clinical axial extension of Class II carious lesions. **METHODS:** The authors examined 51 Class II carious lesions visually, imaged them by means of DIFOTI and radiographed them with D-speed film and a complementary metal oxide silicon (CMOS)-based digital radiographic sensor. They validated axial extension of the lesions clinically. They compared the clinical and radiographic depths of the carious lesion with the size of the lesion on the DIFOTI images. **RESULTS:** The authors detected 84 percent of the lesions with DIFOTI, and 82 percent showed a visible dark shadow under the marginal ridge when examined clinically. DIFOTI correlated significantly with the clinical depth of decay (Pearson  $r = 0.43189$ ). The combination of a CMOS digital sensor and DIFOTI ( $R^2 = 0.7210$ ) provided readings closer to the clinical measures than did the combination of D-speed film and DIFOTI ( $R^2 = 0.6215$ ). **CONCLUSIONS:** DIFOTI images correlated with clinical depth, especially in smaller lesions, and improved the estimation of lesion size when used in conjunction with the CMOS digital sensor and D-speed images. **CLINICAL IMPLICATIONS:** Using radiographs in combination with DIFOTI images could help clinicians determine the presence and, to some extent, the size of proximal caries, especially in smaller lesions. **PMID:** 18832273 [PubMed - indexed for MEDLINE] Free Article Am Dent Assoc. 2008 Oct;139(10):1374-81

### Estudio 2

Caries Res. 1997;31(2):103-10. Assessment of dental caries with Digital Imaging Fiber-Optic Transillumination (DIFOTI): in vitro study. Schneiderman A, Elbaum M, Shultz T, Keem S, Greenebaum M, Driller J. Department of Oral Pathology, Biology, New Jersey Dental School, University of Medicine and Dentistry of New Jersey, Newark 07103, USA.

#### Abstract

This paper describes Digital Imaging Fiber-Optic Transillumination (DIFOTI™), a new method for the reliable detection of dental caries. Images of teeth obtained through visible-light, fiber-optic transillumination (FOTI) are acquired with a digital CCD camera, and sent to a computer for analysis with dedicated algorithms. The algorithms were developed to facilitate the location and diagnosis of carious lesions by the operator in real time, and provide quantitative characterization for monitoring of the lesions. The DIFOTI method has been tested by imaging teeth in vitro. The results suggest the superior sensitivity of DIFOTI for detection of approximal, occlusal and smooth-surface caries vis-à-vis radiological imaging. **PMID:** 9118181 [PubMed - indexed for MEDLINE]

### Estudio 3

From the Divisions of Cariology and Endodontology, Department of Dental Medicine Karolinska Institutet, Stockholm, Sweden

Exploring the boundaries of caries detection  
Álfheiður Ástvaldsdóttir  
Two advanced methods evaluated

#### Abstract

**BACKGROUND:** Caries detection methods require thorough validation. This should include studies which clarify what characteristics of the caries lesions are being measured, the limitations of the method and comparison of performance with conventional caries detection methods. The outcome of validation tests has important clinical implications, such as interpretation of the data at the cut-off points used by the clinician to differentiate between lesions requiring invasive and non-invasive intervention. **CONCLUSIONS:** The LF method can be useful for occlusal caries detection and quantification, under certain conditions: a) no general cut-off threshold can be recommended due to inconsistency between instruments. b) the method does not give information on demineralization or specific information on bacterial content of lesions, but rather responds to the synergistic effect of the caries process. The results of the in vitro investigations suggest that the DIFOTI method may be of value for caries detection and quantification on both approximal and occlusal surfaces. The method shows superior performance to both film and digital radiography, especially for detection of early caries lesions on approximal surfaces. Thus the method shows promise as a means of monitoring early caries lesions and warrants further investigation.

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