Essix A+[®]

Plastic

DIRECTIONS FOR USE
INSTRUCCIONES DE USO
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG

Manufacturer:
Dentsply International Raintree Essix
7290 26th Court East
Sarasota, FL 34243 USA
Tel: 800.883.8733
www.essix.com

Emergo Europe
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

Biostar and Drufomat are not trademarks of
Dentsply Sirona. ©2017 Dentsply Sirona

#RT-DFU-A+ RTE-061-18 Rev.8 3/18

(Instructions for use on Essix® Vacuum Machines only)*

Intended Use:

Essix A+® Plastic is thermoformed to fabricate intra-oral appliances such as bite planes, splints, and snoring appliances. Essix A+® Plastic is suitable for use with non-bruxers.

Indications for Use:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

Contraindications:

This product is contraindicated for patients with a history of allergic reactions to plastics.

⚠ Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

Precautions:

Essix A+® Plastic is not recommended for bruxers.

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only. Risks of reuse include infection and cross-contamination to patient and healthcare provider as well as other wear-related irregularities.

Adverse Reactions: Allergic reactions.

Step by Step Instructions:

- Essix A+® Plastic has a protective film on both sides. REMOVE FILM FROM BOTH SIDES BEFORE THERMOFORMING.
- All models should be trimmed to $\frac{3}{4}$ " high. Full arch models should have the palate trimmed away.
- IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
- Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
- Place the plastic into the frame of the vacuum machine and the model on the base plate.
- After pre-heating, swing the heating element over the plastic.
- Observe the plastic as it heats.
- Essix A+® Plastic should be heated until it sags approximately $\frac{1}{2}$ ".
- Move the heating element away from the plastic.
- Turn on the vacuum and lower the plastic over the model.
- Push the plastic into each undercut with the Essix Accentuator® Tool.
- Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it "shrinks" the material around the cast, realizing a tighter fit.
- Let vacuum continue to run until plastic is cool.
- Trim and finish the appliance as desired.

*Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic "folds" when thermoformed, reduce heating time until material does not fold.

Thickness	Vacuum Heat Time**	Biostar™ Code**	Essix® SelectVac®	Drufomat
0.5mm (.020")	30 seconds	122	1:20	1:00
0.75mm (.030")	35 seconds	132	1:30	1:00
1.0mm (.040")	40 seconds	142	1:35	1:05
1.5mm (.060")	50 seconds	164	2:00	1:10
2.0mm (.080")	70 seconds	205	2:05	1:30
3.0mm (.120")	80 seconds	226	3:20	2:10

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