FICHA TÉCNICA

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TECHNICAL DATA SHEET

PROCLINIC SELF-SEALING STERILIZATION POUCHES		
REFERENCES	<p>79930: 75 x 250 mm 79931: 60 x 100 mm 79932: 90 x 250 mm 79933: 140 x 260 mm 79947: 300 x 450 mm 79954: 90 x 165 mm</p>	
MANUFACTURER	<p>E.C.S. Srl Via Como , 71 23883 Brivio (LC), Italia Tel.: +39.0341.681602 Fax: +39.0341.6054 Email: info @ecssrl.com – www.ecssrl.com</p>	
CLASSIFICATION	<p>Medical Device, Class I, according to Regulation (EU) 2017/745 on medical devices.</p>	
DESCRIPTION	<p>Self-sealing sterilization pouches for single use</p>	
COMPOSITION	<p>Made of medical use laminated transparent film (PET and CPP) and medical paper.</p>	
SPECIFICACIONES	<p>Sterilization with steam: temperature between 121°C-134°C; the indicator colour changes from pink to brown. Sterilization with ethylene oxide: temperature between 30°C-54°C; the indicator colour changes from light blue to yellow/brown.</p>	
PACKAGING	<p>Box with 200 pouches</p>	
SHELF LIFE	<p>5 years from the date of manufacture</p>	
INTENDED USE	<p>Sterilization of medical devices and materials for medical use which are especially suitable for steam or ethylene oxide gas.</p>	
HANDLING AND STORAGE	<p>Store in a cool dry place, away from direct sunlight. Protect from moisture.</p>	

<p>SAFETY INFORMATION</p>	<ul style="list-style-type: none"> • Before starting sterilization process, check that the pouch is filled a maximum of 3/4 of its volume, so that sealings are not obstructed at the edges. • Cover the end and the cutting surfaces of instruments with appropriate protection. • Put the instruments with the handle at the opening of the pouch. • Do not use solvent substances on the pouch. • Possible printings have to be written on a label and sticked on paper or film side; please do not write directly on the pouch, and anyway, use the maximum care to avoid contamination problem during use. • Carefully check the weight of the instruments to be sterilized in order to avoid risks of breaking.
<p>PRESERVATION OF STERILITY</p>	<p>It's essential to grant the correct condition for storage and use of the material after the sterilization process.</p> <p>Preservation of sterility depends on type of packaging and storage conditions. The sterilized material should be stored in closed, clean and dry cupboards, protected from powder, with limited access to closed windows and doors, at a temperature between 18-22°C and humidity between 35-50%.</p> <p>The material should be handled as less as possible and with care. Packages broken, damaged, opened, fallen down or wet must be considered as non-sterile, and must be repacked and sterilized again.</p> <p>Sterility duration for steam depends on the correct and strict observance of the norms and the technical instructions for sterilization.</p> <p>In this context, we suggest a duration of sterility for steam in single pouch (usually used for non-invasive used on laboratory and stay in bed operation unit) of 30 days, and in double pouch (usually used with "no touch" technique for medical devices used for invasive procedures and for medical devices used in Surgery Unit) a duration of 60 days, on condition that in both cases, the conditions suitable for conservation are respected.</p>