Prior to use, carefully read the instructions for use.

EN

GC Fuji® II LC CAPSULE

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE IN CAPSULES
For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Restoration of Class III, V and limited Class I cavities.
- Restoration of primary teeth.
- Core build up.
- Cases where a radiopaque restoration is required.
- Geriatric applications.
- As a base or liner.

CONTRAINDICATIONS

- Pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

DIRECTIONS FOR USE

| | |
|-----------------------------|-------------|
| Powder / Liquid Ratio (g/g) | 0,33 / 0,10 |
| Mixing Time (sec.) | 10" |
| Working Time (min., sec.) | 3'15" |
| Light Curing Time (sec.) | 20" |
| Depth of Cure (A2) (mm) | 1.8 |

Test conditions : Temperature (23 +/-1°C)
Relative humidity (50 +/-10%)

1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Wash and dry but DO NOT DESICCAT. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer (Fig. 1).
- CAPSULE ACTIVATION AND MIXING**
 - Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
 - To activate the capsule, push the plunger until it flush with main body (Fig. 3).
 - Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 4). The capsule is now activated. Note : The capsule should be activated just before mixing and used immediately.
 - Immediately remove the capsule and set it into a mixer (an amalgamator) and mix for 10 seconds at high speed (+/-4.000 RPM) (Fig. 5).

3. RESTORATIVE TECHNIQUE

- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- Make two clicks to prime the capsule then syringe (Fig. 6).
The working time is 3 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time. Note :
 - To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.
- Remove surface moisture but DO NOT DESICCAT.
- Extrude cement directly into preparation (Fig.7). Avoid air bubbles.
- Form the contour and place a matrix if required.
- Light-cure for 20 seconds using a suitable visible light curing device (470nm wavelength) (Fig. 8). Place light source as closely as possible to the cement surface. Note : For cavities deeper than 1.8 mm, use a layering technique.

4. FINISHING

- Finish under water spray using superfine diamond bur, silicone point and polishing strips (Fig. 9).
Note :
Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.)(Fig. 10).

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F).
(Shelf life : 2 years from date of manufacture)

A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, D2 (11 shades based on Vita® shades).

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 pcs each of A4, B2, B3, C2; 5 pcs each of B4, D2).
*Average contents per capsule : 0,33g powder and 0,10g (0,085mL) liquid.
*Minimum net volume of mixed cement per capsule : 0,10mL.
Option : GC Capsule Applier (1 piece).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>
They can also be obtained from your supplier.

Last revised : 02/2018

0086

GC

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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