# GRADIA® PLUS

MANUFACTURED by  
GC DENTAL PRODUCTS CORP.  
2-285 Torrimatsu-cho, Kasugai, Aichi 486-0844, Japan

RESPONSIBLE MANUFACTURER IN CANADA:  
GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED by:  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

 EU: GC EUROPE N.V.  
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,  
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900  
www.gcamerica.com

GC SOUTH AMERICA  
RUA HELIODORA, 399,  
SANTANA - 02022-051  
SÃO PAULO, SP - BRASIL  
TEL: +55-11-2925.0965  
CNPJ: 08.279.999/0001-61  
RESPONSÁVEL TÉCNICA  
MAYARA DE SANTIS RIBEIRO  
CRO/SP: 105.982

GC ASIA DENTAL PTE. LTD.  
11, Tampines Concourse, #03-05 Singapore 528729  
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.  
1753 Botany Rd, Banksmeadow, NSW 2019, Australia  
TEL: +61 2 9301 8200

PRINTED IN BELGIUM

- EN MODULAR COMPOSITE SYSTEM FOR INDIRECT RESTORATIONS
- DE MODULARES KOMPOSITSYSTEM FÜR INDIREKTE RESTAURATIONEN
- FR SYSTEME COMPOSITE MODULAIRE POUR RESTAURATIONS INDIRECTES
- IT SISTEMA MODULARE DI COMPOSITO PER RESTAURI INDIRETTI
- ES SISTEMA DE COMPOSITE MODULAR PARA RESTAURACIONES INDIRECTAS
- NL MODULAIR COMPOSITESTYSTEEM VOOR INDIRECTE RESTAURATIES
- DA MODULOPBYGGET KOMPOSITSYSTEM TIL INDIREKTE RESTAURERINGER
- SV MODULUPPBYGGT KOMPOSITSYSTEM FÖR INDIREKTA RESTAURERINGAR
- PT SISTEMA DE COMPÓSITO MODULAR PARA RESTAURAÇÕES INDIRETAS
- EL ΣΥΣΤΗΜΑ ΣΥΝΘΕΤΗΣ ΡΗΤΙΝΗΣ ΓΙΑ ΕΜΜΕΣΕΣ ΑΠΟΚΑΤΑΣΤΑΣΕΙΣ
- NO MODULÆRT KOMPOSITTSYSTEM FOR INDIREKTE RESTAURERINGER
- FI MODULAARINEN YHDISTELMÄMUOVIJÄRJESTELMÄ EPÄSUORIA TÖITÄ VARTEN
- BG МОДУЛНА КОМПОЗИТНА СИСТЕМА ЗА ИНДИРЕКТНИ ВЪЗСТАНОВЯВАНИЯ
- CS MODULÁRNÍ KOMPOZITNÍ SYSTÉM PRO NEPŘÍMÉ DOSTAVBY
- HR MODULARNI KOMPOZITNI SUSTAV ZA INDIREKTNE RADOVE
- HU MODULÁRIS KOMPOZIT RENDSZER AZ INDIREKT HELYREÁLLÍTÁSOKHOZ
- PL MODUŁOWY SYSTEM KOMPOZYTOWY DO POŚREDNICH UZUPEŁNIENIŃ PROTETYCZNYCH
- RO SISTEM COMPOZIT MODULAR PENTRU RESTAURĂRI INDIRECTE
- RU МОДУЛЬНАЯ КОМПОЗИТНАЯ СИСТЕМА ДЛЯ НЕПРЯМЫХ РЕСТАВРАЦИЙ
- SK MODULÁRNY KOMPOZITNÝ SYSTÉM PRE NEPRIAME NÁHRADY
- SL MODULARNI KOMPOZITNI SISTEM ZA INDIREKTNE OBNOVE
- SR MODULARNI KOMPOZITNI SISTEM ZA INDIREKTNE RADOVE
- UK МОДУЛЬНА КОМПОЗИТНА СИСТЕМА ДЛЯ НЕПРЯМИХ РЕСТАВРАЦІЙ
- TR İNDİREK RESTORASYONLAR İÇİN MODÜLER KOMPOZİT SİSTEMİ
- LV MODULĀRA KOMPOZĪTA SISTĒMA NETIEŠĀM RESTAURĀCIJĀM
- LT KOMPOZITINĖ SISTEMA NETIESIOGINIŲ RESTAURACIJŲ GAMYBAI
- ET MODULAARNE KOMPOSIIDI SÜSTEEM INDIREKTSETELE RESTAURATSIOONIDELE

951941  
X7504

Prior to use, carefully read the instructions.

## GRADIA® PLUS MODULAR COMPOSITE SYSTEM FOR INDIRECT RESTORATIONS

For use only by a dental professional in the recommended indications.

### INDICATIONS

1. Crown and Bridgework (with metal backing)
2. Metal free restorations: Jacket Crowns, Inlays, Onlays and Laminated Veneers
3. Implant superstructures
4. Reproduction of gum tissue for crown restorations

### CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer or alcohol.

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###### 1. Production of Metal Framework

- a) Waxing-up and Casting  
Wax-up framework conventionally, occlusal surface should have a minimum of 1 mm thickness. Apply retention beads to the facing surface. Cast in the usual manner and remove retainers in 0.5 mm width around the marginal veneer surface.
- b) Sandblasting  
Sandblast the metal frame with 50 micron aluminum oxide. Clean and dry with air gun.
- c) Priming  
Immediately apply METAL PRIMER II / METAL PRIMER Z to retention areas. Allow to dry.

###### 2. Build up Veneer

- a) Apply and Light cure the Opaque:  
Immediately apply a layer of OPAQUE over retention area and light cure for 1 minute. Repeat this process until the metal color is masked out.  
Note: For irradiation time by the GC's visible light curing unit, refer to the chart given below.
- b) Forming Dentin:  
Apply Dentin shades\*, building up to the desired thickness. Pre-cure. (In some cases, adjust the shade by Lustre Paint and light cure for 1 minute.)  
*\*A processing aid like MODELING LIQUID can be used. Use a few drops, a brush or adequate instrument to adapt the pastes.*
- c) Cervical application:  
Progressively thinly apply Lustre Paint (A) in 2-3 mm width around the cervical area as required and light cure for 1 minute.
- d) Forming Enamel:  
Apply Enamel shades progressively from the incisal edge. Light cure for 3 minutes.
- e) Adjust and Characterize.  
Adjust contour and characterize surface with diamond and carbide burs. (At this time, scrape off all air inhibition layers on the surface.)

### **3. Glazing (Coating and Polishing)**

#### **3-1 Coating method**

- a) Roughening the surface:  
Roughen the resin surface using a carbide bur or sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion.  
Dry and clean with air gun.
- b) Priming:  
Immediately apply a silane coupling agent (e.g. CERAMIC PRIMER II) to the resin surface and dry.
- c) Coloring, Coating and Final light curing:  
Apply a thin layer (<0.1mm) of Lustre Paint over the resin surface and final-cure.  
To prevent contamination, wash the brush with Lustre Paint Diluting Liquid every time.

#### **3-2 Polishing method**

Use proper instruments and burs to polish and finish the prosthesis, and to confirm, remove any agent or paste for polish and finish in proper way and see the surface has proper shine.

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### **1. Preparation of the working model**

- a) Prepare the die stone in usual manner.
- b) Apply DIE-HARDNER to the die stone surface.
- c) Block out undercuts with wax.
- d) Apply a layer of paraffin wax as die spacer as required.
- e) Apply a thin layer of a resin composite separator.

### **2. Special pre-cure treatment**

- a) In case of a metal core, apply first a thin layer of OPAQUE. Light cure for 1 minute.
- b) In case of a vital abutment tooth, apply first a thin layer of (LB-) Base OD and Pre-cure.

### **3. Build up of Crown**

For following build up on the surface, see instructions BUILDING UP PROCEDURE A. (Section 2. b-e) above.

### **4. Glazing**

See instructions BUILDING UP PROCEDURE A. (Section 3.) above.

### **5. Cleaning, treating restoration for bonding**

Roughen the inner surface of Jacket crown using alumina sandblast (1.5 kg/cm<sup>2</sup>) for obtaining the mechanical adhesion.  
And then clean with air gun.

### **6. Cementation (Chair side)**

- a) Priming  
Apply a thin layer of a silane coupling agent (e.g. CERAMIC PRIMER II) to bonding surface.
- b) Cementation  
Use a dental adhesive resin cement (e.g. G-CEM LinkForce).

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### **1. Preparation of the working model**

Follow the instructions Preparation of the working model (B. Section 1. a-e) above.

### **2. Build up of Inlay**

- a) Build up (LB-) Inlay TD as a dentin color in order to retain the natural tooth color. Pre-cure.  
Note: In the case of a stained tooth, it is recommended to first apply OPAQUE to the cavity floor. Light cure for 1 minute.
- b) Build up (LB-) Inlay E, shape and contour. Final-cure.
- c) Continue BUILDING UP PROCEDURE A. (Section 2. e) above.

### **3. Glazing**

See instructions BUILDING UP PROCEDURE A. (Section 3.) above.

## II. ADDITIONAL APPLICATION AND REPAIR

### A. LABSIDE

1. Roughen the resin surface using a carbide bur or sandblast (1.5 kg/cm<sup>2</sup>).
2. Immediately apply a silane coupling agent (e.g. CERAMIC PRIMER II) to the resin surface and dry.
3. Apply Enamel shades, etc. as required.

### B. CHAIRSIDE

#### MASKING OF EXPOSED METAL SURFACES WHEN REPAIRING INDIRECT RESTORATIONS INTRAORALLY

1. As for PREPARATION of RESTORATION, refer to G-Premio BOND IFU-II.1.  
Note: Use a rubber dam for isolation.
2. Apply GRADIA PLUS OPAQUE to the prepared metal surface only. The layer applied should have a thickness of maximum 0.2 mm. Light-cure GRADIA PLUS OPAQUE for 20 seconds.
3. If the surface is still not sufficiently masked, apply a second layer of OPAQUE and light-cure.

#### Irradiation time

Curing unit	Halogen / LED (700 mW/cm <sup>2</sup> ) & High power LED (more than 1200 mW/cm <sup>2</sup> )
OPAQUE	20 seconds

4. Repair the restoration by applying a chairside composite material (e.g. G-ænial & Essentia) in layers. Light-cure each layer according to the manufacturer's instructions of the respective product
5. Finish and polish using standard techniques.

#### Irradiation time and curing unit

Curing unit	LABOLIGHT DUO		STEPLIGHT SL-I	LABOLIGHT LV-III,II
	Step-mode**	Full-mode	Pre-Cure**	Final cure
OPAQUE	-	1 min.	-	1 min.
PASTE HB, PASTE LB, GUM SHADES LB, GUM SHADES HB	10 sec.*	3 min.	10 sec.*	3 min.
LUSTRE PAINT ***	10 sec.*	90 sec.	10 sec.*	5 min.

\* For one surface of a single crown.

\*\* Distance from light source : 3cm.

\*\*\* Thickness: 0.1mm or less.

**Configuration table: Opaque mixing ratio**

OA1	OA2	OA3	OA3.5	OA4
O-Base:OA 3 : 1	O-BASE:OA 1 : 1	O-BASE:OA 1 : 2	O-BASE:OA 1 : 3	OA
OB1	OB2	OB3		OB4
O-BASE:OB 3 : 1	O-BASE:OB 1 : 1	O-BASE:OB 1 : 3		OB
OC1	OC2	OC3		OC4
O-BASE:OC 3 : 1	O-BASE:OC 1 : 1	O-BASE:OC 1 : 3		OC
	OD2	OD3		OD4
	O-BASE:OD 1 : 1	OA:OD 1 : 2		OD

**Corresponding table (Enamel, Dentin and Lustre Paint)**

Shade	A1	A2	A3	A3.5	A4
Enamel (Paste HB)	EL	EL	ED	ED	ED
Dentin (Paste HB)	DA1	DA2	DA3	DA3.5	DA3.5
Cervical (Lustre Paint)	LP-CLF	LP-CLF	LP-CLF	LP-CLF	LP-A
Shade	B1	B2	B3		B4
Enamel (Paste HB)	EL	EL	ED		ED
Dentin (Paste HB)	DB1	DB1	DB3		DB3
Cervical (Lustre Paint)	LP-CLF	LP-B	LP-CLF		LP-B
Shade	C1	C2	C3		C4
Enamel (Paste HB)	EL	ED	ED		ED
Dentin (Paste HB)	DB1	DB1	DC3		DC3
Cervical (Lustre Paint)	LP-C	LP-C	LP-CLF		LP-C
Shade	D2		D3		D4
Enamel (Paste HB)	EL		ED		ED
Dentin (Paste HB)	DD2		DD2		DD2
Cervical (Lustre Paint)	LP-CLF		LP-A		LP-D

	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	BW	
Opaque*	OA1	OA2	OA3	OA3.5	OA4	OB1	OB2	OB3	OB4	OC1	OC2	OC3	OC4	OD2	OD3	OD4	O-BASE	
Dentin	(HB-)DA1	(HB-)DA2	(HB-)DA3	(HB-)DA3.5		(HB-)DB1		(HB-)DB3		(HB-)DB1		(HB-)DC3			(HB-)DD2		(HB-)DW	
Opacious Dentin	(HB-)ODA					(HB-)ODB			(HB-)ODC				(HB-)ODD		(HB-)ODW			
Cervical (Lustre Paint)	(LP-)CLF, (LP-)A					(LP-)CLF, (LP-)B			(LP-)CLF, (LP-)C				(LP-)CLF, (LP-)D		(LP-)CLF, (LP-)A			
Enamel	(HB-)EL	(HB-)ED				(HB-)EL	(HB-)ED	(HB-)EL	(HB-)ED	(HB-)EL	(HB-)ED	(HB-)EL	(HB-)ED	(HB-)EL	(HB-)ED	(HB-)EL	(HB-)ED	(HB-)EL
Effect 1 (Dentin shade)	(LB-)Base D, (LB-)Base OD, (LB-)DW, (LB-)ODW																	
Effect 2 (Enamel shade)	(HB-)PE, (LB-)Base E, (LB-)Inlay E																	
Effect 3 (Characterize)	(LB-)Base Opal, (LB-)Orange, (LB-)Red, (LB-)Yellow, (LB-)Blue, (LB-)Grey, (LB-)Milky																	
Translucent	(HB-)CLF, (LB-)Base CLF,																	
Cervical Translucent	(LB-)Inlay TD																	
Gum shade	GO-1, GO-2, GLB-1, GLB-2, GLB-3, GLB-CL																	
One body	(LB-)A					(LB-)B				(LB-)C				(LB-)D		(LB-)W		
Effect 4 (Lustre Paint)	(LP-)Blue, (LP-)Cream, (LP-)Grey, (LP-)Lavender, (LP-)CL, (GLP-)Bright red, (GLP-)Violet																	

\*cf. Configuration table: Opaque mixing ratio

## PACKAGES

### 1. Set Packages

- **LAYER SET (19 shades)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 shades)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 shades)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 shades)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 22 Blue
- **ONE BODY SET (11 shades)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **ACCESSORY KIT**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

### 2. Refill Packages

- **OPAQUE, 2.0mL syringe (in 13 shades)**  
O-Base, OA, OB, OC, OD, OA1, OA2, OA3, OA3.5, OB1, OB3, OC3, OD2, Dispensing Tip Needle Type, Light Protective Cover
- **PASTE HB, 3.3mL syringe (in 18 shades)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-DW, HB-ODA, HB-ODB, HB-ODC, HB-ODD, HB-ODW, HB-EL, HB-ED, HB-CLF, HB-PE
- **PASTE LB, 2.0mL syringe (in 15 shades)**  
LB-Base D, LB-Base E, LB-Base CLF, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Red, LB-Yellow, LB-Blue, LB-Grey, LB-Milky, LB-Inlay E, LB-Inlay TD, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **ONE BODY, 2.0mL syringe (in 5 shades)**  
LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Opaque, 2.0mL syringe (in 2 shades)**  
GO-1, GO-2, Dispensing Tip Needle Type, Light Protective Cover
- **GUM Shades HB, 3.3mL syringe (in 4 shades)**  
GHB-1, GHB-2, GHB-3, GHB-CL
- **GUM Shades LB, 2.0mL syringe (in 4 shades)**  
GLB-1, GLB-2, GLB-3, GLB-CL, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Lustre Paint, 0.8mL syringe (in 2 shades)**  
GLP-Bright Red, GLP-Violet, Dispensing Tip Needle Type Small, Light Protective Cover
- **LUSTRE PAINT, 0.8mL syringe (in 10 shades)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-Blue, LP-Cream, LP-Grey, LP-Lavender, LP-CL, Dispensing Tip Needle Type Small, Light Protective Cover

### 3. Accessories

- MODELING LIQUID 3mL
- Lustre Paint DILUTING LIQUID 3mL
- AIR BARRIER 10mL
- SEPARATOR 5mL
- METAL PRIMER II 5mL / METAL PRIMER Z 5mL
- CERAMIC PRIMER II 3mL
- DIE-HARDNER 5mL
- DIAPOLISHER PASTE 2g

## NOTES

- Please review this Instruction for Use carefully before using the product. After reviewing the Instructions for Use, keep them in a convenient place for quick reference if necessary. If this Instruction for use is lost, contact the distributor or supplier shown on the package.
- This product must only be used by dentists or dental technicians.
- To prevent troubles with operating or storing, the product should be used according to the instructions specified in this Instruction for Use. Do not use this product for contraindications given in this Instruction for Use.
- Assure that the crown or inlay is sufficiently thick, using the methods for preparing abutments and cavities described in this manual.
- Wear gloves or take other appropriate measures when using the product.
- After dispensing, use immediately and replace caps on syringe or bottle.
- Avoid separation of the coating (Lustre Paint) from the orthoprosthesis. Sandblast the surface of orthoprosthesis with alumina particles, 50 µm diameter before applying the connector.
- Avoid desorption of the orthoprosthesis. Use a dental adhesive resin cement to lute jacket crowns, veneers, inlays or onlays. Sandblast the inner surface of the jacket crowns, veneers, inlays or onlays using low pressure (1.5 kg/cm<sup>2</sup>).
- The oxygen inhibited layer of the resin surface to be coated should be removed before applying Lustre Paint as it may hinder the proper curing of Lustre Paint.
- For curing Lustre Paint, violet LED (effective wave length: 390-430nm) is required.
- Place the light source as close as possible to the resin surface. Set a 2cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or area in shadow.

## CAUTION

1. In rare cases, this product may cause sensitivity in some people. If such reactions are experienced, discontinue the use of the product and consult a physician.
2. Wear eye protection and plastic or rubber gloves during operation in particular to avoid direct contact with paste and air inhibited resin layers in order to prevent possible sensitivity.
3. LUSTRE PAINT DILUTING LIQUID, SEPARATOR, DIE HARDNER, GC METAL PRIMER II / GC METAL PRIMER Z and GC CERAMIC PRIMER II are highly flammable liquid and vapour. Do not use near naked flame. Keep away from sources of ignition. Keep away from direct sunlight.
4. When adjusting or polishing the veneer surface, use a dust collector and wear a dust mask to avoid inhaling dust.
5. In case of contact with oral tissues or skin, immediately remove with cotton. Flush with water.
6. In case of contact with eyes, immediately wash the eye with copious amounts of water and seek medical attention.
7. This product should not be used for malocclusion, clenching or bruxism.
8. At the time of cementation, sandblast the inner surface of metal free restorations, and then be sure to treat the inner surface with a silane coupling agent (e.g. CERAMIC PRIMER II).
9. Dispose all waste according to local regulation.
10. To prevent cross contamination between patients, do not reuse the dispensing tip.
11. The dispensing tip cannot be sterilized in an autoclave or chemiclave.
12. Wear protective eye glasses during light curing.
13. Lower light intensity may cause insufficient curing.

## STORAGE

Recommended for optimal performance, store in a cool place (4-25°C / 39.2-77.0 °F) away from high temperatures or direct sunlight. Shelf life: 3 years from the date of manufacture.

Some products referenced in the present IFU may be classified as hazardous according to GHS.

Always familiarize yourself with the safety data sheets available at:

<http://www.gceurope.com>

They can also be obtained from your supplier.

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## VERPACKUNGSGRÖSSEN

### 1. Sets

- **LAYER SET (19 Farbtöne)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 Farbtöne)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Gelb, LB-Rot, LB-Grau, LB-Blau, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 Farbtöne)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blau, LP-Grau, LP-Creme, LP-Lavendel, Lustre Paint Verdunnungsflüssigkeit, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 Farbtöne)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Hellrot, GLP-Lila, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 22 Blue
- **ONE BODY SET (11 Farbtöne)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **ACCESSORY KIT (Zubehör Set)**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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- **OPAQUE, 2,0-ml-Spritze (in 13 Farbtönen)**  
O-Base, OA, OB, OC, OD, OA1, OA2, OA3, OA3.5, OB1, OB3, OC3, OD2, Dispensing Tip Needle Type, Light Protective Cover
- **PASTE HB, 3,3-ml-Spritze (in 18 Farbtönen)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-DW, HB-ODA, HB-ODB, HB-ODC, HB-ODD, HB-ODW, HB-EL, HB-ED, HB-CLF, HB-PE
- **PASTE LB, 2,0-ml-Spritze (in 15 Farbtönen)**  
LB-Base D, LB-Base E, LB-Base CLF, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Rot, LB-Gelb, LB-Blau, LB-Grau, LB-Milky, LB-Inlay E, LB-Inlay TD, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **ONE BODY, 2,0-ml-Spritze (in 5 Farbtönen)**  
LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Opaque, 2,0-ml-Spritze (in 2 Farbtönen)**  
GO-1, GO-2, Dispensing Tip Needle Type, Light Protective Cover
- **GUM Shades HB, 3,3-ml-Spritze (in 4 Farbtönen)**  
GHB-1, GHB-2, GHB-3, GHB-CL
- **GUM Shades LB, 2,0-ml-Spritze (in 4 Farbtönen)**  
GLB-1, GLB-2, GLB-3, GLB-CL, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Lustre Paint, 0,8-ml-Spritze (in 2 Farbtönen)**  
GLP-Bright Red, GLP-Violet, Dispensing Tip Needle Type Small, Light Protective Cover
- **LUSTRE PAINT, 0,8-ml-Spritze (in 10 Farbtönen)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-Blau, LP-Creme, LP-Grau, LP-Lavendel, LP-CL, Dispensing Tip Needle Type Small, Light Protective Cover

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## PRESENTACIONES

### 1. Set

- **LAYER SET (19 colores)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 colores)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 colores)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 colores)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (11 colores)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **ACCESSORY KIT**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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- **OPAQUE, 2.0mL jeringa (en 13 colores)**  
O-Base, OA, OB, OC, OD, OA1, OA2, OA3, OA3.5, OB1, OB3, OC3, OD2, Dispensing Tip Needle Type, Light Protective Cover
- **PASTE HB, 3.3mL jeringa (en 18 colores)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-DW, HB-ODA, HB-ODB, HB-ODC, HB-ODD, HB-ODW, HB-EL, HB-ED, HB-CLF, HB-PE
- **PASTE LB, 2.0mL jeringa (en 15 colores)**  
LB-Base D, LB-Base E, LB-Base CLF, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Red, LB-Yellow, LB-Blue, LB-Grey, LB-Milky, LB-Inlay E, LB-Inlay TD, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **ONE BODY, 2.0mL jeringa (en 5 colores)**  
LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Opaque, 2.0mL jeringa (en 2 colores)**  
GO-1, GO-2, Dispensing Tip Needle Type, Light Protective Cover
- **GUM Shades HB, 3.3mL jeringa (en 4 colores)**  
GHB-1, GHB-2, GHB-3, GHB-CL
- **GUM Shades LB, 2.0mL jeringa (en 4 colores)**  
GLB-1, GLB-2, GLB-3, GLB-CL, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Lustre Paint, 0.8mL jeringa (en 2 colores)**  
GLP-Bright Red, GLP-Violet, Dispensing Tip Needle Type Small, Light Protective Cover
- **LUSTRE PAINT, 0.8mL jeringa (en 10 colores)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-Blue, LP-Cream, LP-Grey, LP-Lavender, LP-CL, Dispensing Tip Needle Type Small, Light Protective Cover

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\*Jmf. Konfigurationstabell: Opaker, blandningsrelation

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### 1. Set fp.

- **LAYER SET (19 pastor)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 pastor)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 färger)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 färger)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (11 pastor)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **TILLBEHÖRSKIT**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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## EMBALAGENS

### 1. Conjuntos

- **Conjunto LAYER SET (19 tonalidades)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **Conjunto LAYER PRO SET (19 tonalidades)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **Conjunto PAINT SET (10 tonalidades)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid (liquido de diluicao), Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **Conjunto GUM SHADES SET (11 tonalidades)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid (liquido de diluicao), Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **Conjunto ONE BODY SET (11 tonalidades)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **KIT DE ACESSÓRIOS**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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- **OPAQUE, 2,0mL seringa (em 13 tonalidade)**  
O-Base, OA, OB, OC, OD, OA1, OA2, OA3, OA3.5, OB1, OB3, OC3, OD2, Dispensing Tip Needle Type, Light Protective Cover
- **PASTE HB, 3,3mL seringa (em 18 tonalidade)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-DW, HB-ODA, HB-ODB, HB-ODC, HB-ODD, HB-ODW, HB-EL, HB-ED, HB-CLF, HB-PE
- **PASTE LB, 2,0mL seringa (em 15 tonalidade)**  
LB-Base D, LB-Base E, LB-Base CLF, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Red, LB-Yellow, LB-Blue, LB-Grey, LB-Milky, LB-Inlay E, LB-Inlay TD, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **ONE BODY, 2,0mL seringa (em 5 tonalidade)**  
LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Opaque, 2,0mL seringa (em 2 tonalidade)**  
GO-1, GO-2, Dispensing Tip Needle Type, Light Protective Cover
- **GUM Shades HB, 3.3mL seringa (em 4 tonalidade)**  
GHB-1, GHB-2, GHB-3, GHB-CL
- **GUM Shades LB, 2,0mL seringa (em 4 tonalidade)**  
GLB-1, GLB-2, GLB-3, GLB-CL, Dispensing Tip Plastic Type Wide, Light Protective Cover
- **GUM Shades Lustre Paint, 0,8mL seringa (em 2 tonalidade)**  
GLP-Bright Red, GLP-Violet, Dispensing Tip Needle Type Small, Light Protective Cover
- **LUSTRE PAINT, 0,8mL seringa (em 10 tonalidade)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-Blue, LP-Cream, LP-Grey, LP-Lavender, LP-CL, Dispensing Tip Needle Type Small, Light Protective Cover

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## AMBALARE

### 1. Seturi

- **LAYER SET (19 nuanțe)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 nuanțe)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 nuanțe)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 nuanțe)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (11 nuanțe)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **KIT ACCESORII**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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## PAKIRANJE

### 1. Seti pakiranj

- **LAYER SET (19 barv)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 barv)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad
- **PAINT SET (10 barv)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid (razredčilo), Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 barv)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid (razredčilo), Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (11 barv)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **SET DODATKOV**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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## PAKOVANJE

### 1. Komplet pakovanje

- **LAYER SET (komplet za slojevanje) (19 boja)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (komplet za slojevanje) (19 boja)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange (narandžasta), LB-Yellow (žuta), LB-Red (crvena), LB-Grey (siva), LB-Blue (plava), LB-Milky (mlečno bela), LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (komplet za bojenje) (10 boja)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue (plava), LP-Grey (siva), LP-Cream (krem), LP-Lavender (lavanda), Lustre Paint Diluting Liquid (tečnost za razređivanje), Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (komplet boja za gingivu) (11 boja)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red (svjetlo crvena), GLP-Violet (ljubičasta), GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid (tečnost za razređivanje), Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (komplet jednokomponentnih boja) (11 boja)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **KOMPLET PRIBORA**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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- **LAYER SET (19 відтінків)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **LAYER PRO SET (19 відтінків)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAINT SET (10 відтінків)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid, Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **GUM SHADES SET (11 відтінків)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY SET (11 відтінків)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **ACCESSORY KIT**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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## PAKUOTĖS

### 1. Rinkinys

- **SLUOKSNIIVIMO RINKINYS (19 spalvų)**  
HB-DA1, HB-DA2, HB-DA3, HB-DA3.5, HB-DB1, HB-DB3, HB-DC3, HB-DD2, HB-EL, HB-ED, HB-CLF, HB-PE, HB-ODA, HB-ODB, O-Base, OA, OB, OC, OD, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PAPILDOMO SLUOKSNIIVIMO RINKINYS (19 spalvų)**  
HB-ODC, HB-ODD, HB-ODW, HB-DW, LB-Base E, LB-Base CLF, LB-Base D, LB-Base OD, LB-Base Opal, LB-DW, LB-ODW, LB-Orange, LB-Yellow, LB-Red, LB-Grey, LB-Blue, LB-Milky, LB-Inlay E, LB-Inlay TD, Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **DAŽYMO RINKINYS (10 spalvų)**  
LP-A, LP-B, LP-C, LP-D, LP-CLF, LP-CL, LP-Blue, LP-Grey, LP-Cream, LP-Lavender, Lustre Paint Diluting Liquid (skiediklis), Dispensing Tip Needle Type Small & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B
- **DANTENŲ MODELIAVIMO RINKINYS (11 spalvų)**  
O-Base, GO-1, GO-2, LP-CL, GLP-Bright Red, GLP-Violet, GLB-1, GLB-2, GLB-3, GLB-CL, GHB-2, Lustre Paint Diluting Liquid (skiediklis), Mixotip, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type Small & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Disposable Palette, Mixing Pad No. 14B, Plastic Spatula No. 2 Blue
- **ONE BODY RINKINYS (11 spalvų)**  
O-Base, OA, OB, OC, OD, LB-Base OD, LB-A, LB-B, LB-C, LB-D, LB-W, Dispensing Tip Plastic Type Wide & Light Protective Cover, Dispensing Tip Needle Type & Light Protective Cover, Mixing Pad No. 22, Plastic Spatula No. 2 Blue
- **PRIEDAI**  
AIR BARRIER, SEPARATOR, DIE-HARDNER, CERAMIC PRIMER II, METAL PRIMER Z, Acrylic Primer, DIAPOLISHER PASTE, MODELLING LIQUID, Brush Round No. 1, Brush Flat No. 1, Brush Holder (ivory & white handle), Brush No. 7, Shade Guide Kit, Mixing Pad No. 22

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