# A-Silicone for Laboratory Use

## DE GEBRAUCHSINFORMATION

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EN INSTRUCTIONS FOR USE
-------------------------

### HIGH-PRECISION ADDITION CURED SILICONES (POLYVINYLSILOXANES) FOR THE DENTAL LABORATORY

**Proclinic A-Silicone** – Rigid setting (85 Shore-A) lab putty mixed in a 1:1 ratio characterised by an excellent balance between elasticity and final hardness.

#### Use

- Matrix for constructing provisional prostheses.
- Addition or replacement of elements on full dentures or metal frame works.
- Matrix for use when making a removable gingival mask.
- Partial counter die in flask for partial or full-full dentures.
- Blocking in the repair of dentures.
- Occlusal regulation key for positioning in an articulator
- Duplication of full dentures and gypsum models without using insulators.
- Support for working with light-curing composite resins.
- Control matrix for designing a metal framework.
- Ceramic moulded on articulator technique.

**Note:** The decision to use Proclinic A-Silicone is at the operator's discretion, whether he believes it more appropriate to use an extremely rigid silicone or a rigid silicone with some elasticity.

#### INSTRUCTIONS FOR USE

##### Dosage:

Take equal parts of the two components from the tubs using exclusively the dosing spoons provided (1). **Do not swap over the dosing spoons.**

##### Mixing / use:

Mix the Base and Catalyst vigorously for about 30 sec. until a uniformly coloured putty free of streaks is obtained (2). Then start to model the putty as required. The working time, including the mixing phase, is about 2 min. at room temperature (3). The material sets starting from a minimum of about 4 min. \*\*The prosthesis may then be removed (4).

#### WARNINGS / PRECAUTIONS

The 1:1 proportions must be strictly respected to obtain the correct times and avoid altering the products' characteristics. A higher dose of the catalyst does not reduce the setting time. Close the tubs immediately after use, take care not to swap the lids, or dosing spoons. The setting of addition silicones may be inhibited by latex gloves. We recommend the use of vinyl gloves or thorough hand washing. The surfaces with which the silicone comes into contact must be cleaned, degreased and dried.

#### STORAGE

Shelf life of 24 months if stored at a temperature between 5° and 27°C (41° - 80°F).

TECHNICAL DATA	Proclinic A-Silicone
Mixing time	30"
Working time*	2'
Setting time*	8'
Removal time**	6'
Shore-A hardness	85

\*The times mentioned above are intended from the start of the mixing phase at 23°C (73°F).

\*\*The times indicated may vary according to room temperature and mixing technique.

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**Proclinic S.A.**