Prior to use, carefully read the instructions for use.

OSTRON 100

ACRYLIC TRAY AND BASEPLATE RESIN

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

Construction of individual impression trays and baseplates.

WARNING

Contains methyl methacrylate. In susceptible individuals, sensitisation to the product cannot be excluded. Use of the product should be discontinued if allergic reactions are observed.

DIRECTIONS FOR USE

A. INDIVIDUAL TRAY

1. PREPARATION
 - a) Take an impression for making a study model.
 - b) To create space for the impression material and for covering undercuts, apply a suitable spacer to the model.
2. MIXING
 - a) Standard powder/liquid ratio is 10 g/4.5 g and is obtained by one scoop powder (13 g) and a full measure of liquid (5.8 g). Two scoops of powder are sufficient for an upper arch, and one and a half scoops for lower.
 - b) Mix them in the bowl for 30 seconds.
3. FORMING
 - a) Roll out the mixture uniformly to a thickness of 1.5 to 2 mm.
 - b) Place the flattened mixture over the spacer, form a tray and trim off excess.
 - c) Make a handle with residual mixture before hardening. This can be done by applying liquid to both tray and handle.
 - d) When retention holes are needed, the use of amalgam-carrier before hardening or round bur after hardening is recommended.
4. TRIMMING
 - a) After Ostron becomes hard, remove it from the model and brush the inside clean. If a wax spacer is used, it must be washed out with hot water containing detergent.
 - b) Smooth the edges using a suitable instrument.
5. SETTING
 - a) Apply the impression material adhesive to the inside and edge of the tray, and leave until it dries.
 - b) Load the impression material into the tray and set it into the mouth.

B. BASE PLATE

- a) Soak the model in water or apply a separator (GC Acro Sep) to the mucosal surface of the model.
- b) After mixing, roll out the mixture and form the base plate on the model.
- c) Trim off the excess on the model using a wax-carver before hardening.

STORAGE

Store in a cool place away from direct sunlight. (Shelf life: 3 years from date of manufacture)

PACKAGES

1. Powder : 500 g tin, 1 kg tin, 5 kg tin, 10 kg tin (Europe : 1 kg only).
2. Liquid : 250 g (270 mL) bottle, 2.5 kg (2.7 L) tin (Europe : 250 g only).

SHADES

Transparent blue, Blue, Pink, White, Clear (Europe : transparent blue, clear only).

WARNING

1. The liquid is volatile and highly flammable. Close the container immediately after use and keep away from sources of ignition.
2. Avoid direct contact between skin and uncured resin. It is recommended to wear rubber or plastic gloves.
3. Avoid contact with skin or eyes. In case of contact with eyes, flush immediately with water and seek medical attention.

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