Biosil® f

Instructions for Use	English
Mode d'emploi.....	French
Instrucciones para el uso	Italian
Verarbeitungsanleitung	German
Instrucciones de uso	Spanish
Gebruiksaanwijzing	Dutch
Naudojimo instrukcija.....	Lithuanian
Návod k použití.....	Czech
Návod na použitie.....	Slovakian
Használati utasítás.....	Hungarian
Οδηγίες χρήσης	Greek

Instructions for use

Biosil[®] f

Gebrauchsanweisung
Mode d'emploi
Instrucciones para el uso
Istruzioni per l'uso

 Dentsply
Sirona

Biosil® f

en	Instructions for use	04 – 08
de	Gebrauchsanweisung	9 – 13
fr	Mode d'emploi	14 – 18
it	Istruzioni per l'uso	19 – 23
es	Instrucciones para el uso	24 – 28

1 Indications for use

en

Biosil® f is a CoCr alloy for the fabrication of partial model cast dentures.

For dental use only.

1.1 Technical data

Composition in mass %

Co	Cr	Mo	Si	Mn	C
64.8	28.5	5.3	0.5	0.5	0.4

Biosil® f is free from nickel and beryllium and may contain Fe.

Color	white
Melting range	1320 – 1380 °C
Vickers hardness	400 HV 10
0.2% yield strength	700 MPa
Tensile strength	900 MPa
Modulus of elasticity	220 GPa
Elongation at fracture	5%
Density	8.4 g/cm ³
Preheating temperature	1000 °C
Casting temperature	1500 °C

It is in accordance with the international standard ISO 22674.

2 Contraindications

Do not use in case of hypersensitivity to one or more of the metal constituents of the alloy.

3 Warnings

In the mechanical processing of dental alloys, always work with local aspiration plus face or respiratory protection.

Cobalt - vapour is released during processing.

May cause cancer.

May release metal vapours when melting.

Do not breathe Dusts and vapours.

Ensure suitable suction / aeration at the work place and with operational machinery.

4 Precautions

Avoid occlusal and approximal contact of different alloy types.

5 Adverse reaction

Allergies to metals contained in the alloy as well as electrochemically conditioned dysesthesia are possible.

In individual cases systemic side effects of metals contained in the alloy have been reported.

6 Step-by-Step instructions

6.1 Sprues

2 – 4 sprues – depending on the size of the frame work should be attached to each object. In order to facilitate flow of way and controlled solidification, wax sprue wire with a diameter of 3 – 4 mm should be attached to the thickest points of the casting tapering. The length of the sprues must be determined individually by the selection and the position of the sprue former.

6.2 Investing

Apply a thin coat of Waxit[®] wetting agent to the wax object or spray with Waxit[®] and then blow it dry carefully (do not use compressed air). Invest the pattern with Optivest[®] or Biosint[®] investment. Observe the working instructions of the investment that is used. A fine investment can be used. However, it is not required if a fine-grain investment compound like e.g. Optivest[®] is used.

6.3 Preheating

Preheat slowly with the suction system activated. Adhere to the heating rates and holding times given by the manufacturer of the respective investment. After the final temperature of 1000 °C has been reached, preheat the casting ring for 60 min. In the case of particularly filigree constructions the preheating temperature must be raised to 1050 °C. The holding time must be extended correspondingly, if a larger number of casting rings is used.

6.4 Casting

Due to their shape the ingots can be placed easily into the ceramic crucibles. Use ceramic crucibles only.

a) HF-casting unit (e.g. Degutron®): Do not use graphite or carbon crucibles. After a time of further heating of 5 – 7 s, depending on the quantity of alloy, casting is carried out before the oxide skin tears up.

b) Flame-melting (e.g. Motorcast® compact): melt with the reducing flame of a propane-oxygen mixture. After exceeding the liquidus temperature and a period of continued heating – 5 – 15 s depending on the alloy quantity – begin casting.

c) Arc (e.g. Degumat®): melt with the arc on step 3 – 4. Begin the casting process after a period of further heating – 2 – 7 s depending on the alloy quantity. Since melting in the Degumat® unit is carried out under protective gas (argon), no oxide film is formed on the melt.

6.5 Devesting

After casting let the casting ring cool down at room temperature for at least 20 min. Quenching in water is not recommended because of the risk of deformation of the casting. After devesting blast with blasting agent. In order not to destroy the passivation layer on the surface, CoCr alloys **should not be pickled**.

6.6 Finishing

All finishing instruments that are suitable for cast dentures can be used. Use covering varnish to protect filigree elements, e.g. clasps, against abrasion in case of electro-polishing.

6.7 Soldering

Use flux Oxynon[®] for soldering of Biosil[®] f to precious metal alloys. The solder that is selected depends on the precious metal alloy to be soldered. For soldering Biosil[®] f we recommend the Degudent[®]-Lot (solder) N1W.

6.8 Reusability

The fabrication of high-quality dentures requires the use of extremely pure materials. Therefore, to avoid risks only original Biosil[®] f ingots should be used.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

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