# FREEPRINT®

## crown

3D PRINTING MATERIAL

MED RESIN

385 nm

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Composite for dental 3D printing

## **INDICATION**

Permanent single crowns (anterior and posterior areas), permanent inlays and onlays, veneers, long-term provisional crowns, bridges and permanent denture teeth.

## **PATIENT TARGET GROUP**

Persons being treated in the context of a dental procedure.

## **INTENDED USERS**

Dentist, dental technician

## **SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE**

see "Annex 1" (enclosed separately)

## **PROCESSING**

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Do not use thermal methods for disinfection or sterilization. This might deform the workpiece.
- ▶ **FREEPRINT® crown** is used for the fabrication of permanent single crowns in the anterior and posterior region, denture teeth, inlays and onlays. It is also used for long-term provisional crowns and bridges having up to 3 pontics with an intermediate pontic and veneers. The minimum occlusal thickness is 1.5 mm, circularly 1.5 mm and cervically 1.0 mm. The connection cross-section for intermediate pontics of a bridge must be at least 16 mm<sup>2</sup>. Failure to observe this information may lead to an undesirable result.
- ▶ A layer thickness of 50 µm is suitable for the printing process.
- ▶ Mechanically polish the surface. Pre-polishing is done with rotating brushes and pre-polishing paste, a high-gloss polish with buffing and high-gloss polishing agents for resin composite.
- ▶ The restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.
- ▶ **Attachment of permanent denture teeth**  
Once the denture base (see additional IFU **FREEPRINT® denture**) and the denture teeth have been removed from the respective build platforms and the corresponding cleaning processes completed, leave the printed components in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again. Apply liquid **FREEPRINT® denture** to the receptacles with a small brush or a syringe, and connect the repositioned denture teeth by exposure in a xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2 × 1800 flashes under inert gas conditions (nitrogen).
- ▶ Processing temperature 23 °C ± 2 °C.

## **SAFETY INFORMATION**

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (e.g. protective gloves, goggles) when handling the uncured material.

- ▶ For finishing the cured material, suitable personal protective equipment must be used (protective gloves, goggles, face mask).
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is ensured only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

### NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet!

### For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent public authority of the Member State in which the user and/or patient resides.

### Brief Report on Safety and Clinical Performance (SSCP):

The Summary Safety and Clinical Performance Report (SSCP) for BASIS UDI-DI ++EDET0996K (freeprint® crown)

is stored in the European Database for Medical Devices (EUDAMED) under the internet link

<https://ec.europa.eu/tools/eudamed>.

### STORAGE

- ▶ Store **FREEPRINT® crown** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

### CONTRAINDICATION

Contains (meth)acrylates and phosphine oxides.

Ingredients of **FREEPRINT® crown** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **FREEPRINT® crown** intraorally only when fully polymerized.

### SIDE EFFECTS

Product may cause allergic reactions.

### DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

#### Storage:

#### Processing:

at 23 °C ± 2 °C

# Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

## Building process

Creation of a print job in compliance with the machine and material parameters, see table

## Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

## Pre-cleaning

Remove components from the platform and remove excess resin with a gentle stream of (pressurized) air.

## Cleaning

see "Annex 1, CLEANING EQUIPMENT"

Then clean cavities and surfaces additionally with compressed air and turn the component over. Remove support structures from the component.

## Main cleaning

see "Annex 1, CLEANING EQUIPMENT"

Dry cavities and surfaces with compressed air. The surface must no longer be shiny, otherwise resin residues will be present.

## Drying

Store the component protected from light at room temperature for 30 minutes before post-exposure.

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DETAX GmbH