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Light-curing resin based on (meth)acrylate, for cast parts in dental precision casting, burn-out for DLP printers with UV-LED 385 nm

Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm
MAX UV Pro2 UV Pico2 UV	MiiCraft y-Serie	DII-Serie

Printers may only be operated using material parameters authorised by DETAX.

### Processing

- ▶ Phosphate-bonded investment material can be used for fabrication of the casting mould. A preheating temperature of approx. 800 °C, is to be maintained for approx. 30-45 minutes, depending on the size of the casting mould, regardless of the final casting temperature.
- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Processing temperature 23 °C ± 2 °C.
- ▶ A correction or repair can be done with easyform gel LC after printing.

### Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of FREEPRINT® cast 2.0 may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

### Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

### Storage

FREEPRINT® cast 2.0 is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

### Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations..

### Area of application:

Fabrication of residue-free combustible cast objects for precision casting technology employing the DLP process.

### Processing:

at 23 °C ± 2 °C

### Storage

### Ordering information:

**FREEPRINT® cast 2.0 385**  
 500 g, red **02548**  
 1.000 g, red **02632**

**FREEPRINT® model 2.0 385**  
 1.000 g  
 sand **02128**  
 light-grey **02099**  
 grey **02177**  
 caramel **02850**  
 white **02148**

500 g  
 caramel **02844**

**FREEPRINT® gingiva 385**  
 500 g **02820**  
 1.000 g **02843**

## Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

### Construction process

Generation of a Print Job complying with machine and material parameters

### Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

### Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

### Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

### Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to post-curing, check the openings and additional bore holes for residues.

### Post-curing

Post-curing is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.