# dima Print Stone

## dima Print Stone

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### GB Instructions for use

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### SE Bruksanvisning

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### NO Bruksanvisning

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**KULZER**  
MITSUBI CHEMICALS GROUP

**Manufacturer:**  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)  
+49 800 4372522

**Made in Germany**

Distributed in USA/Canada exclusively by:  
Kulzer, LLC  
4315 South Lafayette Blvd.  
South Bend, IN 46614-2517  
1-800-431-1785

## GB Instructions for use dima Print Stone

**1. Intended purpose**  
Photopolymer for the fabrication of dental models. For use by dental health care professionals only.

**1.1** Indications  
For 3D fabrication of dental models

**1.2** Contraindications  
Not for use in the mouth. In case of known or suspected allergy to components of the product, its use is contraindicated.

**1.3** Potential side effects  
Hypersensitivities to the product or its components cannot be excluded in individual cases.

**2. Product features**

**2.1** Package sizes/Supply Forms  
1000g bottle  
Available shades: beige / gray

**2.2** Accessories  
The material can be printed using printers with a wavelength of 405 or 385 nm (currently released devices – see table below).  
- Stone beige: 405 and 385nm wavelength.  
- Stone gray: 385nm wavelength only.

The cleaning takes place in the cleaning device, e.g. cara Print Clean *pro* or ultrasonic bath. Post-curing of the material must be carried out with a suitable exposure device.

dima Print Stone beige			
3D printer / Manufacturer		Post-curing unit / Manufacturer	
• cara Print 4.0 • cara Print 4.0 <i>pro</i>	Kulzer	• HiLite power 3D • cara Print LEDCure	Kulzer
• MAX UV • PRO 2 UV • PRO 4K	Asiga	• HiLite power 3D • cara Print LEDCure	Kulzer

dima Print Stone gray			
3D printer / Manufacturer		Post-curing unit / Manufacturer	
cara Print 4.0 <i>pro</i>	Kulzer	• HiLite power 3D • cara Print LEDCure	Kulzer
MAX UV	Asiga	• HiLite power 3D • cara Print LEDCure	Kulzer

Always use the latest firmware version of your 3D printer.  
Layer thickness: 50 µm / 70 µm / 100 µm

**2.3** Composition  
Contents: (Meth-)acrylate containing monomers, photoinitiators, stabilizers, (stone gray only: polymer bead), pigments

### 3. Product use

Materials		dima Print Stone beige / gray	
Homogenization		Shake for five minutes	
Cleaning solution		Isopropanol (pure)	
Pre-cleaning		3 min.	
Post-cleaning (fresh rinse)		2 min.	
Post-curing		• HiLite power 3D 2 x 5 min (1x top, 1x bottom) • cara Print LEDCure in accordance with the exposure program	

Processing temperatur: 18–30°C (64–86°F) at 30-80% humidity.

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## Explanation of symbols on labelling

**REF** Reorder no.

**LOT** Batch code

**Manufacturer**

**Date of manufacture**

**Keep dry**

**Keep away from sunlight**

**Consult instructions for use**

**Use-by date**

**Storage temperature**

Please, carefully observe the following safety instructions and all processing notes in other sections.

This product is intended to be used only in accordance with this Instructions for Use. Any use of this product inconsistent with these instructions is at the discretion and sole responsibility of the dental health care professional.

#### 4.1 Warnings and safety precautions

**Causes skin irritation.** Causes serious eye damage.  
**May cause an allergic skin reaction.** Harmful to aquatic life with long lasting effects. Wash thoroughly after handling. Avoid release to the environment. Wear protective gloves / eye protection / face protection. If on skin: Wash with plenty of water. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention.

#### 4.2 Storage conditions

Do not use the material after the expiration date. Do not store above room temperature of 25°C (77°F). Keep away from sunlight. Store in a dry place.

**4.3 Disposal Information**  
Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

**5. Compliance notice**  
Please quote batch number and the article number in all correspondence about the product. Article number and batch number: See note on packaging.

Additional information can be found on the download page at: **www.kulzer.com/cara-print**

<p><b>4.2</b> Application subtheme</p> <p>When filling the vat or when refilling the liquid into the bottle, it must not be mixed with other photopolymers, as otherwise negative changes in the product properties (colour, mechanical and chemical properties) may occur.</p> <p>a. Shake vigorously in the original container for five minutes before use.</p> <p>Inadequate shaking may result in color variations and printing errors.</p> <p>b. Print the objects according to the instruction manual supplied with the 3D printer. Select the print program dima Print Stone beige or dima Print Stone gray.</p> <p>c. Thoroughly remove remaining liquids using compressed air.</p> <p>d. Pre-clean the printed objects for three minutes in isopropanol using the appropriate supplies (see Item 2.2).</p> <p>Note: Thoroughly rinse the inside of the structures with the spray bottle (isopropanol). After cleaning, make sure the printed objects are dry and free from solvent residues.</p> <p>e. Repeat step d using pure isopropanol for two minutes.</p> <p>Note: The objects should not remain in the isopropanol for more than five minutes, otherwise they will be damaged.</p> <p>f. Thoroughly rinse the inside of the structures again with the spray bottle (isopropanol). After cleaning, make sure the printed objects are dry and free from solvent residues.</p> <p>g. Place the printed object for final polymerization in the photo polymerization unit according to the table above. Note: The final properties depend on the adherence to the recommended post-curing process.</p> <p><b>4.3</b> General information</p> <p>Please, carefully observe the following safety instructions and all processing notes in all other sections. This product is intended to be used only in accordance with this Instructions for Use. Any use of this product inconsistent with these instructions is at the discretion and sole responsibility of the dental health care professional.</p> <p><b>4.1</b> Warnings and safety precautions</p> <p>Causes skin irritation. Causes serious eye damage.  May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Wash thoroughly after handling. Avoid release to the environment. Wear protective gloves / eye protection / face protection. If on skin: Wash with plenty of water. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention.</p> <p><b>4.2</b> Storage conditions</p> <p>Do not use the material after the expiration date. Do not store above room temperature of 25°C (77°F). Keep away from sunlight. Store in a dry place.</p> <p><b>4.3</b> Disposal Information</p> Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances. <p><b>5. Compliance notice</b></p> Please quote batch number and the article number in all correspondence about the product. Article number and batch number: See note on packaging. <p>Additional information can be found on the download page at: <b>www.kulzer.com/cara-print</b></p> <p>Date: 2021-07</p>	
--	--

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