# FREEPRINT®

## IBT


3D PRINTING MATERIAL

MED RESIN

385 nm

**indirect bonding tray**

elastic, tear resistant  
biocompatible, Class I

# DETAX

<b>DE</b>	Gebrauchsanweisung.....	3	<b>NO</b>	Bruksanvisning.....	36
<b>EN</b>	Instructions for use.....	6	<b>PL</b>	Instrukcja stosowania.....	39
<b>FR</b>	Mode d'emploi.....	9	<b>PT</b>	Manual de instruções.....	42
<b>ES</b>	Modo de empleo.....	12	<b>RO</b>	Instrucțiuni de utilizare.....	45
<b>IT</b>	Istruzioni per l'uso.....	15	<b>SV</b>	Bruksanvisning.....	48
<b>DA</b>	Brugsanvisning.....	18	<b>SK</b>	Návod na použitie.....	51
<b>EL</b>	Οδηγίες χρήσης.....	21	<b>SL</b>	Navodila za uporabo.....	54
<b>HR</b>	Upute za uporabu.....	24	<b>CS</b>	Návod k použití.....	57
<b>LV</b>	Lietošanas instrukcija.....	27	<b>HU</b>	Használati útmutató.....	60
<b>LT</b>	Naudojimo instrukcija.....	30	<b>RU</b>	Инструкция по применению.....	63
<b>NL</b>	Gebruiksaanwijzing.....	33	<b>TR</b>	Kullanım kılavuzu.....	66

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Resin for dental 3D printing

## INDICATION

Othodontic transfer trays for positioning brackets

## PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

## INTENDED USERS

Orthodontist, dentist, dental technician

## SUITABLE FOR THE FOLLOWING DLP-PRINTER

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV	MiiCraft y-Serie	DII-Serie	Solflex Series

Printers may only be operated using material parameters authorised by DETAX!

## PROCESSING

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 8).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ Maximum curing depth\* at direct post-exposure:  
clear: 6 mm  
\*In case of large objects and exposure on both sides, the material thickness can be up to 12 mm.
- ▶ Processing temperature 23 °C ± 2 °C.

## SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.

- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.
- ▶ Read and understand the safety data sheet!

## For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.de](mailto:incident@detax.de) and to the competent authority of the Member State in which the user and/or patient is established.

## STORAGE

**FREEPRINT® IBT** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

## CONTRAINDICATION

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® IBT** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® IBT** only insert intraorally in completely polymerised state.

## ADVERSE EFFECTS

Product may cause allergic reactions.

## DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

## Symbols:

 Medical Device

## Storage:

## Application:

At 23 °C ± 2 °C

# Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to post-exposure, check the openings and additional bore holes for residues. Then blow off with compressed air.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 3 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between. Then clean the platforms fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. and blow off with compressed air.

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GmbH & Co. KG