# dimaprint Mouth Guard

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Distributed in USA/  
Canada exclusively by:  
Kulzer, LLC  
4315 South Lafayette Blvd.  
South Bend, IN 46614-2517  
1-800-431-1785

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**3.2** Application substeps  
When filling the vat or when refilling the liquid into the bottle, it must not be mixed with other photopolymers, as otherwise negative changes in the product properties (colour, mechanical and chemical properties) may occur.

**a.** Shake vigorously in the original container for five minutes before use. Inadequate shaking may result in color variations and printing errors.

**b.** Print the objects according to the instruction manual supplied with the 3D printer. Select printing parameters for dimaprint Mouth Guard (white or blue).

**c.** Thoroughly remove remaining liquids using compressed air. **d.** Pre-clean the printed objects for three minutes in Isopropanol using the appropriate squeegee (see Item 2.2).

**e.** Repeat step **d** using pure Isopropanol for two minutes. **Note:** The objects should not remain in the Isopropanol for more than five minutes, otherwise they will be damaged. **Note:** Thoroughly rinse the inside of the structures again with the spray bottle (Isopropanol). After cleaning, make sure the printed objects are dry and free from solvent residues.

**f.** Place the printed object for final polymerization in the photo polymerization unit according to the table above using the spray bottle (Isopropanol). After cleaning, make sure the printed objects are dry and free from solvent residues.

**g.** Additional cleaning after postcuring: leave printed objects in fresh Isopropanol for 3min, after that in water for 10 min.

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Materials	dimaprint Mouth Guard				
Homogenization	Shake for five minutes				
Cleaning solution	Isopropanol (pure)				
Pre-cleaning	3 min.				
Post-cleaning (fresh rinse)	2 min.				
Post-curing	<div><ul style="list-style-type: none"><li>• HiLite power 3D 2 x 5 min (1x top, 1x bottom)</li><li>• cara Print LEDcure in accordance with the exposure program</li></ul></div>				
Additional cleaning after postcuring	3 min. in fresh Isopropanol and 10 min in water				

**3.3** Application notes  
Store the Mouth Guard in an adequate ventilated boxes to avoid dima.

**Note:** The cold dima Print Mouth Guard can be too stiff, this will disappear after a few minutes of wearing. The cumulative wearing time should not exceed 30 days (or 30 x 24 hours).

**4. General information**  
Please, carefully observe the following safety instructions and all processing notes in other sections. This product is intended to be used only in accordance with this Instructions for Use. Any use of this product inconsistent with these instructions is at the discretion and sole responsibility of the dental health care professional.

**4.1** Warnings and safety precautions  
Avoid contact with skin and eyes. May cause skin and eye irritation. Skin contact may cause sensitization. Wear protective gloves/protective clothing/eye protection/face protection. In case of contact of the product with skin, immediately wash with plenty of water and soap. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing if symptoms persist, consult an eye specialist. Toxic to aquatic life with long lasting effects.

**4.2** Storage conditions  
Do not use the material after the expiration date. Store at room temperature, not above 25°C (77°F). Avoid direct sunlight. Store in a dry place.

**4.3** Disposal information  
Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially empty packaging in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

**5. Complaint notice**  
Please quote batch number and article number in all correspondence about the product.

Additional information can be found on the download page at: [www.kulzer.com/cara-print](http://www.kulzer.com/cara-print)

## Instructions for use dimaprint Mouth Guard

**1. Intended purpose**  
Photopolymer for the fabrication of sport Mouth Guards via additive manufacturing.

**1.1** Limitations for use  
Not for very high impact sports where extreme hard or small area hits are expected. In case of known or suspected allergy to components of the product, its use is contraindicated.

**1.2** Caution  
Hypersensitivites to the product or its components cannot be excluded in individual cases. In case of suspicion please contact the manufacturer for ingredients information.

### 2. Product features

**2.1** Package sizes/Delivery forms  
1000g bottle

### 2.2

The material can be printed using printers with a wavelength of 405 or 385 nm (currently released devices – see table below). The cleaning takes place in the cleaning device, e.g. cara Print Clean pro or ultrasonic bath. Postcuring of the material must be carried out with a suitable exposure device.

3D printer / Manufacturer	Post-curing unit / Manufacturer				
cara Print 4.0	Kulzer	<div><ul style="list-style-type: none"><li>• HiLite power 3D</li><li>• cara Print LEDcure</li></ul></div>	Kulzer		
cara Print 4.0 pro	Kulzer	<div><ul style="list-style-type: none"><li>• HiLite power 3D</li><li>• cara Print LEDcure</li></ul></div>	Kulzer		
Layer thickness: 100 µm					

Always use the latest firmware version of your 3D printer.

### 2.3

**Contents:** (Meth)acrylate-containing monomers, photo initiators, pigments

### 3. Product use

#### 3.1

Lead the bottle at around 25°C (77°F) for 2 hours prior to 3D printing, and thoroughly shake before use. In case that the temperature of ink at printing is not normal room temperature, accuracy of printing may get influenced. Recommended minimum wall thickness 2,5 mm

Materials	dimaprint Mouth Guard				
Homogenization	Shake for five minutes				
Cleaning solution	Isopropanol (pure)				
Pre-cleaning	3 min.				
Post-cleaning (fresh rinse)	2 min.				
Post-curing	<div><ul style="list-style-type: none"><li>• HiLite power 3D 2 x 5 min (1 superior, 1 inferior)</li><li>• cara Print LEDcure según el programa de exposición</li></ul></div>				
Limpeza adicional tras el post-curado	3 min en isopropanol fresco y 10 min en agua				

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