dima[®] Print C&B temp

Manufacturer:
DENTCA, INC.
357 VAN NESS WAY, STE 250
TORRANCE, CA 90501
USA

CE 0044

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GB Instructions for use dima[®] Print C&B temp

Indications for Use

Intended Use
dima Print C&B temp is a light-curable polymerizable resin to fabricate, by additive manufacturing, temporary crowns or bridges. The fabricated temporary crowns or bridges are an alternative to preformed temporary crowns or bridges.
The fabrication of temporary crowns or bridges with dima Print C&B temp requires digital models of crowns or bridges, a stereolithographic additive printer, and curing light equipment.

Intended Use for Canadian customers
dima Print C&B temp is a light-curable polymerizable resin to fabricate, by additive manufacturing, temporary crowns or bridges to be used for less than 30 days.
* The fabricated temporary crowns or bridges are an alternative to preformed temporary crowns or bridges.
The fabrication of temporary crowns or bridges with dima Print C&B temp requires digital models of crowns or bridges, a stereolithographic additive printer, and curing light equipment.
* The device is registered in Canada as class II which is used for less than 30 days.

Requirements
1. Digital crown or bridge model file with minimum thickness 2 mm around the occlusal and shoulder area; STL format file
2. Additive printer and its operation software:

Additive Printer	Operation Software	Provider
cara Print 4.0 or cara Print 4.0 pro	cara Print CAM or cara Print CAM 2.0	Kulzer
Asiga Max, Pro2 or Pro 4K	Asiga Composer	Asiga

3. Curing light equipment

Curing Equipment	Provider
HiLite power 3D / LEDcure	Kulzer

Specific Manufacturing Considerations

1. Digital crown or bridge model file
 - 1.1 File format: STL file
 - 1.2 File size less than 100 MB of each file
 - 1.3 File size should be uploadable in the 3D printer operation software.
2. Additive printer

- 2.1 Hardware
 - a. Wavelength: 385 nm or 405 nm
 - b. Light source
 - Stereolithographic (SLA) method; laser with 25 mW < X < 250 mW
 - Digital Light Processing (DLP) method; high power LED or lasers
 - c. Build Volume: > 103 x 58 x 130 mm (Least fit one arch)
 - d. Build Path: line drawing path or surface layer drawing path
- 2.2 Features of Operation Software
 - a. STL file import
 - b. Automatic rotation and placement
 - c. Layer slicer for path inspection
 - d. Auto and manual generation of supports
- 2.3 Printing Parameter

Printer Model	Layer Thickness (micron)	Recommended orientation angle (degree)	Support point size (mm)	Support density (mm)
cara Print 4.0 or cara Print 4.0 pro	50, 100	Occlusal to build platform	1–2	min. 1
Asiga Max, Pro2 or Pro 4K	50–100	20–90	0.7–1.5	3.0–5.0

2.4 Environmental Conditions

- a. Temperature: 18–30°C (64–86°F)
- b. Relative Humidity: 30–90%
- 2.5 Cleaning Kit
 - Rinse bath and tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, Scraper

3. Recommended Curing light Equipment (Post curing units)

Manufacturer / Model	Curing Chamber	Supply voltage (V / Hz)	Lamp power (W)	Light intensity (mW/cm ²)	Lamp wavelength (nm)	Curing time (min)
Kulzer / HiLite power 3D	Required	100,115,230 / 50–60	200	N/A	390–540	20
Kulzer / LEDcure	Required	100–240 / 50, 60	150	N/A	370–470	Select the program dima C&B temp

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Stand: 2023-06

3.2 Accessories

- a. USP Grade glycerin
- b. Transparent glass container and 2 glass plates
- c. Heat-protective gloves and silicone coated stainless steel tong
- d. Thermocouple

4. Notification

- 4.1 The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users should follow this document to use the device.
- 4.2 If there is a serious incident (death or permanent damage to a patient) that has occurred in relation with this device please report to DENTCA (info@dentca.com) or your local authority of medical device.
- 4.3 When you received damaged or unintentionally open bottles before use, or if the packaging is exposed to environmental conditions outside of the specified in the label, please inform to DENTCA (info@dentca.com).

Warnings:

1. The dima Print C&B temp resin contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If contact with skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
2. Avoid inhalation or ingestion. High vapor concentration can cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the material may cause more serious health effects. Monitor air quality per OSHA standards. Eye Contact: Immediately flush eyes with plenty of clean water for at least 20 minutes, and consult a physician. Wash the contacted area thoroughly with soap and water. Inhalation: In case of exposure to a high concentration of vapor or mist, remove person to fresh air. Give oxygen or artificial respiration as required. Ingestion: Contact your regional poison control center immediately. BURN HAZARD: GLYCERIN BATH CAN REACH TEMPERATURES OF 90°C (-200°F) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.

1. When washing the printed crown and bridge with isopropanol or grinding the crown and bridge, it should be in a properly ventilated environment with proper protective masks and gloves.
2. Store dima Print C&B temp resin at 15–25°C (60–77°F) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date.
3. Expired or unused dima Print C&B temp should be completely cured or polymerized prior disposal.

Adverse Reactions:

1. Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
2. Proper ventilation and personal protective equipment should be used when grinding printed crown and bridge as the particles generated during grinding may cause respiratory, skin and eye irritation.

Procedure to Fabricate the Temporary Crowns and Bridges

1. Printing Preparation
 - a. Select the shade based on prescription. (Recommend to use a different resin tank or tray for the different shade.)
 - b. Open the 3D printer cover and fill the resin tank or tray of the printer with dima Print C&B temp resin up to the required filling line by manufacturer. (When filling the resin into the resin tank or tray, gloves and mask should be used.)
 - c. Close the printer cover.
2. Printing
 - a. Load the crown or bridge model file to be printed in printer operation software which printer manufacturer recommended.
 - b. Use the printer operation software to load the model in order to locate the model in proper position to the build platform.
 - c. For the crown or bridge model, rotate the tip of teeth to face the build platform.
 - d. Generate support sticks on the model using the recommended setting by printer provider. If the support is not enough, add supports on the model. (Avoid the support structures on the connection area between teeth.)
 - e. Start printing.
3. Cleaning
 - a. Detach the printed model from the build platform.
 - b. Use a small flush cutter to remove the support sticks from the printed model.
 - c. Wash the printed object with isopropyl alcohol.
 - d. Use air blowing to dry the printed object or dry it at room temperature under ventilation system or open area.
 - e. Post curing of printed object and finishing
 - a. Smooth the support marks using a bur after washing the printed object with water and drying.
 - b. The printed one should be cured by soaking into the glycerin container for the required curing time under recommended post curing unit. (Glycerin temperature should be greater than 60°C/140°F and it is recommended the glycerin to replace every 80 hrs running or every three months whichever first.)
 - c. Take out the printed one from curing oven using coated tong. (Be careful hot glycerin!).
 - d.

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