Gebrauchsanweisung | Instructions for use |  
Mode d'emploi | Modo de empleo | Modalità d'uso

3D Isoliermittel  
3D separating medium  
Isolant 3D  
Barniz aislante 3D  
Isolante 3D

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## Dear customer

Thank you for choosing a quality product from Dentaaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application.

In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaaurum.com](http://www.dentaaurum.com), even though you may frequently use the same product.

## 1. Manufacturer

Dentaaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

## 2. Product description

The 3D separating medium forms a separating layer between the 3D-printed plastic model and Orthocryl®/Orthocryl® LC.

To fabricate dental models, we recommend using high quality plastic models that were made using DLP (Digital Light Processing). e. g. prime4me® base model (REF 430-110-10). Besides being highly accurate and therefore ensuring very precise fitting, these models allow for the usual process in the laboratory.


Due to the large variety of materials for model printing, it is essential to test the suitability of the 3D separating medium before its use.

### 3. Processing instructions


**3.1** If necessary, block out undercuts in the area of the teeth/jaw during fabrication using the 3D printing process.

**3.2** Scrape the plastic model and, if required, create holes/slits for the expansion screw.


**3.3** Place the expansion screw in position and bend the wire elements.

**3.4** Apply the 3D separating medium evenly with a brush over the entire model before the expansion screw (or screws) and wire elements are finally fixed in place and allow the separating medium to dry thoroughly. Repeat the process in order to achieve an even separating layer everywhere.


**3.5** If required, block out undercuts in the tooth/jaw area and fix the expansion screw(s) and wire elements as usual on the dried plastic model. Use Thermo wax (REF 120-170-00) to avoid a milky effect in the plastic.

**3.6** If necessary, use a wax strip as a boundary for the plastic.

**3.7** Process the plastic as usual. Please follow the corresponding Instruction for use.

**3.8** Residues from the separating medium can easily be removed completely with a steam cleaner or with water and a brush.

## 4. Storage and shelf life

Protect the 3D separating medium from contamination and store in a closed container. Do not use after the expiry date.

## 5. Availability

3D separating medium, 100 ml, REF 162-801-00

## 6. Quality information

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum Implants' own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

## 7. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at [www.dentaurum.com](http://www.dentaurum.com) (Explanation of symbols on labels REF 989-313-00).

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