IPS Ivocolor®

[en] Instructions for Use

Stains and glazes

Rx ONLY

CE 0123

Date information prepared:
2023-04-03 / WW / Rev.0
web

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Characterization of ceramic restorations from Ivoclar Vivadent AG in the anterior and posterior region

Patient target group

Patients with permanent teeth

Intended users / Special training

- Dentists (clinical procedure)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

No special training required.

Use

For dental use only.

Description

IPS Ivocolor® is a universal range of stains and glazes for the individual staining and glazing of ceramic materials. The products are coordinated with the layering ceramics, press ceramics and CAD ceramics of Ivoclar Vivadent AG.

They are compatible with ceramics in the CTE range of $9.0\text{--}18.0 \times 10^{-6}/\text{K}$ (according to ISO 6872:2015).

Product name	Product description
IPS Ivocolor Essence IPS Ivocolor Shades	Characterization of ceramic restorations from Ivoclar Vivadent AG in the anterior and posterior region
IPS Ivocolor Essence Fluid	Mixing of IPS Ivocolor Essence Powder for ceramic restorations in the anterior and posterior region
IPS Ivocolor Glaze Paste IPS Ivocolor Glaze Powder	Glazing of ceramic restorations from Ivoclar Vivadent AG in the anterior and posterior region
IPS Ivocolor Glaze Paste FLUO IPS Ivocolor Glaze Powder FLUO	Glazing of full-contour ceramic restorations from Ivoclar Vivadent AG in the anterior and posterior region
IPS Ivocolor Mixing Liquid Allround IPS Ivocolor Mixing Liquid LongLife	Mixing of IPS Ivocolor Powders and Pastes for ceramic restorations in the anterior and posterior region

Technical data

Characteristics	Unit	Specifications	
Flexural strength	MPa	$\geq 50^*$	
Chemical solubility	$\mu\text{g}/\text{cm}^2$	$\leq 100^*$	
Coefficient of thermal expansion ($25\text{--}T_G$ °C)	$10^{-4}/\text{K}^1$	Glaze	$9.3 \pm 0.5^*$
		Glaze FLUO / Essence / Shade Pigment content $\leq 4\%$	$9.2 \pm 0.5^*$
		Essence / Shade Pigment content > 4 to $\leq 12\%$	$9.1 \pm 0.5^*$
		Essence / Shade Pigment content $> 12\%$	$8.8 \pm 0.5^*$
Glass transformation temperature T_G	°C	$460 \pm 20^*$	
Radioactivity U^{238}	Bq/g	$\leq 1^*$	

* The performance criteria defined in EN ISO 6872:2015 Dentistry – Ceramic materials (ISO 6872:2015) are met.
Classification: Dental ceramics Type I / Class 1

Firing table

Stain and glaze firing when using the staining technique (monolithic)

B [°C]	S* [min]	t \nearrow [°C/min]	T1 [°C]	H [min]	V1 [°C]	V2 [°C]	L** [°C]	tl [°C/min]
IPS e.max® Press								
403	IRT/6	60	710	1	450	709	0	0
IPS e.max® CAD								
403	IRT/6	60	710	1	450	709	0	0
IPS Empress® CAD								
403	IRT/6	60	710	1	450	709	0	0
IPS e.max® ZirCAD								
403	IRT/6	45	710	1	450	709	0	0

B = Stand-by temperature, S = Closing time, t \nearrow = Temperature increase rate, T1 = Firing temperature, H = Holding time, V1 = Vacuum 1, V2 = Vacuum 2, L = Long-term cooling, tl = Cooling rate

* IRT standard mode ** Long-term cooling for the last firing

Note: If several IPS e.max ZirCAD restorations (multi-unit bridges with large pontics or several full-contour restorations) are fired in one firing cycle, heating of the objects to be fired may be delayed. This can be compensated by reducing the heating rate or extending the holding time once the final temperature is achieved.

Stain and glaze firing when using the layering technique (full-coverage veneers, partial veneers)

B [°C]	S* [min]	t \nearrow [°C/min]	T1 [°C]	H [min]	V1 [°C]	V2 [°C]	L** [°C]	tl [°C/f/min]
IPS e.max® Ceram/IPS e.max® ZirCAD								
403	IRT/6	60	710	1	450	709	450	0
IPS e.max® Ceram/IPS e.max® Press								
403	IRT/6	60	710	1	450	709	0	0
IPS e.max® Ceram/IPS e.max® CAD								
403	IRT/6	60	710	1	450	709	0	0
IPS Style® Ceram								
403	IRT/6	60	750	1	450	749	0	0
IPS InLine®								
403	IRT/6	60	830	1	450	829	0	0

B = Stand-by temperature, S = Closing time, t \nearrow = Temperature increase rate, T1 = Firing temperature, H = Holding time, V1 = Vacuum 1, V2 = Vacuum 2, L = Long-term cooling, tl = Cooling rate

* IRT standard mode ** Long-term cooling for the last firing

- The firing parameters listed represent standard values and apply to the Ivoclar Programat® furnaces P310, P510, P710, P300, P500, P700, EP 3010, EP 5010, EP 3000, EP 5000. The temperatures indicated also apply to furnaces of older generations (e.g. P80, P90, P95, P100, X1). If one of these older furnaces is used, however, the temperatures may deviate by ± 10 °C/18 °F depending on the age of the heating muffle.

- If a non-Ivoclar furnace is used, temperature adjustments may be necessary.
- Regional differences in the power supply or the operation of several electronic devices by means of the same circuit may render adjustments of the temperatures necessary.

Indications

None

Areas of application

Staining and/or glazing of ceramic materials from Ivoclar Vivadent AG:

- Veneering ceramics (all-ceramics & metal ceramics)
- Press ceramics
- Milling ceramics
- Zirconium oxide

Contraindication

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

IPS Ivocolor is not suitable for the staining and glazing of partially crystallized IPS e.max® CAD restorations (in their blue state).

Processing restrictions

The product is not intended to be reprocessed or reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

Restoration of esthetics

Composition

- **IPS Ivocolor Shade**
Components: Alkali aluminosilicate glass, 1,3-butandiol, glycerine, water, pigments
- **IPS Ivocolor Essence**
Components: Alkali aluminosilicate glass, pigments
- **IPS Ivocolor Glaze Powder/FLUO**
Components: Alkali aluminosilicate glass
- **IPS Ivocolor Glaze Paste/FLUO**
Components: Alkali aluminosilicate glass, 1,3-butandiol, glycerine, water
- **IPS Ivocolor Mixing Liquids/Essence Fluid**
Components: 1,3-butandiol, 1,5-pentandiol, glycerine, water

2 Application

- **IPS Ivocolor Shade**

IPS Ivocolor Shade Dentin and Shade Incisal are intended for the superficial shading of the restoration to achieve the desired A-D tooth shade.

The ready-to-use pastes have a gel-like consistency and need to be mixed / diluted with IPS Ivocolor Mixing Liquids if required. Slightly wetting the restoration surface with IPS Mixing Liquid is recommended prior to applying the shades, as this will make application easier. Apply the shades in a thin coat.

- **IPS Ivocolor Essence**

IPS Ivocolor Essence materials are intended for the superficial individual characterization of restorations as well as for the internal characterization of layering ceramics.

The Essence powders can either be used by themselves or they can be mixed with each other as desired. For initial mixing, IPS Ivocolor Essence Fluid can be used. Apply the Essence materials in a thin coat.

- **IPS Ivocolor Glaze Powder/FLUO**

IPS Ivocolor Glaze Powder/FLUO materials are used for the glazing of restorations.

Mix the powders to the desired consistency using IPS Ivocolor Mixing Liquids. For this purpose, add the mixing liquid dropwise to the glaze powder and mix until a homogeneous paste is achieved. If the glaze is correctly mixed, its consistency should not be too thin, so as to ensure an appropriate lustre of the fired restoration. Apply the glaze in a thin coat. For zirconium oxide restorations, we recommend applying the glaze in a more viscous, i.e. less diluted, consistency.

- **IPS Ivocolor Glaze Paste/FLUO**

IPS Ivocolor Glaze Paste/FLUO materials are used for the glazing of restorations.

The gel-like pastes are mixed / diluted to the desired consistency using IPS Ivocolor Mixing Liquids. If the glaze is correctly mixed, its consistency should not be too thin, so as to ensure an appropriate lustre of the fired restoration. Apply the glaze in a thin coat. For zirconium oxide restorations, we recommend using the glaze in a more viscous, i.e. less diluted, consistency.

- **IPS Ivocolor Mixing Liquid allround**

Liquid for mixing all IPS Ivocolor powders and pastes

- **IPS Ivocolor Mixing Liquid longlife**

Liquid for mixing all IPS Ivocolor powders and pastes. Due to the viscous consistency and slow evaporation of the liquid, the mixed powders and pastes feature a higher viscosity and offer longer working time.

- **IPS Ivocolor Essence Fluid**

Liquid for the mixing of IPS Ivocolor Essence powders ensuring a paste-like consistency.

IPS Ivocolor Essence Fluid is only suitable for the initial mixing of IPS Ivocolor Essence powders.

Notes

- IPS Ivocolor is only approved for use in conjunction with Ivoclar IPS ceramics. If IPS Ivocolor is used for other ceramics, the sole responsibility lies with the user.
-  The proper preparation of the surface is dependent on the material used for the restoration. Proceed according to the instructions for use of the material in use.
- Before applying stains and glazes, the restoration must be free of grease and dirt.
- In order to remove any residues, clean the restoration with a steam cleaner or in an ultrasonic bath.
- Use the IPS Ivocolor components as indicated in the combination table.
- Avoid pooling and the application of inadequately thick layers.
- The consistency used affects the application behaviour and the firing result. After correct application of the glaze material, the restoration will show a uniform lustre and surface texture.
- More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers.
- Apply IPS Ivocolor stains and glazes with a brush to achieve the desired effect. The products can be used for staining and glazing, for glazing only or for the internal characterization of layering ceramics. After applying the stains and glazes, a glaze and stain firing cycle must be conducted
- The IPS Ivocolor Shade Guide is used to select the correct Ivocolor Shade and Essence colour.

After finalization, the dentist will incorporate the restoration in the mouth of the patient based on the instructions of the manufacturer.

Components Combination Table

Components	Initial mixing			Touching up dried out pastes		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

IPS Ivocolor® Shade

Combination table for the shading of IPS veneering ceramics and the all-ceramic materials of the IPS e.max® system*

A-D shade	Shade 0	Shade 1	Shade 2	Shade 3	Shade 4	Shade 5	Shade 6	Shade 7	Shade 8
D4							●		
D3								●	
D2								●	
C4							●		
C3							●		
C2							●		
C1							●		
B4						●			
B3						●			
B2					●				
B1					●				
A4				●					
A3.5			●						
A3			●						
A2			●						
A1		●							
BL4	●								
BL3	●								
BL2	●								
BL1	●								

* When using IPS e.max ZirCAD LT discs in LT group shades, other shade combinations apply.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/elfu
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
Basic UDI-DI: 76152082ACERA009F8

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.

4 Shelf life and storage

- Storage temperature:
 - Shade/Paste/Mixing Liquids/Essence Fluid: 2–28°C (36–82 °F).
 - Powder: No special storage conditions are required.
- Immediately seal the container tightly after each use to prevent any adverse environmental effects on the product.
- Do not use the product after the indicated expiry date.
- Expiry date: See note on the packages.

5 Additional information

Keep material out of the reach of children.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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- IPS Ivocolor Mixing Liquid allround

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- IPS Ivocolor Mixing Liquid longlife

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Componenti	Prima miscelazione			Ritocco di paste asciugate		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Componentes	Initial mixing			Touching up dried out pastes		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Komponenter	Initial blandning			Förbättring av torra pastor		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Färg	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-

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Komponentkombinationstabel

Komponenter	Første blanding			Opfriskning af udtørrede pastaer		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Komponentit	Alkusekoitus			Kuivuneiden tahnojen notkistaminen		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Komponenter	Første gangs blanding			Oppfrisking av inntørkede materialer			
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-	
IPS Ivocolor Shade	✓	✓	-	✓	✓	-	
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-	
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-	

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Combinatietabel componenten

Componenten	Initieel mengen			Uitgedroogde pasta's bijwerken		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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IPS Ivocolor Mixing Liquid allround

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IPS Ivocolor Mixing Liquid longlife

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Bileşenler	İlk karıştırma			Kurumuş macunlara rötüş		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Компоненты	Первичное замешивание			Касание высохших паст		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	–
IPS Ivocolor Shade	✓	✓	–	✓	✓	–
IPS Ivocolor Glaze Powder / FLUO	✓	✓	–	✓	✓	–
IPS Ivocolor Glaze Paste / FLUO	✓	✓	–	✓	✓	–

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Tabela kombinacij komponent

Komponente	Začetno mešanje			Popravljanje izsušenih past		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Tablica s kombinacijama komponenti

Komponente	Početno miješanje			Dotjerivanje osušenih pasta		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Компоненти	Първоначално смесване			Освежаване на засъхнали пасти		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Komponentët	Përzjerja fillestare			Rregullimi i pastave të thara		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO.	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Componente	Primul amestec			Retușarea umidității pastelor uscate		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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Komponendid	Algne segamine			Kuivanud pasta retuššimine		
	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid	IPS Ivocolor Mixing Liquid allround	IPS Ivocolor Mixing Liquid longlife	IPS Ivocolor Essence Fluid
IPS Ivocolor Essence	✓	✓	✓	✓	✓	-
IPS Ivocolor Shade	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Powder/FLUO	✓	✓	-	✓	✓	-
IPS Ivocolor Glaze Paste/FLUO	✓	✓	-	✓	✓	-

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