Molloplast® B

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Instructions for use
with cleaning recommendations for the patient

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**permanently soft
+ anti-plaque effect**

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**Product description:**

Heat curing one-component silicone, permanently soft relining material, suitable for denture resins based on (meth-)acrylate.

Indications for use:

Sharp-ridged alveolar processes with thin tegument. Cushioning of dentures against the sharply marked linea mylohyoidea, foramen mentale or other prominent points of the alveolar process such as trigonum retromolare or tubera. Cushioning against the sutura mediana. For post-damming of the A-line and to steady the position of the upper and lower dentures by increased adhesion. For the construction of clasless dentures and for covering major surgical defects with obturators.

Contra-indication:

Do not apply to patients suffering from a wellknown hypersensitivity to organic peroxides.

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Silicona de componente único polimerizada por calor, material de rebase permanentemente blando, apto para resinas para dentaduras a base de (met-)acrilato.

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Molloplast® B

Instructions for use

Relining of new dentures with Molloplast® B

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Spacer

is fabricated with a thermo-vacuum apparatus. Thickness of the spacer is to be determined by the dentist and should be between 1 - 3 mm. If a thermo-vacuum unit is not available, use cold curing acrylic and trim to the requested thickness (Fig. 2). Set up teeth and try-in with spacer in patient's mouth (Fig. 3).

Investing

After the try-in invest the wax pattern with spacer into the flask with counter. An insulation of the teeth by means of FLEXISTONE® or FLEXISTONE® *Plus* is recommended. Place flask in boiling water for approx. 4 min., then open flask and remove all traces of wax.

Important:

Wax residues must not penetrate or adhere to plaster.

Insulation

Apply an alginate-based insulating liquid and reposition the spacer.

Acrylic mixing

Acrylic should be mixed according to instructions and packed into a lukewarm flask. Prior to pressing, place a polyethylene (PE-)foil between acrylic and spacer. (Fig. 4).

Pressing

Pre-press in intervals of 10-15 min. at 100 kp. Open the flask and remove PE-foil and excess acrylic. Close flask with spacer and PE-foil and put it into a clamp. Screw clamp tight, place flask in cold water and bring it to the boil. Let boil for approx. 30 min., then cool down flask. This pre-polymerisation prevents possible reactions with the acrylic monomer.

Packing of Molloplast® B

Remove flask from the clamp, open and remove spacer and PE-foil.

Important:

The acrylic must show adequate stiffness to avoid deformation by the Molloplast® B material during pressing (Fig. 5). Molloplast® B needs no mixing. Take Molloplast® B with a clean plastic spatula from the jar and pack it on the pre-pressed acrylic in the flask. Place a new PE-foil between Molloplast® B and counter of the flask. Close flask and pre-press in intervals at 100 kp. Open flask and remove PE-foil and excess Molloplast® B material. Close flask again and press for approx. 10 - 15 min. at 100 kp. For polymerisation place flask in the clamp.

Polymerisation

Place flask in cold water and heat up slowly to 100 °C / 212 °F.

Polymerisation in boiling water at 100 °C/212 °F for approx. 2 hours.

Cool down flask slowly. Do not "shocktreat" by rinsing cold water over the flask!

Polymerisation of Molloplast® B by microwaving

Curing of Molloplast® B is recommended in a commercially available microwave flask at 850 watts for 10 min., using microwave suitable dental plaster for investing. Investing and polymerisation must be carried out within 8 hours. Ovens without rotating plate may require a gradual, all-side polymerisation by turning the flask. De-flask after self-cooling to room temperature.

Trimming of the denture

Use normal procedure for trimming the acrylic. The relined denture base and the functional borders can be trimmed with Molloplast® special cutters or grinding sleeves (15 - 20 thousand r.p.m.), without problems and heat build-up. Use Molloplast® pre-polishing discs for smoothing of rough spots (Fig. 6).

Polishing

Use standard methods for pre-polishing with pumice, brush or buff and polish with commercially available materials. Mechanical high gloss polishing of Molloplast® B relined dentures is not possible -for high gloss polishing Lustrol gloss varnish is recommended, mixing ratio: 1:1. Prior to the application of Lustrol gloss varnish base and functional borders of the

relined denture have to be absolutely clean and dry. Apply Lustrol with a brush in a thin layer and let it dry for approx. 10 - 15 min. - resulting in a high lustre polish (Fig. 7). The dentist can now insert the soft-lined denture into the patient's mouth.

First relining of existing dentures with Molloplast® B

For relining of an existing denture the dentist takes an impression with the existing denture, using a medium or light bodied impression material (Fig. 8).

Plaster cast

A master cast of dental stone class III is poured. Let the master cast set and invest it into a flask with counter. For easier removal of the denture from the plaster counter it is recommended to coat the teeth and also a part of the denture surface with silicone (FLEXISTONE® or FLEXISTONE®Plus), prior to making the flask counter. Allow plaster to set completely, open the flask and remove the impression material from the denture.

Denture base

Reduce acrylic material for approx. 1 - 3 mm by means of a coarse grinding stone.

Important:

Provide sufficient retentions to ensure optimum bonding with Molloplast® B. For a smooth surface between acrylic and Molloplast® B we recommend to cut grooves evenly into the functional borders by means of a wheel or fissure bur. Afterwards brush Primo adhesive 1 - 2 times onto the entire denture base. Let Primo air-dry for 60 - 90 min. prior to applying Molloplast® B (Fig. 9).

Insulation

Insulate master cast with an alginate-based insulating liquid.

Packing of Molloplast® B

Take ready-to-use Molloplast® B with a clean plastic spatula from the jar and apply it evenly

onto the denture base. Cover with a PE-foil, close the flask and place it under the press. The flask must be lukewarm.

Pressing

Press for approx. 4 min., open the flask and remove PE-foil and Molloplast® B excess material. Control and if necessary apply more Molloplast® B material. Close the flask and place it under the press for approx. 15 min. at 100. Place flask in the clamp and tighten.

Polymerisation

Place flask in cold water and heat up slowly to 100 °C / 212 °F.

Polymerisation in boiling water at 100 °C / 212 °F for approx. 2 hours.

Cool down flask slowly. Do not "shocktreat" by

rinsing cold water over the flask!

Polymerisation of Molloplast® B by microwaving.

Curing of Molloplast® B is recommended in a commercially available microwave flask at 850 watts for 10 min., using microwave suitable dental plaster for investing. Investing and polymerisation must be carried out within 8 hours. Ovens without rotating plate may require a gradual, all-side polymerisation by turning the flask. De-flask after self-cooling to room temperature.

Trimming of the denture

Use normal procedure for trimming the acrylic. The relined denture base and the functional borders can be trimmed with Molloplast® special cutters or grinding sleeves (15 - 20 thousand r.p.m.), without problems and heat build-up.

Use Molloplast® pre-polishing discs for smoothing of rough spots (Fig. 6).

Polishing

Use standard methods for pre-polishing with pumice, brush or buff, and polish with commercially available materials. Mechanical high gloss polishing of Molloplast® B relined dentures is not possible - for high gloss polishing Lustrol gloss varnish is recommended, mixing ratio: 1:1. Prior to the application of Lustrol gloss varnish base and functional borders of the relined denture have to be absolutely clean and dry. Apply Lustrol with a brush in a thin layer and let it dry for approx. 10 - 15 min. resulting in a high lustre polish (Fig. 7).

Primo adhesive

Primo adhesive ensures optimum bonding of Molloplast® B to the acrylic denture base. However, the adhesive is only applied onto already polymerised denture base material. Prior to the application of Primo adhesive the denture base is prepared according to Molloplast® B instructions for use. Then Primo adhesive is applied uniformly with a brush 1 - 2 times onto the clean and dried denture base. Please note: Primo coating must cover the entire surface of denture acrylic in contact with Molloplast® B. After drying for approx. 60 - 90 min. Molloplast® B is pressed onto the denture and polymerised according to instructions.

General information

Wear protective gloves!

Only for dental use by qualified personnel.

Primo adhesive: avoid contact with eyes, mucous membranes, and skin.

Lustrol gloss varnish: avoid contact of catalyst with eyes and skin, this can lead to irritation. In case of accidental contact, rinse eyes immediately with running water for at least 15 minutes, consult an ophthalmologist.

Please follow the instructions of the safety data sheets!

Special recommendations to avoid incorrect processing

1. Air inclusions in polymerised dentures:

- Not enough pressure while pressing.
- Flask is placed into hot water.
- Packing of Molloplast® B before acrylic is stiff.
- Wetting of the acrylic with monomer liquid prior to packing Molloplast® B.

2. Insufficient or no bonding between Molloplast® B and acrylic:

- Not enough retentions have been cut into the existing denture base or Primo adhesive has not been applied.

- b) The acrylic denture has been wetted with alginate-based insulating liquid.
- c) Hands and material were not free of greasy spots.
- d) Excessive heat build-up during trimming may cause separation of Molloplast® B from functional borders. Use Molloplast® cutters and prepolishers only.

Recommendations for repairs of Molloplast® B

Broken dentures relined with Molloplast® B can be repaired with Molloplast® B. Partial supplementations of Molloplast® B relinings are also possible. Invest the denture as usual. Grind off approx. 1 mm of the old relining material, using Molloplast®-cutters. In case of breakage open the area to be repaired V-shaped, beginning from the base. Apply new Molloplast® B supplementary material on the repairing area of the denture invested in the flask and put it under the press. (Test pressing with foil is possible). Let polymerise according to instructions.

Do not use Primo adhesive!

Hygienic care and cleaning of your permanently soft relined denture

For the patient
Please make a copy and give it to the patient

Your new Molloplast® B relining has many advantages: It resembles natural gingiva, is comfortable to wear and is resistant to mechanical pressure. Please read this hand-out carefully and follow the guidelines to maximise the use of your relined denture:

- Remove the denture after each meal and rinse it well using tap water
- Rinse also your mouth carefully, especially after drinking fruit or vegetable juices.

- At least once daily the denture must be thoroughly cleaned, removing all plaque. This is best done using a liquid soap and a soft tooth-brush.
- Also commercially available denture cleaning liquid may be used - however do not immerse the denture for more than 10-15 min. to avoid discoloration.
- At night the relined denture should be removed, cleaned and stored in a dry place.

- Frequent consumption of tea, coffee or wine may cause discoloration of the denture. Special note to smokers: Tobacco-containing products may cause irreversible discoloration of the relined denture!

Additional care with “good morning spray”
good morning spray improves the oral hygiene. The appearance of plaque from oral-specific bacteria is delayed. After thorough morning cleaning of the denture, moisten the relined denture base by

applying 2-3 spray strokes, then insert the denture immediately.

Please note: the spray application does not replace any cleaning procedure of the relined denture.

We advise you to clean your teeth plaque-free according to your dentist's instructions - it is for your own benefit.

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