HeraCeram Zirkonia 750

(DE)	Gebrauchsanweisung	2
(GB)	Instructions for use	4
(FR)	Mode d'emploi	6
(ES)	Instrucciones de uso	8
(IT)	Istruzioni per l'uso	10
(PT)	Instruções de uso	12
(NL)	Gebruiksaanwijzing	14
(SE)	Bruksanvisning	16
(DK)	Brugervejledning	18
(NO)	Bruksanvisning	20
(FI)	Käyttöohjeet	22
(CZ)	Návod k použití	24
(HU)	Használati utasítás	26
(LV)	Lietošanas instrukcija	28
(LT)	Naudojimo instrukcija	30
(PL)	Instrukcja obsługi	32
(HR)	Upute za uporabu	34
(RO)	Instrucțiuni de utilizare	36
(RU)	Инструкция по применению	38
(SK)	Návod na použitie	40
(BG)	Инструкции за употреба	42
(UA)	Інструкція по застосуванню	44
(EE)	Kasutusjuhend	46

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HeraCeram Zirkonia 750 is the ideal veneering material for ceramic frameworks manufactured from zirconium dioxide and lithium di-silicate, with a CTE ($\alpha_{25-500^{\circ}\text{C}}$) of 10.2–10.5 $\mu\text{m/mK}$.

Veneering of zirconium dioxide frameworks

Framework preparation: A very thin/even layer of Adhesive 750 is applied and subsequently fired in the ceramic furnace under vacuum. (See firing procedures at the table below). This will ensure optimal adhesion between the zirconium dioxide framework and HeraCeram Zirkonia 750.

Increaser: For further shade adaptation, the veneering surface is covered with a layer of increaser using the appropriate colour.

Dentine/Incisal layers: The ceramic can be applied to suit your build up technique depending upon preference. For example a simplified two-layer technique or, customised/individualised approach.

Note: Wear a face mask and eye protection when grinding ceramic materials and work with a suction device. Avoid breathing in the ceramic dust.

Glaze firing: The degree of glazing is defined by the final firing temperature and the amount of holding time. For this reason, the specifications for the glaze firing temperature and holding time are for orientation only and may need to be adjusted to achieve the required result. Glazing is often down to personal preference therefore is also possible to polish HeraCeram Zirkonia 750 mechanically. Our Signum® HP pastes have proven themselves ideal for final polishing.

Veneering of lithium di-silicate frameworks

Framework preparation: Do not blast the lithium di-silicate frameworks after finishing; simply clean them carefully (e.g., steam cleaner).

Adhesive: A thin layer of Adhesive 750 is then applied to the surface of the framework and fired.

Dentine/Incisal layering: See above. Please pay attention to the firing procedure

Glaze firing: See above. Please pay attention to the firing procedure

Colour mapping table

Shade	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Increaser	IN A1	IN A2	IN A3	IN A3,5	IN A4	IN B1	IN B2	IN B3	IN B4	IN C1	IN C2	IN C3	IN C4	IN D2	IN D3	IN D4
Dentine	D A1	D A2	D A3	D A3,5	D A4	D B1	D B2	D B3	D B4	D C1	D C2	D C3	D C4	D D2	D D3	D D4
Chromadentine	CD A1	CD A2	CD A3	CD A3,5	CD A4	CD B1	CD B2	CD B3	CD B4	CD C1	CD C2	CD C3	CD C4	CD D2	CD D3	CD D4
Incisal	S1	S1	S2	S2	S4	S1	S1	S2	S4	S1	S3	S3	S3	S1	S2	S2
Enhancer*	EHA	EHA	EHA	EHA	EHA	EHB	EHB	EHB	EHB	EHC	EHC	EHC	EHC	EHA	EHA	EHA
MD/SD	MD1	MD1	SD2	SD2	SD2	MD2	MD2	MD3	MD3	MD2	SD1	SD2	SD2	MD1	MD3	SD1
Value	VL1	VL2	VL3	VL4	VL4	VL1	VL2	VL3	VL4	VL1	VL2	VL3	VL4	VL2	VL3	VL4
Opal incisal*	OS1	OS1	OS2	OS2	OS4	OS1	OS1	OS2	OS4	OS1	OS3	OS3	OS3	OS1	OS2	OS2
Opaltranspa*	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT	OT
Margin HM/LM	1	1	2	2	6	3	3	4	4	1	1	6	6	1	2	4

* The use of the individual materials should always reflect the patient requirements and the specifications in this table are to be understood as non-binding recommendations.

Firing procedure

IMPORTANT NOTICE: The firing temperatures given are intended as guidelines. Deviations are possible due to varying furnace power and types and may need to be adapted.

General firing table zirconium dioxide		Adhesive 750	Primary dentine	Secondary dentine	Glaze	Correction material
Start temperature	°C	500	500	500	500*	500
Drying	min	5	5	5	4	4
Pre-heating	min	1	1	1	1	1
Temperature increase	°C/min	75	60	60	60	60
Final temperature	°C	800	750	745	725	720
Holding time	min	2	1	1	1	1
Vacuum start	°C	500	500	500	–	500
Vacuum stop	°C	800	750	745	–	720

* Start temperature 600°C for stain shade and glaze pastes!

General firing table lithium di-silicate		Adhesive 750	Primary dentine	Secondary dentine	Glaze	Correction material
Start temperature	°C	400	400	400	400*	400
Drying	min	5	5	5	4	4
Pre-heating	min	1	1	1	1	1
Temperature increase	°C/min	50	50	50	60	50
Final temperature	°C	750	750	745	725	715
Holding time	min	1	1	1	1	1
Vacuum start	°C	400	400	400	–	400
Vacuum stop	°C	750	750	745	–	715

* Start temperature 600°C for stain shade and glaze pastes!

Safety instructions

The data given are furnished to the best of our knowledge and are based on in-house test results. To be used only by dental technicians for its intended purpose. This product or one of its components may, in particular cases, cause hypersensitive reactions. If any reaction is suspected, information about the ingredients can be obtained from the manufacturer. The use of the product is contraindicated in cases of known or suspected allergies to components of the ceramic. In the case of a suspected allergy, we recommend taking an allergy test prior to use. Never swallow or ingest this product! Avoid contact with eyes. In case of contact with eyes, rinse immediately and thoroughly with plenty of cool water and seek medical advice.

Please observe the information in the corresponding safety data sheet.

Please quote the batch number in all correspondence about the product.

Batch number LOT: See note on packaging.

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Kulzer GmbH
Kulzer, LLC