BIO-C[®] REPAIR

CHOOSE YOUR LANGUAGE

ENGLISH.....	3	HRVATSKI	73
ESPAÑOL.....	8	SLOVENŠČINA	79
PORTUGUÊS.....	14	SLOVENSKY	84
FRANÇAIS	20	POLSKI	90
ITALIANO	26	EESTI KEEL.....	96
ROMÂNĂ	32	LIETUVIŠKAI	102
NEDERLANDS.....	38	ČESKY	107
DEUTSCH	44	MAGYAR	113
SVENSKA	50	TÜRKÇE	119
NORSK.....	56	БЪЛГАРСКИ	125
DANSKI	61	СРПСКИ.....	131
SUOMI.....	67	ΕΛΛΗΝΙΚΑ	137

ENGLISH

DEFINITION

Endodontic bioceramic reparative ready-to-use, intended to be used for patients requiring endodontic procedures, with the function of repairing endodontic complications in the tooth's root canal.

INDICATIONS

- A. Treatment of perforation of root or furcation by canal;
- B. Treatment of perforation of root or furcation surgically;
- C. Treatment of internal resorption surgically or by root canal;
- D. Treatment of external resorption;
- E. Retrofilling in paraendodontic surgery;
- F. Direct and indirect pulp capping;
- G. Apexification;
- H. Apicogenesis and Pulpotomy.

INTENDED PATIENTS' POPULATION

BIO-C® REPAIR should be utilized in cases involving patients who need treatment for lateral and furcation root canal perforations, internal resorption, reverse root filling, pulp capping, apexification, and apicogenesis.

INTENDED USER

BIO-C® REPAIR is indicated for use by dentists who are trained to perform endodontic procedures.

INTENDED ENVIROMENT

Dental offices.

CONTRAINDICATIONS

BIO-C® REPAIR contains components that may cause reactions in susceptible individuals. Do not use on patients who are allergic to the product's components.

COMPOSITION

Calcium silicate, calcium aluminate, calcium oxide, zirconium oxide, iron oxide, silicon dioxide and dispersing agent.

TECHNICAL INFORMATION

- Working time: ready-to-use product, no mixing required.
- Setting time: 30 - 120 minutes, depending on local humidity.

Important: dry with paper cones, but without excessively drying the cavity. The humidity released from the tissues will initiate the material's setting reaction.

- High alkalinity.
- High radiopacity.

WARNING

- Avoid contact with eyes or skin. In case of contact, wash immediately with water.
- Avoid contact with the oral mucosa. In case of contact, wash with water and avoid ingestion of the product. If any sensitivity persists, seek medical attention immediately.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority health in which the user and/or patient is established

PRECAUTIONS

- Wear protective eyewear, mask and gloves when handling the product. In case of contact with eyes or skin, wash immediately with water.
- Use absolute isolation in all indications for use via root canal.
- Before use, check if the product has hardened on the syringe tip and, if so, discard the first portion.
- Use the product only after controlling the infection or the material will not set due to the acidic pH of the environment.
- The cavities where the product will be inserted do not need to be completely dry, but they should not be excessively damp.
- Do not use the product to fill root canals, as its flowability is inadequate for this purpose.
- To avoid product hardening, close the syringe immediately after applying the material using adequate pressure.
- Never store the product in a refrigerator, as it is sensitive to humidity.
- Do not use the product after its expiration date.
- Do not use the product if its packaging is damaged.
- The safety and efficacy of this product have not been established in pregnant, lactating women or for children and adolescents.
- Keep out of reach of children.

SYRINGE

After use, clean and disinfect the syringe for future use and confirm that the silica sachet inside the envelope remains blue. After this procedure, return the properly closed syringe together with the silica sachet inside the package, which must be fully closed and stored correctly.

CAPSULE

The BIO-C® REPAIR single dose capsule is not reusable. After use, dispose of the capsule and envelope according to the "DISPOSAL METHOD" section.

TECHNIQUE FOR USE

CAPSULE

1. Before applying BIO-C® REPAIR single dose, check the product for any damage or defects.
2. Insert the capsule into the gun and remove the capsule cap.
3. Squeeze the gun with a slow, even pressure to extrude the correct amount of material from the capsule.
4. Depending on the application site, BIO-C® REPAIR single dose can be applied with the gun itself or with an appropriate instrument.
5. Insert the material into the treated area.
6. Complete the treatment according to the appropriate methods.

Make sure the capsule is at the end of the gun for easier removal. Remove the capsule from the gun and dispose of it after use in accordance with national laws and local authority recommendations for proper disposal. Do not reuse the capsule on another patient. Sterilize the gun using the appropriate method.

A. TREATMENT OF PERFORATION OF ROOT OR FURCATION BY CANAL

1. Anesthetize, install the absolute isolation, execute the biomechanical preparation of the duct and cleaning of the existing perforation;
2. Execute hemostatic control and maintain the cavity slightly moist;
3. Insert BIO-C® REPAIR and, if there is overflow, remove the excess material;
4. X-ray to check the correct seating of the material in the place applied;
5. Fill in the remainder of the duct with gutta-percha cones and a filling cement;
6. Execute the coronary sealing and restoration.

Important:

- In perforations of furcation apply a layer of glass ionomer as a mechanical barrier followed by restoration with final material;
- Do not use composite resin directly on BIO-C® REPAIR before its final setting;
- X-ray and execute preservation for at least two years.

B. TREATMENT OF PERFORATION OF ROOT OR FURCATION SURGICALLY

1. Anesthetize, raise the gingival flap and execute osteotomy to locate the perforation;
2. Prepare the perforation with an appropriate instrument to facilitate the insertion of the material;
3. Execute hemostatic control and maintain the cavity slightly moist;

4. Insert BIO-C® REPAIR in the cavity, adapt with condensers and remove excess if overflowing;
5. Reposition the gingival flap and sew it;
6. X-ray and execute preservation for at least two years.

C. TREATMENT OF INTERNAL RESORPTION SURGICALLY OR BY ROOT CANAL

1. Anesthetize and install the absolute isolation;
2. Remove the granulation tissue of the reabsorption area;
3. Neutralization of the medium with calcium hydroxide paste;
4. Remove the calcium hydroxide paste in the following session and execute endodontical treatment in the conventional way until the height of the reabsorption;
5. Apply BIO-C® REPAIR in the place of reabsorption, condensing it against the walls with specific instruments or a sterile ball of cotton slightly moistened;
6. Execute the coronary sealing with glass ionomer and restoration with the material preferred;
7. X-ray and execute preservation for at least two years.

D. TREATMENT OF EXTERNAL RESORPTION

1. Anesthetize, raise the gingival flap and execute osteotomy to locate the perforation;
2. Remove all the cement and dentin affected, removing the cells inducing reabsorption;
3. Execute hemostatic control and dry the cavity. It is advisable to condition the surface of the root affected, preferably with citric acid;
4. Insert BIO-C® REPAIR in the cavity which must be slightly moist. Adapt the material with condensers and remove the excess if overflowing;
5. Reposition the gingival flap and sew it;
6. X-ray and execute preservation for at least two years.

Important: If there is external reabsorption above the bone crest it is not advisable to use BIO-C® REPAIR.

E. RETROFILLING IN PARAENDODONTIC SURGERY

1. Anesthetize, raise the gingival flap and execute osteotomy to expose the dental apex;
2. Cut 3 mm of the root apex which contains a large number of collateral canals (apical delta);
3. Execute retro preparation with appropriate instrument;
4. Execute hemostatic control and maintain the cavity slightly moist;
5. Insert the BIO-C® REPAIR in the retrocavity, adapt it with appropriate condensers and remove any excess extravasated material;
6. Return the gingival flap to its position and sew it;
7. X-ray and execute preservation for at least two years.

F. DIRECT AND INDIRECT PULP CAPPING

1. Anesthetize and install the absolute isolation;
2. Remove the caries; If there is pulp exposure execute hemostasis;
3. Apply BIO-C® REPAIR on the exposed pulp with suitable instruments and tamp it with a slightly moistened sterile cotton ball;
4. Recover BIO-C® REPAIR with a glass ionomer cement and wait 4 to 6 weeks to execute the final restoration. The glass ionomer can be used as the basis of the restoration.

Warning: for Young permanent or deciduous teeth with a great exposure, execute the pulpotomy and then do steps 3 and 4 above.

G. APEXIFICATION

First Session

1. Anesthetize and install the absolute isolation;
2. Execute the preparation of the duct by the conventional technique;
3. Apply an intracanal dressing with calcium hydroxide paste for a week.

Second Session

1. Remove the intracanal dressing with files and irrigating solution;
2. Dry the duct with paper cones;
3. Make a mechanical barrier by applying a collagen sponge in the apical region, avoiding any BIO-C® REPAIR leakage;
4. Apply BIO-C® REPAIR condensing it to form an apical cover of 3 to 4 mm;
5. X-ray to check the correct filling of the canal;
6. Place a sterile cotton ball slightly moistened at the entrance of the canal followed by provisional restoration.

Third Session

1. Remove the provisional restoration and fill the canal to the limit of the apical cover;
2. Restore with glass ionomer;
3. Execute the clinical and radiographic control of 3 to 6 months until the radiographic check of apical barrier of hard tissue.

H. APICIGENESIS AND PULPOTOMY

1. Anesthetize, install the absolute isolation and then make the access to the pulp chamber;
2. Remove the roof of the pulp chamber and execute excision of the coronary pulp with sharpened curettes or drills;

3. Cut the pulp 0.5 mm below the inlet of the root canal;
4. Wash the surgical wound abundantly with distilled water, saline solution or other appropriate sterile liquid used for the same purpose;
5. Execute hemostasis and maintain the cavity slightly moist;
6. Apply BIO-C® REPAIR on the pulp stump and adapt it with a sterile cotton ball slightly moistened;
7. X-ray to check the correct filling of the canal;
8. Restore the cavity provisionally with glass ionomer;
9. Execute the clinical and radiographic control of 3 to 6 months until the radiographic check of radicular formation;
10. After the end of the radicular formation choose between conventional endodontic treat or only the restoration of the coronary endodontic cavity.

STORAGE

- Keep the product in a dry and well ventilated location between 15 and 30°C with relative humidity below 60%;
- Do not store in a refrigerator;
- Do not store the product near ammonia, ammonium nitrate and products containing chlorine. Avoid the use of disinfecting solutions containing any of these ingredients.

DISPOSAL METHOD

The product must be disposed of in accordance with local regulations.

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